USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

1. Technical data

Parameter description	Parameter value				
Product name	Platform scale				
Model	SBS-BW-3T	SBS-BW-5T	SBS-BW-1T	SBS-BW-3T/1KG	SBS-BW-5T/2KG
Measuring range [kg]	up to 3000	up to 5000	up to 1000	up to 3000	up to 5000
Division [kg]	1	2	0.5	1	2
Dead weight [kg]	119	151	79	57	82
Weighing surface [cm]	150 150		100 100		
Functions	Hold				
	Tare				
	Zero				
	Parts counting				
	Change the weight units (kg/lb)				
	PACK				
Display	LCD				

Parameter description	Parameter value		
Product name	Floor scale		
Model	SBS-BW-1T/0.5A	SBS-BW-3T/1A	SBS-BW-5T/2A
Measuring range [kg]	up to 1000	up to 3000	up to 5000
Division [kg]	0.5	1	2
Dead weight [kg]	52	99	98.5
Weighing surface [cm]	100 100	150 151.5	
Functions	Hold		
	Tare		
	Zero		
	Parts counting		
	Change the weight units (kg/lb)		
	PACK		

	PACK
Display	LCD

Parameter description	Parameter value	
Product name	Animal scale	
Model	SBS-PS-100	SBS-PS-200
Measuring range [kg]	up to 1500	up to 5000
Division [kg]	2	
Dead weight [kg]	76	146
Weighing surface [cm]	100x200	100x220
Functions	Hold	
	Tare	
	Zero	
	Parts counting	
	Change the weight units (kg/lb)	
	PACK	
Display	LCD	

2. General description

The user manual is designed to assist in the safe and trouble-free use of the device. The product is designed and manufactured in accordance with strict technical guidelines, using state-of-the-art technologies and components. Additionally, it is produced in compliance with the most stringent quality standards.

DO NOT USE THE DEVICE UNLESS YOU HAVE THOROUGHLY READ AND UNDERSTOOD THIS USER MANUAL.

To increase the product life of the device and to ensure trouble-free operation, use it in accordance with this user manual and regularly perform maintenance tasks. The technical data and specifications in this user manual are up to date. The manufacturer reserves the right to make changes associated with quality improvement. The device is designed to reduce noise emission risks to a minimum, taking into account technological progress and noise reduction opportunities.

2.1. Legend

Icon	Description
	The product satisfies the relevant safety standards.
	Read instructions before use.
	The product must be recycled.

	WARNING! or CAUTION! or REMEMBER! Applicable to the given situation. (general warning sign)
	ATTENTION! Electric shock warning!
	Class II protection device with double insulation.

PLEASE NOTE! DRAWINGS IN THIS MANUAL ARE FOR ILLUSTRATION PURPOSES ONLY AND IN SOME DETAILS MAY DIFFER FROM THE ACTUAL PRODUCT.

3. Usage safety

ATTENTION! READ ALL SAFETY WARNINGS AND ALL INSTRUCTIONS. FAILURE TO FOLLOW THE WARNINGS AND INSTRUCTIONS MAY RESULT IN ELECTRIC SHOCK, FIRE AND/OR SERIOUS INJURY OR EVEN DEATH.

The terms "device" or "product" are used in the warnings and instructions to refer to:

Platform scale / Floor scale / Animal scale

3.1. Electrical safety

For the battery charger

- a) The plug must fit the socket. Do not modify the plug in any way. Using original plugs and matching sockets reduces the risk of electric shock.
- b) Avoid touching earthed elements such as pipes, heaters, boilers and refrigerators. There is an increased risk of electric shock if the earthed device is exposed to rain, comes into direct contact with a wet surface or is operating in a damp environment. Water getting into the device increases the risk of damage to the device and of electric shock.
- c) Do not touch the device with wet or damp hands.
- d) Use the cable only for its designated use. Never use it to carry the device or to pull the plug out of a socket. Keep the cable away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cables increase the risk of electric shock.
- e) If using the device in a damp environment cannot be avoided, a residual current device (RCD) should be applied. The use of an RCD reduces the risk of electric shock.
- f) Do not use the device if the power cord is damaged or shows obvious signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service centre.

- g) ATTENTION! DANGER TO LIFE! While cleaning, never immerse the device in water or other liquids.

3.2. Safety in the workplace

- a) Make sure the workplace is orderly and well lit. A messy or poorly-lit workplace may lead to accidents. Try to think ahead, observe what is going on and use common sense when working with the device.
- b) If there are any doubts as to the correct operation of the device, contact the manufacturer's support service.
- c) Only the manufacturer's service point may repair the device. Do not attempt any repairs independently!
- d) If a fire starts, use solely powder or carbon dioxide (CO₂) fire extinguishers suitable for use on live devices to put it out.
- e) Please keep this manual available for future reference. If this device is passed on to third parties, the manual must be passed on with it.

REMEMBER! WHEN USING THE DEVICE, PROTECT CHILDREN AND OTHER BYSTANDERS.

3.3. Personal safety

- a) Do not use the device when tired, ill or under the influence of alcohol, narcotics or medication which can significantly impair the ability to operate the device.
- b) The device is not designed to be handled by persons (including children) with limited mental and sensory functions or persons lacking relevant experience and/or knowledge unless they are supervised by a person responsible for their safety or they have received instructions on how to operate the device.

3.4. Safe device use

- a) Do not overload the device. Use appropriate tools for the given task. A correctly-selected device will perform the task for which it was designed better and in a safer manner. Platform scales have a zone of overload. It must be observed to avoid damages of the sensor and to prevent a loss of warranty (the balance gives a warning signal up to 135% overload).
- b) Do not use the device if the ON/OFF switch does not function properly (does not switch the device on and off). Devices which cannot be switched on and off using the ON/OFF switch are hazardous, should not be operated and must be repaired.
- c) When not in use, store in a safe place, away from children and people not familiar with the device, who have not read the user manual. The device may pose a hazard in the hands of inexperienced users.
- d) Keep the device out of the reach of children.

- e) Device repair or maintenance should be carried out by qualified persons, only using original spare parts. This will ensure safe use.
- f) To ensure the operational integrity of the device, do not remove factory-fitted guards and do not loosen any screws.
- g) Clean the device regularly to avoid a durable encrustation of dirt.
- h) Ensure that the weight placed on the scales is stable.
- i) Check that all screws, rods and brackets are secure.
- j) Do not weigh in case of extreme humidity and rain.
- k) Avoid working at high ambient temperatures.
- l) Allow the balance to recharge every 2-3 months for at least 24 hours to ensure maximum capacity and longevity.

ATTENTION! DESPITE THE SAFE DESIGN OF THE DEVICE AND ITS PROTECTIVE FEATURES, AND DESPITE THE USE OF ADDITIONAL ELEMENTS PROTECTING THE OPERATOR, THERE IS STILL A SLIGHT RISK OF ACCIDENT OR INJURY WHEN USING THE DEVICE. STAY ALERT AND USE COMMON SENSE WHEN USING THE DEVICE.

4. Use guidelines

The device is designed for performing control measurements and readings of the weight and/or the price of items/animals placed on the weighing platform.


This scale is not suitable for commercial transactions, legal purposes, medical uses, pharmaceutical applications, price determination in sales.

The user is liable for any damage resulting from unintended use of the device.

4.1. Device description

4.1.1. Control panel

For model SBS-BW-1T, SBS-BW-3T, SBS-BW-5T, SBS-BW-3T/1KG, SBS-BW-5T/2KG, SBS-BW-1T/0.5A, SBS-BW-3T/1A, SBS-BW-5T/2A

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4.2. Device use

4.2.1. Key functions

Warning signals and state Battery Level (Ib)

If your balance shows **LB** and produces a beep at the same time or the display gradually goes out, the battery is discharged. Connect the scale on the sector. The **LB** display and the beep tone stop, and the AC LED will turn on. To fully charge the battery, let the scale load for 4h.

4.2.2. X display security

For model SBS-PS-100, SBS-PS-200

Status "NO"

When the scale displays "NO", the circuit does not work. The balance must be repaired.

4.2.3. Buttons and options

With the switch I/O on the back of the device, you can turn the scale on. A self-test starts. It begins with a test of the screen and the diodes. The screen should display "0.0.0.0.0." and all the LEDs should light. Then the information on the current calibration values is displayed (from the factory: "0150.05": "0150" stands for 150kg and "0.5" for the measurement interval), as well as the programming of the calculator ($U_{er_3, 2}$ for ver3,2) and the state of the battery (e.g.,

dC_U5, 0 for 5.0 V) After the test, the scale is ready to use when it displays "0.00". Press the button I / O again to turn off the scale.

For the SBS-BW-3T/1KG scale model, the connection between the scale head and the scale platform is wireless. The communication between the head and the platform is being established right after switching on the scale head and the scale platform. The scale platform is equipped with three control lights:

- AC – means that the storage battery of the platform is charging (AC power is connected)
- Low Voltage – means that the platform storage battery needs to be charged
- SIGNAL – after proper wireless connection between the scale head and the scale platform this light will be blinking.

If the scale display will show Err24, it means that there is wireless connection error (the platform is off or the platform storage battery needs to be charged).

4.2.4. Error messages

For model SBS-BW-3T, SBS-BW-5T, SBS-BW-1T, SBS-BW-3T/1KG, SBS-BW-5T/2KG

- "Ad_Err" – The weight cannot be measured by the sensor. Necessary repairs.
- "Lo_bat" – Charge the battery
- "OL Overload" – Too much weight on the scale
- "Seterra" – The internal data are not correct. Necessary repairs.

4.2.5. Tare (zero)

Before weighing, make sure that ZERO and STABLE indicators are on the left and the desired unit is on the right. This means the optimum state of operation of the device. If the ZERO indicator is not turned on or the display is not showing zero after taking off the load from the scale, then the key **[ZERO/TARE]** must be pressed. The scale will confirm the tare showing "0.00". The weight window will display the net weight value.

Do not press the button **[ZERO/TARE]** until the weight is stable, (i.e., until the LED STABLE is lit) to avoid measurement errors.

To tare, place an empty container on the platform when the STABLE LED is lit. Then press **[ZERO/TARE]** for the scale to ignore the weight of the container.

CAUTION! After performing tare, the maximum weighing capacity is reduced by the weight of the container. For example:

- Container weight (b) = 90 kg
- Normal maximum capacity (s) = 150 kg
- Max. capacity left = $s - b = 60$ kg

If the scale is overloaded, immediately remove the weight. When overload occurs, the display will show: "OL ---- ()"

4.2.6. Hold (hold)

If the weight is superior to 0.05 kg in weighing mode, the save mode can be activated by pressing the **[HOLD]** button. At the same time, the HOLD LED lights up as a confirmation. Press the **[HOLD]** key again to exit the hold mode and to return to the weighing mode.

4.2.7. Unit change (unit/kg/lb, lb/kg)

With the **[UNIT]** key, you can change units between lb and kg. Press the button **[UNIT]** to choose the right unit.

4.2.8. Pack function

There are 2 selections under PACK function: IN & OUT. The default setting is “NO”, which means ALARM and PACK functions in non-work status; When it is under IN, it will alert between UPPER and LOWER limit; When it is under OUT setting, it will alarm that the weight exceeds the UPPER limit or below LOWER limit. It shows HHHHHH when exceeding UPPER limit and LLLLLL when below LOWER limit.

4.2.9. How to set in

Put the preset weight (e.g., 1000.0g), keep pressing **[PACK]** and the indicator shows SET-AL; then the scale enters into choosing function mode. Press +/- for enter into mode IN. Press **[ZERO/TARE]** to confirm and then enter into next setting; The indicator shows AL-UP, then enter into LOWER limit set (e.g., input 995.0g). Press **[ZERO/TARE]** to confirm and go for next step.

The indicator shows AL-DN, then enter into UPPER limit (e.g. input 1005.0g). Press **[ZERO/TARE]** for confirmation and the indicator shows IN. It means the setting is finished when the indicator keeps alarming.

Note: IN range between 1000.0g±5.0g is an example for customer’s references.

4.2.10. How to set upper/lower limit

[UNIT] and **[HOLD]** could be used to set UPPER & LOWER limit by inputting the numbers 0-9 into the setting. Press **[UNIT]** to increase and **[HOLD]** to decrease. Press **[ZERO/TARE]** to confirm and quit setting.

4.2.11. Setting of alarm

When indicator shows IN, it will alert the weight is between UPPER and LOWER limit. When window shows OUT, it will alert the load exceeds the UPPER limit or below LOWER limit. When window shows NO, it means it will not alarm without any load. Press **[ZERO/TARE]** to confirm the current selection.

4.2.12. How to set out

Put the preset weight (e.g., 1000.0g), keep pressing **[PACK]** and the indicator shows SET-AL; then the scale enter into choosing function mode. Press +/- for enter into mode OUT. Press **[ZERO/TARE]** to confirm and then enter into next setting; The indicator shows AL-UP, then enter into UPPER limit set(e.g., input 1005.0g). Press **[ZERO/TARE]** to confirm and go for next step. The display will show AL-DN and the LOWER limit should be set (e.g., input 995.0 kg). Press **[ZERO/TARE]** to confirm and then the display will show “1000.0 kg”. This means that settings are finished. Exceeding the LOWER limit (995.0 kg) is indicated by a beep and displaying LLLLLL on the display. Exceeding the UPPER limit is indicated by a beep and displaying HHHHHH on the display.

4.2.13. Parts counting

Setting the average weight of the sample.

Turn on the scale and wait for the weighing mode. Put the samples on the platform - the number of the samples must be 10/20/50/100/200/500/1000. Hold the **[PCS]** for a while until the display will show the "P-XX", then press **[UNIT]** or **[HOLD]** in order to set the number of samples – the number must be the same as the amount of pieces on the platform. Press **[PCS]** and the display will show the number of the pieces which are currently on the platform.

4.2.14. Counting

Put the parts on the platform and the scale will show the actual amount of these parts.

4.2.15. Calibration

For model SBS-BW-1T, SBS-BW-3T, SBS-BW-5T, SBS-BW-3T/1KG, SBS-BW-5T/2KG, SBS-BW-1T/0.5A, SBS-BW-3T/1A, SBS-BW-5T/2A SBS-PS-100, SBS-PS-200

Calibration:

- 1) Press and hold **[ZERO/TARE]** while turning on the scale. The display will show "-----".
- 2) Press **[PCS] + [UNIT] + [PCS] + [ZERO/TARE]** to enter the calibration mode.
- 3) Input the calibration weight that you will put on the platform:
 - Press **[PACK]** and **[PCS]** to move the cursor to your desired position.
 - Press **[HOLD]** and **[UNIT]** to increase or decrease the value.
 - Example: If you set display to 030.000, you need to put on a 30 kg weight.
- 4) Press **[ZERO/TARE]** to confirm.
- 5) Place the weight on the platform, then press **[ZERO/TARE]** again.
- 6) If the displayed value matches the weight, calibration is complete.

Parameter setting:

- 1) Press and hold **[ZERO/TARE]** while turning on the scale. The display will show "-----".
- 2) Press **[UNIT] + [PCS] + [UNIT] + [ZERO/TARE]**. The display will show "FUN-1".
 - **FUN-1: Maximum Capacity**
 - Press **[ZERO/TARE]** to enter the setting mode.
 - Press **[PACK]** and **[PCS]** to move the cursor to your desired position.
 - Press **[HOLD]** and **[UNIT]** to adjust the value.
 - Press **[ZERO/TARE]** to confirm. The display will then show "FUN-2".
 - **FUN-2: Zero Tracking Range**
 - Press **[ZERO/TARE]** to display the range of zero tracking.
 - Press **[ZERO/TARE]** again to confirm the range, and the display will then show


“FUN-3”.

- **FUN-3: Division**
 - Press [ZERO/TARE] to display the current division.
 - Press [HOLD] and [UNIT] to set the required division.
 - Press [ZERO/TARE] to confirm. The display will then show “FUN-4”.
 - **FUN-4: Decimal Point Position**
 - Press [ZERO/TARE] to display available options of the digit of point (e.g., 0, 0.0, 0.00, 0.000).
 - Press [HOLD] and [UNIT] to change the decimal position.
 - Press [ZERO/TARE] to confirm. The display will then show “FUN-5”.
 - Press [ZERO/TARE] twice to enter **FUN-6**.
 - **FUN-6: Backlight Setting (Only applicable to scales equipped with an LCD with backlight.)**
 - Press [ZERO/TARE] again to adjust the backlight.
 - Press [HOLD] and [UNIT] to cycle through options: “ON / OFF / AUTO”.
 - Choose your desired option and press [ZERO/TARE] to confirm.
- 3) When finished, the display will show “AND”. Press [ZERO/TARE] to enter weighing mode. Now you can start weighing.


4.2.16. Wireless pairing procedure

For model SBS-BW-3T/1KG | SBS- PS-100 | SBS-PS-200

- A) Press the button I / O of the platform base (Transmitter), the indicator [SIGNAL] will flash.
- B) Press and hold the [PCS] button continuously, then turn on the ON/OFF switch on the control panel (Receiver). Now the display will show [idSEt].

- C) Press the button [pair switch] on the platform base (Transmitter), the display will show [PASS], which means the pairing is completed.

- D) After approximately 1 second, return to the settings interface. Now the base can turn on again.

4.3. Cleaning and maintenance

4.3.1. General instructions

- a) Use only non-corrosive cleaners to clean the surface.
- b) After cleaning the device, all parts should be dried completely before using it again.
- c) Store the unit in a dry, cool place, free from moisture and direct exposure to sunlight.
- d) Do not spray the device with a water jet or submerge it in water.
- e) Do not allow water to get inside the device through vents in the housing of the device.
- f) The device must be regularly inspected to check its technical efficiency and spot any damage.
- g) Use a soft cloth for cleaning.
- h) Do not leave the battery in the device if it will not be used for a longer period of time.
- i) Do not use sharp and/or metal objects for cleaning (e.g., a wire brush or a metal spatula) because they may damage the surface material of the appliance.

- j) Do not clean the device with an acidic substance, agents of medical purposes, thinners, fuel, oils or other chemical substances because it may damage the device.

4.3.2. Safe removal of batteries and rechargeable batteries

Batteries are installed in the devices.

Remove used batteries from the device using the same procedure by which you installed them.

Recycle batteries with the appropriate organisation or company.

4.3.3. Disposing of used devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

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