USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

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DE	Produktname	Magnetrührer mit Heizplatte
EN	Product name	Pellet Press
PL	Nazwa produktu	Prasa do peletu
CZ	Název výrobku	Lis na pelety
FR	Nom du produit	Presse à granulés
IT	Nome del prodotto	Pressa per pellet
ES	Nombre del producto	Prensa de pellets
HU	Termék neve	Pelletprés
DA	Produktnavn	Pillepresse
FI	Tuotteen nimi	Pelletinpuristin
NL	Productnaam	Pelletpers
NO	Produktnavn	Pellets-presse
SE	Produktnamn	Pelletspress
PT	Nome do produto	Prensa de pellets
SK	Názov produktu	Lis na pelety
BG	Име на продукта	Преса за пелети
EL	Όνομα προϊόντος	Πρέσα πέλλετ
HR	Naziv proizvoda	Preša za pelete
LT	Produkto pavadinimas	Granulių presas
RO	Numele produsului	Presă pentru peleți
SL	Ime izdelka	Stiskalnica za pelete
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-MR-2000
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

1. Technical data

Table 1: Technical data of the product

Parameter description	Parameter value
Product name	Ceramic Magnetic Stirrer
Model	Model SBS-MR-2000
Rated voltage [V~]/ Frequency [Hz]	230/50
Rated power [W]	600
Fuse	T5AL, 250V
Stirring speed [rpm]	100-2000
The maximum volume of liquid [ml]	5000
Maximum temperature [°C]	0 - 350
Work surface dimensions [mm]	190x190
Device dimensions [mm]	Device dimensions [mm] 345 with power cord plugx190x400
Protection rating IP	20
Protection class	Protection class
Weight [kg]	4.65

2. General description

The user manual is designed to assist in the safe and trouble-free use of the device. The product is designed and manufactured in accordance with strict technical guidelines, using state-of-the-art technologies and components. Additionally, it is produced in compliance with the most stringent quality standards.

DO NOT USE THE DEVICE UNLESS YOU HAVE THOROUGHLY READ AND UNDERSTOOD THIS USER MANUAL.

To increase the product life of the device and to ensure trouble-free operation, use it in accordance with this user manual and regularly perform maintenance tasks. The technical data and specifications

in this user manual are up to date. The manufacturer reserves the right to make changes associated with quality improvement. The device is designed to reduce noise emission risks to a minimum, taking into account technological progress and noise reduction opportunities.

2.1. Legend

Icon	Description
	The product satisfies the relevant safety standards.
	Read instructions before use.
	The product must be recycled.
	ATTENTION! Electric shock warning!
	Do not use the device to experiment with flammable, explosive, poisonous and highly caustic substance.
	ATTENTION! Hot surface, risk of burns!

PLEASE NOTE! DRAWINGS IN THIS MANUAL ARE FOR ILLUSTRATION PURPOSES ONLY AND IN SOME DETAILS MAY DIFFER FROM THE ACTUAL PRODUCT.

3. Usage safety

ATTENTION! READ ALL SAFETY WARNINGS AND ALL INSTRUCTIONS. FAILURE TO FOLLOW THE WARNINGS AND INSTRUCTIONS MAY RESULT IN ELECTRIC SHOCK, FIRE AND/OR SERIOUS INJURY OR EVEN DEATH.

The terms "device" or "product" are used in the warnings and instructions to refer to: CERAMIC MAGNETIC STIRRER.

3.1. Electrical safety

a) The plug must fit the socket. Do not modify the plug in any way. Using original plugs and matching sockets reduces the risk of electric shock.

b) Avoid touching earthed elements such as pipes, heaters, boilers and refrigerators. There is an increased risk of electric shock if the earthed device is exposed to rain, comes into direct contact with a wet surface or is operating in a damp environment. Water getting into the device increases the risk of damage to the device and of electric shock.

c) Do not touch the device with wet or damp hands.

- d) Use the cable only for its designated use. Never use it to carry the device or to pull the plug out of a socket. Keep the cable away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cables increase the risk of electric shock.
- e) If using the device in a damp environment cannot be avoided, a residual current device (RCD) should be applied. The use of an RCD reduces the risk of electric shock.
- f) Do not use the device if the power cord is damaged or shows obvious signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service centre.
- g) To avoid electric shock, do not immerse the cord, plug or device in water or other liquids. Do not use the device on wet surfaces.

3.2. Safety in the workplace

- a) Make sure the workplace is clean and well lit. A messy or poorly lit workplace may lead to accidents. Try to think ahead, observe what is going on and use common sense when working with the device.
- b) Do not use the device in a potentially explosive environment, for example in the presence of flammable liquids, gases or dust. The device generates sparks which may ignite dust or fumes.
- c) If you discover damage or irregular operation, immediately switch the device off and report it to a supervisor without delay.
- d) If there are any doubts as to the correct operation of the device, contact the manufacturer's support service.
- e) Only the manufacturer's service point may repair the device. Do not attempt any repairs independently!
- f) In case of fire, use a powder or carbon dioxide (CO₂) fire extinguisher (one intended for use on live electrical devices) to put it out.
- g) Please keep this manual available for future reference. If this device is passed on to a third party, the manual must be passed on with it.
- h) Keep packaging elements and small assembly parts in a place not available to children.

REMEMBER! WHEN USING THE DEVICE, PROTECT CHILDREN AND OTHER BYSTANDERS.

3.3. Personal safety

- a) Do not use the device when tired, ill or under the influence of alcohol, narcotics or medication which can significantly impair the ability to operate the device.
- b) The device is not designed to be handled by persons (including children) with limited mental and sensory functions or persons lacking relevant experience and/or knowledge unless they are supervised by a person responsible for their safety or they have received instruction on how to operate the device.
- c) When working with the device, wear personal protective equipment that prevents:

splashing or contact with evaporating liquid,
contact with toxic or flammable gases.

d) To prevent accidental starting, make sure the knobs are in off position before connecting to the power source.

e) The device is not a toy. Children must be supervised to ensure that they do not play with the device.

3.4. Safe device use

a) Do not use the device if the knobs or control buttons do not function properly. Devices that cannot be controlled by the controls are dangerous, may not be operated and must be repaired.

b) Disconnect the device from the power supply before commencement of adjustment, cleaning and maintenance. Such a preventive measure reduces the risk of accidental activation.

c) When not in use, store in a safe place, away from children and people not familiar with the device who have not read the user manual. The device may pose a hazard in the hands of inexperienced users.

d) Keep the device in perfect technical condition. Before each use check for general damage and especially check for cracked parts or elements and for any other conditions which may impact the safe operation of the device. If damage is discovered, hand over the device for repair before use.

e) Keep the device out of the reach of children.

f) Device repair or maintenance should be carried out by qualified persons, only using original spare parts. This will ensure safe use.

g) To ensure the operational integrity of the device, do not remove factory-fitted guards and do not loosen any screws.

h) Do not move, adjust or rotate the device in the course of work.

i) Do not leave this appliance unattended while it is in use.

j) Clean the device regularly to prevent stubborn grime from accumulating.

k) The device is not a toy. Cleaning and maintenance may not be carried out by children without supervision by an adult person.

l) It is forbidden to interfere with the structure of the device in order to change its parameters or construction.

m) Keep the device away from sources of fire and heat.

n) NOTE! When the heating function is being used, the workplate warms up to a high temperature. Do not touch the hot parts of the device as they may cause burns.

o) Heated liquid can cause burns. It is advisable to use gloves when handling a vessel filled with hot liquid. Take special care not to spill the liquid.

- p) When working with flammable materials, pay special attention that the set heating temperature is not too high.
- q) Make sure that the liquid container can withstand the set heating temperature. A melted container may cause a fire.
- r) Never heat an empty container.
- s) Control the amount of liquid in containers. The liquid may evaporate.
- t) Spilling the flammable contents of the container onto the heating elements of the device may cause a fire.
- u) The set heating temperature must be at least 25°C lower than the ignition temperature of the stirred substance.
- v) Special attention must be paid to the hazards arising from:
- Flammable materials and substances with low boiling point
 - Overfilled containers
 - Containers made of brittle material or of low melting point
- w) Pathogenic materials must be stirred only in closed containers.
- x) Do not cover the device while it is in operation.
- y) Make sure that the power cord does not touch any heated parts of the device.
- z) Keep the device away from strong magnetic fields.
- aa) Do not use metal vessels or vessels with a thick bottom. This will make work with the device more difficult.
- bb) When mixing viscous liquids reduce the rotational speed and mix a small volume of liquid only.
- cc) The device can work continuously for 8 hours at moderate speed and 4 hours at high mixing speed.



ATTENTION! DESPITE THE SAFE DESIGN OF THE DEVICE AND ITS PROTECTIVE FEATURES, AND DESPITE THE USE OF ADDITIONAL ELEMENTS PROTECTING THE OPERATOR, THERE IS STILL A SLIGHT RISK OF ACCIDENT OR INJURY WHEN USING THE DEVICE. STAY ALERT AND USE COMMON SENSE WHEN USING THE DEVICE.

4. Use guidelines

The device is designed for use in laboratories, schools, universities, and similar institutions. It is used for non-contact stirring and for heating liquids in laboratory vessels.

The user is liable for any damage resulting from unintended use of the device.

4.1. Device description


Workplate

2. Control panel

3. Sensor pole

4. Sensor cable

5. Sensor
6. The part for adjusting the height and extension of the clamp
7. Sensor cable
8. Sensor cable socket
9. Sensor pole holder
10. Power cord socket
11. ON / OFF button

- A. TEMP – the button activates the heating temperature setting mode.
- B. SPEED – the button activates the mixing speed setting mode.
- C. TIME – the button activates the device operation time setting mode.
- D. SHIFT – the button increases the value in the parameter setting mode / Press and hold for 5 seconds to switch between readings from temperature sensors.
- E. AT – the button reduces the value in the parameter setting mode. Press and hold for 6 seconds to activate the auto-tuning function.
- F. STIRRING - mixing switch.
- G. Display.

4.2. Preparing for use

4.2.1. Appliance location

Ambient conditions:

The device should be used in a working environment with the following parameters:

- Ambient temperature : 0-40°C
- Ambient humidity : < 85%RH

Ensure good ventilation in the room in which the device is being used. There should be at least 10 cm distance between each side of the device and the wall or other objects. The device should always be used when positioned on an even, stable, clean, fireproof and dry surface, and be out of the reach of children and persons with limited mental and sensory functions. Position the device such that you always have access to the power plug. The power cord connected to the appliance must be properly grounded and correspond to the technical details on the product label.

4.3. Assembling the device

Before starting to use the device, install the temperature probe correctly.

- a) Place the pole in the holder. The hole overhang should enter the groove in the bottom part of the post.
- b) Fit the adjusting element and tighten the clamping knob.
- c) Place the probe clamp into the adjusting element and tighten the clamping knob.
- d) Place the probe in the clamp and make sure it is correctly clamped. The probe should be oriented with the cable upwards.
- e) Connect the probe cable to the suitable socket in the housing of the device.
- f) The device is ready to use.

4.4. Device use

- a) Pour the liquid into the container and insert the magnetic stirrer into the liquid.
- b) Place the vessel with liquid on the workplate. Dip the tip of the temperature sensor into the liquid to a depth of at least 20 mm. The recommended depth is $\frac{2}{3}$ of the height of the liquid. When using a metal container, do not place temperature probe on the bottom of the container. Placing the temperature probe on the bottom of the container may result in an incorrect measurement reading. The end of the measuring probe must be at least 5 mm from the bottom of the container; a distance of 10 mm is the best.
- c) Plug the power cord into the device and into the electrical outlet.
- d) Activate the heating temperature setting mode using the TEMP [A] button. Increase or decrease the temperature with [D] or [E] buttons. Confirm the selection and return to the start screen with the TEMP [A] button.
- e) Switch on the speed setting mode using the SPEED button [B]. Increase or decrease the speed with [D] or [E] buttons. Confirm the selection and return to the start screen with the TEMP [A] button. To start mixing, press the [F] button.
- f) Activate the device operation time setting mode with the TIME [C] button. Increase or decrease the operation time with the [D] or [E] buttons. Confirm the selection and return to the start screen with the TIME [A] button.
 - If the setpoint is "0", the heating and the set temperature will be maintained continuously until the time is changed or the device is turned off. If the mixing function is

activated, it will be active until the time is changed, mixing is turned off by pressing [F] or the device is turned off.

- If the set value is ≥ 1 , the heating will continue for the set time. If the mixing is switched on during heating, it will be active until the set time has elapsed, the mixing is switched off with the [F] key or the device is switched off.
- After the set operation time has elapsed, the message "End" will be displayed (the heating and mixing will turn off) and the device will beep. To turn off the beep, press any key. **IMPORTANT!** Do not turn on the heat under an empty vessel.

g) After work with the device, switch off the mixing, set the heating temperature to 0, switch off the device with the ON / OFF button and disconnect it from the power supply.

TEMPERATURE MEASUREMENT POINT

The device enables switching between readings from temperature sensors after pressing and holding the "SHIFT" button [D] for 5 seconds.

"A" PT100 sensor for measuring the temperature of liquid

"B" K-type sensor for measuring the temperature of the workplate.

OVER TEMPERATURE ALARM

If the set temperature is exceeded, an over temperature alarm will come on. A beep is heard and the "ALM" alarm message lights up. To turn off the beep, press any key. If the over temperature alarm occurs as a result of changing the setpoint to a lower value than the displayed measurement, the alarm message will light up but there will be no beep. When the over temperature alarm is active, the controller automatically disconnects the heating power.

TEMPERATURE MEASUREMENT ERROR ALARM

"E-X" on the device's display means that the temperature sensor has failed, the temperature exceeds the measuring range or the controller itself is faulty. The controller will automatically disconnect the heating power, a beep will be heard and the "ALM" alarm message will light up. In this case, the technical condition of the temperature sensor and cabling should be carefully checked.

E-1: indicates the temperature error of the first channel (resistance thermometer PT100);

E-2: indicates the temperature error of the second channel (K-type thermocouple);

E-3: indicates the ambient temperature error.

To turn off the beep, press any key.

AUTOTUNING

The "Autotuning" function enables automatic tuning of PID settings to the most optimal ones for the device operation.

After pressing and holding the "AT" [E] button for 6 seconds, a window will open asking you to confirm the start of the "autotuning" mode. To start the mode, change the value from 0 to 1 (using the D and E buttons) and confirm the selection with the "TEMP" [A] key. The device will start working in the autotuning mode. When the "AT" icon disappears from the display, the

autotuning operation is finished. To turn off the autotuning mode manually, press and hold the "AT" [E] key for 6 seconds. The autotuning function will not activate when the time display is "End".

During the autotuning process:

If an over temperature alarm occurs, the beep will not be heard, the "ALM!" warning message will not light up but the heating will be turned off.

The "TEMP" button is inactive.

SETTING THE CONTROLLER PARAMETERS

To manually set the controller parameters, press the "TEMP" [A] key for 5 seconds. "Lc" will appear on the display and a password will be required. Enter the password using the [D] and [E] keys, confirm with the "TEMP" [A] key. Depending on what password is given and if it is correct, the device will go to the appropriate mode for selecting the controller parameter settings.

List of parameters I

Message	Function name	Function description	Factory value a (setting range)
Lc	Password	Lc = 3	0
P	Proportionality range	Adjustment of the proportional function.	(0.1~300.0°C) 30.0
I	Integral time	Adjustment of the integral function.	(1~1000s) 200
d	Differentiation time	Adjustment of the differential function.	(0~1000s) 150
T	Time to reach the set temperature	The time at which the set temperature is to be reached.	(1~60s) 5
doT1	Accuracy of the first sensor	0: SP without a decimal place, 1: SP with a decimal place.	(0~1) 0
AH1	Over temperature alarm in the first sensor	When selecting the first sensor, if "SV> (SP + AH1)" the "ALM!" message appears. A beep is heard and the heating switches off.	(0~100.0°C) 20.0
Pb1	Zero point adjustment in the first sensor	When choosing the first sensor for correction of errors generated during low temperature measurement. Pb1 = actual value - measured value	(-50.0 ~50.0°C) 0.0
PL1	Maximum point adjustment in the first sensor	When choosing the first sensor for correction of errors generated during high temperature measurement. PL1 = 1000 × (actual value - measured value) / measured value	(-999~999) 0
SPH	Maximum setpoint	Maximum temperature setting.	(0~400.0°C) 400.0

List of parameters II

Message	Function name	Function description	Factory value a (setting range)
Lc	Password	Lc = 9	0
FoP	Temperature at which the fan turns on.	If ambient temperature > FoP", the fan turns on.	(0~80°C) 40
FcP	Temperature at which the fan turns off.	If "Ambient temperature <FcP", the fan turns off.	(0~80°C) 30
ndT	Timer mode	0: Display of the operating time when the measured temperature reaches the set value. 1: The working time is always displayed.	(0~1) 0
Hn	Timer unit	0: minute 1: hour	(0~1) 0
SPd	Constant temperature deviation	When SP> = (SV - SPd), the controller enters the temperature maintaining state.	(0.1 ~100.0°C) 0.5
EST	The beep sound time after the timer expires.	When the set time elapses, the device will beep for EST seconds. Note: if EST = 9999, it means that the beep will be sounded continuously.	(0~9999s) 60
PoT	Time after which the mixing switch will start.	The mixing switch will start after PoT seconds when it is pressed and held.	(0~10s) 2
AH2	Over temperature alarm in the second sensor.	When the second sensor is selected, if "SV> (SP + AH1)", the "ALM!" message will be displayed. A beep is heard and the heating output switches off.	(0~100°C) 20
Pb2	Zero point adjustment in the second sensor.	When the second sensor is selected for the correction of errors generated when measuring the low temperature. Pb2 = actual value - measured value	(-50~50°C) 0
PL2	Adjusting the maximum point in the second sensor.	When selecting the second sensor for the correction of errors generated when measuring the high temperature. PL2 = 1000 × (actual value - measured value) / measured value	(-999 ~999) 0

List of parameters III

Message	Function name	Function description	Factory value a (setting range)
Lc	Password	Lc = 27	0
Fc	Temperature unit	0: degree Celsius; 1: degree Fahrenheit.	(0~1) 0

List of parameters IV

Message	Function name	Function description	(setting range) Factory value a
Lc	Password	Lc = 67	0
rST	Restoring default values	0: cancel to reset to the default value; 1: confirm to restore the default value.	(0~1) 0

To manually set the controller parameters, press the "SPEED" key for 5 seconds. "Lc" will appear on the display and a password will be required. Enter the password using the [D] and [E] keys, confirm with the "SPEED" [B] key. After confirming the correct password, the device will go to the mode appropriate for selecting the controller parameter settings.

List of parameters V

Message	Function name	Function description	(setting range) Factory value a
Lc	Password	Lc = 3	0
Pd	Proportionality range	Adjustment of the proportional function.	(1~99) 40
Id	Integral time	Adjustment of the integral function.	(1~99) 2
InT	Mixing speed increase time	The time required to increase the speed from minimum to maximum.	(5~60) 10
dET	Mixing speed reduction time.	The time required to reduce the speed from maximum to minimum.	(5~60) 10
SdL	Minimum setpoint	Minimum speed setting.	(60~SdH) 200
SdH	Maximum setpoint	Maximum speed setting.	(SdL~6000) 2000
PoL	Pairs of poles	Motor pole pairs	(1~32) 1
db	Speed measurement	Permissible value of speed measurement	(0~99) 5

4.4.1. Cleaning and maintenance

- a) Unplug the mains plug and allow the device to cool completely before each cleaning, adjustment or replacement of accessories, or if the device is not being used.
- b) Use only non-corrosive cleaners to clean the surface.
- c) After cleaning the device, all parts should be dried completely before using it again.
- d) Store the unit in a dry, cool place, free from moisture and direct exposure to sunlight.
- e) Do not spray the device with a water jet or submerge it in water.
- f) Do not allow water to get inside the device through vents in the housing of the device.

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