USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

ELECTRIC CENTRIFUGE

DE	Produktname	Laborzentrifuge
EN	Product name	Electric centrifuge
PL	Nazwa produktu	Wirówka laboratoryjna
CZ	Název výrobku	Laboratorní odstředivka
FR	Nom du produit	Centrifugeuse électrique
IT	Nome del prodotto	Centrifuga elettrica
ES	Nombre del producto	Centrífuga eléctrica
HU	Termék neve	Laboratóriumi centrifugál
DA	Produktnavn	Elektrisk centrifuge
FI	Tuotteen nimi	Sähkökäyttöinen laboratoriokeskittäjä
NL	Productnaam	Elektrische laboratoriumcentrifuge
NO	Produktnavn	Elektrisk sentrifuge
SE	Produktnamn	Elektrisk centrifug
PT	Nome do produto	Centrífuga eléctrica
SK	Názov produktu	Laboratórna centrifuga
BG	Име на продукта	Лабораторна центрофуга
EL	Όνομα προϊόντος	Ηλεκτρική φυγόκεντρος
HR	Naziv proizvoda	Laboratorijska centrifuga
LT	Produktu pavadinimas	Laboratorinė centrifuga
RO	Numele produsului	Centrifugă electrică
SL	Ime izdelka	Laboratorijska centrifuga
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LZ-1000SLS
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value
Product name	Electric centrifuge
Model	SBS-LZ-1000SLS
Rated voltage [V \sim] / Frequency [Hz]	230 / 50
Rated power [W]	100
Rotation speed [rpm]	4000
Speed accuracy [rpm]	± 20
Maximum capacity [ml]	8x15
Maximum Relative Centrifugal Force (RCF) [g]	1880
Time range [min]	0-99
Protection class	I
Protection rating IP	IPX0
Dimensions [Width x Depth x Height; mm]	310 x 245 x 205
Weight [kg]	6.4

1. General description

The user manual is designed to assist in the safe and trouble-free use of the device. The product is designed and manufactured in accordance with strict technical guidelines, using state-of-the-art technologies and components. Additionally, it is produced in compliance with the most stringent quality standards.

**DO NOT USE THE DEVICE UNLESS YOU HAVE
THOROUGHLY READ AND UNDERSTOOD THIS USER
MANUAL.**

To increase the product life of the device and to ensure trouble-free operation, use it in accordance with this user manual and regularly perform maintenance tasks. The technical data and specifications in this user manual are up to date. The manufacturer reserves the right to make changes associated with quality improvement. The device is designed to reduce noise emission risks to a minimum, taking into account technological progress and noise reduction opportunities.

Legend

	The product satisfies the relevant safety standards.
	Read instructions before use.
	The product must be recycled.
	WARNING! or CAUTION! or REMEMBER! Applicable to the given situation. (general warning sign)
	Wear protective goggles.
	Wear protective gloves.
	Wear a protective apron.
	ATTENTION! Electric shock warning!
	ATTENTION! Rotating parts, entanglement hazard!

	ATTENTION! Fire hazard - flammable materials!
	WARNING! Toxic substances, danger of poisoning!
	Only use indoors.

PLEASE NOTE! Drawings in this manual are for illustration purposes only and in some details may differ from the actual product.

2. Usage safety

ATTENTION! Read all safety warnings and all instructions. Failure to follow the warnings and instructions may result in electric shock, fire and/or serious injury or even death.

The terms “device” or “product” are used in the warnings and instructions to refer to < *ELECTRIC CENTRIFUGE* >.

2.1. Electrical safety

- a) The plug must fit the socket. Do not modify the plug in any way. Using original plugs and matching sockets reduces the risk of electric shock.
- b) Avoid touching earthed elements such as pipes, heaters, boilers and refrigerators. There is an increased risk of electric shock if the earthed device is exposed to rain, comes into direct contact with a wet surface or is operating in a damp environment. Water getting into the device increases the risk of damage to the device and of electric shock.
- c) Do not touch the device with wet or damp hands.
- d) Use the cable only for its designated use. Never use it to carry the device or to pull the plug out of a socket. Keep the cable away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cables increase the risk of electric shock.

- e) If using the device in a damp environment cannot be avoided, a residual current device (RCD) should be applied. The use of an RCD reduces the risk of electric shock.
- f) Do not use the device if the power cord is damaged or shows obvious signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service centre.
- g) To avoid electric shock, do not immerse the cord, plug or device in water or other liquids. Do not use the device on wet surfaces.
- h) ATTENTION! DANGER TO LIFE! While cleaning, never immerse the device in water or other liquids.
- i) Do not use in very humid environments or in the direct vicinity of water tanks.
- j) Prevent the device from getting wet. Risk of electric shock!

2.2. Safety in the workplace

- a) Make sure the workplace is clean and well lit. A messy or poorly lit workplace may lead to accidents. Try to think ahead, observe what is going on and use common sense when working with the device.
- b) If you discover damage or irregular operation, immediately switch the device off and report it to a supervisor without delay.
- c) If there are any doubts as to the correct operation of the device, contact the manufacturer's support service.
- d) Only the manufacturer's service point may repair the device. Do not attempt any repairs independently!
- e) In case of fire, use a powder or carbon dioxide (CO₂) fire extinguisher (one intended for use on live electrical devices) to put it out.
- f) Children or unauthorised persons are forbidden to enter a work station. (A distraction may result in loss of control over the device).
- g) Regularly inspect the condition of the safety labels. If the labels are illegible, they must be replaced.
- h) Please keep this manual available for future reference. If this device is passed on to a third party, the manual must be passed on with it.

-
- i) Keep packaging elements and small assembly parts in a place not available to children.
 - j) Keep the device away from children and animals.

Remember! When using the device, protect children and other bystanders.

2.3. Personal safety

- a) Do not use the device when tired, ill or under the influence of alcohol, narcotics or medication which can significantly impair the ability to operate the device.
- b) The device is not designed to be handled by persons (including children) with limited mental and sensory functions or persons lacking relevant experience and/or knowledge unless they are supervised by a person responsible for their safety or they have received instruction on how to operate the device.
- c) The device can be handled only by physically fit persons who are capable of handling it, properly trained, familiar with this manual and trained within the scope of occupational health and safety.
- d) When working with the device, use common sense and stay alert. Temporary loss of concentration while using the device may lead to serious injuries.
- e) Use personal protective equipment as required for working with the device, specified in section 1 (Legend). The use of correct and approved personal protective equipment reduces the risk of injury.
- f) To prevent the device from accidentally switching on, make sure the switch is on the OFF position before connecting to a power source.
- g) Do not wear loose clothing or jewellery. Keep hair, clothes and gloves away from moving parts. Loose clothing, jewellery or long hair may get caught in moving parts.
- h) The device is not a toy. Children must be supervised to ensure that they do not play with the device.
- i) Do not put your hands or other items inside the device while it is in use!

2.4. Safe device use

- a) Do not use the device if the ON/OFF switch does not function properly (does not switch the device on and off). Devices which cannot be switched on and off using the ON/OFF switch are hazardous, should not be operated and must be repaired.
- b) Disconnect the device from the power supply before commencement of adjustment, cleaning and maintenance. Such a preventive measure reduces the risk of accidental activation.
- c) When not in use, store in a safe place, away from children and people not familiar with the device who have not read the user manual. The device may pose a hazard in the hands of inexperienced users.
- d) Keep the device in perfect technical condition. Before each use check for general damage and especially check for cracked parts or elements and for any other conditions which may impact the safe operation of the device. If damage is discovered, hand over the device for repair before use.
- e) Keep the device out of the reach of children.
- f) Device repair or maintenance should be carried out by qualified persons, only using original spare parts. This will ensure safe use.
- g) To ensure the operational integrity of the device, do not remove factory-fitted guards and do not loosen any screws.
- h) Avoid situations where the device stops working during use due to excessive loading. This may result in overheating of the drive elements and damage to the device.
- i) Do not touch articulated parts or accessories unless the device has been disconnected from the power source.
- j) Do not move, adjust or rotate the device in the course of work.
- k) Do not leave this appliance unattended while it is in use.
- l) Clean the device regularly to prevent stubborn grime from accumulating.
- m) The device is not a toy. Cleaning and maintenance may not be carried out by children without supervision by an adult person.
- n) It is forbidden to interfere with the structure of the device in order to change its parameters or construction.
- o) Keep the device away from sources of fire and heat.

- p) Do not overload the device.
- q) Do not cover the ventilation openings!
- r) A flammable, explosive or corrosive test sample may damage or shorten the life of the device. Any damage caused by the above-mentioned type test samples are not covered by the warranty.
- s) The maximum rotor speed must not be exceeded.
- t) The condition of the rotor should be frequently checked (for cracks, corrosion, damaged seals, scratches, etc.). Do not use the device if the rotor is not in good technical condition.
- u) You must never open the device cover while it is working! The lid can be opened with the device switched off and the rotor at a standstill.
- v) Uneven alignment or use of damaged test tubes can cause severe vibrations and damage the device.

ATTENTION! Despite the safe design of the device and its protective features, and despite the use of additional elements protecting the operator, there is still a slight risk of accident or injury when using the device. Stay alert and use common sense when using the device.

3. Use guidelines

The device is used to divide various mixtures, suspensions and body fluids into components of different densities.

The product is outside of the scope of MDD 93/42/EEC Medical Devices regulations. It is not a medical device and should not be used for medical research and tests in hospitals and other health facilities.

The user is liable for any damage resulting from unintended use of the device.

3.1. Device description


3.1.1 Product overview

1. Lid
2. Rotor
3. Feet
4. Control panel

Please note: The ON/OFF button and the power socket are located on the back of the device housing (not shown in the figure above).

3.1.2 Control panel description

1. Display
2. START button
3. Program selection button
4. Displayed parameter selection button (timer, rotational speed, maximum relative centrifugal force (RCF))
 5. "Up" button
 6. "Down" button
7. Device operation LED
8. Device stop LED
9. STOP button

3.2. Preparing for use

3.2.1 Appliance location

The ambient temperature must not exceed 45°C and the relative humidity should not exceed 85%. Only use the device in well-ventilated areas. Do not obstruct the air outlet of the device. Keep the device away from any hot surfaces. The device should always be used when positioned on an even, stable, clean, fireproof and dry surface, and be out of the reach of children and persons with limited mental and sensory functions. The valve that disables the compressed air supply should be immediately accessible from the area where the device is being used. Make sure

that the pressure of the air supplying the device does not exceed the recommended range!

3.2.2 Assembling the device

Place the device on a flat, stable surface so that all four rubber feet are firmly set on the surface.

Connect the power cord to the device and its other end to a 230V AC socket. Please note: The power socket must be grounded. This will reduce the risk of electric shock.

3.3. Device use

1. Check that the nut on the rotor has not loosened, if so, tighten it.
2. Get the test tubes with a preparation ready for centrifugation and place them in the device so that they meet the following premises:
 - The test tube should be filled to a maximum of 75% of its total volume.
 - Place the test tubes symmetrically in the rotor.
 - The sample-filled tubes placed opposite each other must be of equal weight.
 - Always verify the condition of the test tube used. Do not insert cracked or deformed test tubes in the device!
3. Close the device lid.
4. Switch the device connected to the power supply on, using the ON/OFF button at the back of the device.
5. Set the device parameters.
 - **Speed setting**

Press and hold the (4) button on the control panel until “A1 ✕ 3100” flashes on the display. Then, use the “Up” (5)/“Down” (6) buttons to set the preferred speed value. The maximum speed must not be exceeded.

Note: If the density of the sample intended for centrifugation exceeds 1.2 g/ml, the maximum speed value should be recalculated acc. to the formula below:

$$N = n\sqrt{1.2 / S}$$

Where,

- N - new maximum rotation speed
- n- original maximum rotation speed
- S- density of the tested sample

- **Time setting**

Press and hold the (4) button on the control panel until “04:39” flashes on the display. Then, use the “Up” (5)/”Down” (6) buttons to set the preferred time value for the device operation.

- **Acceleration and deceleration values setting**

Press and hold the (4) button on the control panel until “2” flashes on the display. Then use the “Up” (5)/”Down” (6) buttons to set the preferred value.

Please note: The device is factory set to “b5” . There is no need for the user to set this parameter themselves.

6. To start the process, press the START button (2), and the power LED (7) lights up green. The display will show the current speed value (the set value will be reached in about 1-2 minutes), and the countdown time to the end of spinning.
7. To stop the device, press the STOP (9) button or wait until the timer counts down to the value 0. The LED (8) will turn red. The device will not be fully stopped until the actual speed value is 0 on the display.
8. Remove the test tubes with centrifuged samples.
9. Switch the device off using the On/Off button. Even after the device is turned off, it is still live. Disconnect the device from the power source when not in use.

3.4. Cleaning and maintenance

3.4.1 General instructions

- a) Always unplug the device before cleaning it.
- b) Use only non-corrosive cleaners to clean the surface.
- c) After cleaning the device, all parts should be dried completely before using it again.

-
- d) Store the unit in a dry, cool place, free from moisture and direct exposure to sunlight.
 - e) Do not spray the device with a water jet or submerge it in water.
 - f) Do not allow water to get inside the device through vents in the housing of the device.
 - g) Clean the vents with a brush and compressed air.
 - h) The device must be regularly inspected to check its technical efficiency and spot any damage.
 - i) Use a soft, damp cloth for cleaning.
 - j) Do not use sharp and/or metal objects for cleaning (e.g. a wire brush or a metal spatula) because they may damage the surface material of the appliance.
 - k) Do not clean the device with an acidic substance, agents of medical purposes, thinners, fuel, oils or other chemical substances because it may damage the device.
 - l) If there is the concern that the device is contaminated with toxic, pathogenic, radioactive substances or infectious test tubes with blood that are harmful to health promptly sterilize the equipment in accordance to valid laboratory methods.
 - m) Periodically lubricate the central part of the rotor.

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