USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

VACUUM DRYING OVEN

DE	Produktname	Vakuumtrockenschrank
EN	Product name	Vacuum drying oven
PL	Nazwa produktu	Suszarka próżniowa
CZ	Název výrobku	Vakuová sušicí pec
FR	Nom du produit	Étuve de séchage sous vide
IT	Nome del prodotto	Essiccatore a vuoto
ES	Nombre del producto	Horno de secado al vacío
HU	Termék neve	Vákuumszáritó kemence
DA	Produktnavn	Vakuumtørreovn
FI	Tuotteen nimi	Höyrytyskuivaustaso
NL	Productnaam	Vacuümdroogkast
NO	Produktnavn	Vakuumtørkeovn
SE	Produktnamn	Vakuumtorkugn
PT	Nome do produto	Estufa de secagem a vácuo
SK	Názov produktu	Vakuová sušiča pec
BG	Име на продукта	Вакуумна сушилна пещ
EL	Όνομα προϊόντος	Κλίβανος ξήρανσης κενού
HR	Naziv proizvoda	Vakuumsko sušilnica
LT	Produktu pavadinimas	Vakuuminė džiovinimo spinta
RO	Numele produsului	Etuvă de uscare în vid
SL	Ime izdelka	Vakuumsko sušilno omaro
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-ADO-2500 SBS-VDO-92
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical Data

Description of the parameter	Value of the parameter	
Product name	Vacuum drying oven	
Model	SBS-ADO-2500	SBS-VDO-92
Supply voltage [V~] / Frequency [Hz]	230/50	
Rated power [W].	850	2100
Safety class		
Inner dimensions (W x D x H) [mm]	300 x 300 x 180	450 x 450 x 350
Dimensions [mm]	480 x 480 x 605	640 x 590 x 780
Weight [kg]	45.3	83
Drying chamber volume [l]	25	92
Temperature setting range [°C]	0-250	
Temperature accuracy [°C]	+/- 1	
Number of shelves [pcs].	2	
Distance between shelves [mm]	100	185
Maximum shelf load [kg]	15	
Length of supply cable [m]	1.5	
Maximum chamber leakage value [kPa/h]	0.3	
Maximum vacuum value [Pa]	133	
Time to reach temperature [min.]	80	120

1. General Description

The instruction manual is intended to assist in safe and reliable use. The product is designed and manufactured strictly according to technical specifications using the latest technology and components and maintaining the highest quality standards.

CAREFULLY READ AND UNDERSTAND THIS MANUAL BEFORE PROCEEDING.

To ensure long and reliable operation of the unit, make sure to operate and maintain it properly in accordance with the guidelines in this instruction manual. The technical data and specifications contained in this instruction manual are up to date. The manufacturer reserves the right to make changes in order to improve the quality. Taking the technical progress and the possibility of reducing noise into account, the unit is designed and built in such a way so that risks resulting from noise emissions are reduced to the lowest possible level.

Explanation of symbols

	The product complies with applicable safety standards.
	Please read the instructions before use.
	Recyclable product.
	CAUTION! or WARNING! or REMINDER! describing a situation. (general warning sign).
	Wear protective gloves.
	CAUTION! Warning of electric shock!
	Caution! Hot surface can cause burns!
	For indoor use only.

CAUTION! Danger of fire!

CAUTION! The illustrations in this instruction manual are for reference only and may differ from the actual product in some details.

The original instruction manual is in the German language version. Other language versions are translations from German.

2. Safety of use

CAUTION! Read all safety warnings and all instructions. Failure to follow the warnings and instructions may result in electric shock, fire and/or severe personal injury or death.

The term "unit" or "product" in the warnings and in the description of the instructions refers to the VACUUM DRYING OVEN.

2.1. Electrical safety

- a) The plug of this unit must fit into the outlet. Do not modify the plug in any way. Original plugs and matching outlets reduce the risk of electric shock.
- b) Avoid touching grounded parts, such as pipes, heaters, ovens, and refrigerators. There is an increased risk of electric shock if your body is grounded and touches the unit exposed to direct rain, wet pavement, and operation in a humid environment. If water enters the unit, there is an increased risk of damage to the unit and electric shock.
- c) Do not touch the unit with wet or damp hands.
- d) Do not use the cord in an unintended manner. Never use it to carry the unit or to pull the plug out of the socket. Keep the cord away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cords increase the risk of electric shock.
- e) If you cannot avoid using the unit in a wet environment, use a residual current unit (RCD). Using an RCD reduces the risk of electric shock.

-
- f) Do not use the unit if the power cord is damaged or shows signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service department.
 - g) To avoid electric shock, do not immerse the cable, plug, or unit itself in water or other liquid. Do not use the unit on wet surfaces.

2.2. Safety in the workplace

- a) Keep the work area tidy and well lit.
- b) If you have any doubts as to whether the unit is working properly or if it is damaged, contact the manufacturer's service department.
- c) Only the manufacturer's service department can repair the unit. Do not carry out repairs yourself!
- d) In case of open flames or fire, use only dry powder or snow (CO₂) fire extinguishers to extinguish the live equipment.
- e) No children or unauthorized persons are allowed in the work area.
- f) Use the unit in a well-ventilated area.
- g) Connecting and disconnecting the vacuum pump should be done with the air valve closed.
- h) Check the condition of the safety stickers regularly. Replace them if they are illegible.
- i) Keep these instructions for use for future reference. If the unit is to be passed on to a third party, the operating instructions must also be handed over together with the unit.
- j) Keep the packaging and small assembly parts out of the reach of children.
- k) Keep the unit away from children and animals.
- l) When using this unit together with other units, also follow the other instructions for use.

Please note! Keep children and other bystanders safe while operating the unit.

2.3. Personal safety

- a) Do not operate this unit if you are tired, ill or under the influence of alcohol, drugs or medication that could impair your ability to operate the unit.
- b) The unit is not intended to be used by persons (including children) with reduced mental, sensory or intellectual functions or persons who lack experience and/or knowledge unless they are supervised or have been instructed by a person responsible for their safety on how to operate the unit.

-
- c) The unit may be operated by persons who are physically fit, capable of operating it and appropriately trained, and who have read this instruction manual and have been trained in occupational safety and health.
 - d) Use caution and common sense when operating this unit. A moment's inattention during operation may result in serious personal injury.
 - e) The unit is not a toy. Children should be watched to ensure that they do not play with the unit.

2.4. Safe use of the unit

- a) Disconnect the unit from the power supply and the vacuum pump before adjusting, cleaning, and maintenance.
- b) Keep unused equipment out of the reach of children and out of the reach of anyone unfamiliar with the unit or this instruction manual. These units is dangerous in the hands of inexperienced users.
- c) Keep the unit in good working condition.
- d) Keep the unit out of the reach of children.
- e) Repairs and maintenance should be carried out by qualified personnel using only original spare parts. This will ensure the safety of use.
- f) To ensure the designed operational integrity of the unit, do not remove factory-installed covers or loosen screws.
- g) When transporting or moving the unit from storage to the place of use, observe the health and safety rules for manual handling applicable in the country where the unit is used.
- h) Do not move, shift, or rotate the machine while in operation.
- i) Clean the unit regularly to prevent permanent dirt build-up.
- j) The unit is not a toy. Cleaning and maintenance must not be performed by children without adult supervision.
- k) Do not tamper with the unit to alter its performance or design.
- l) Keep the unit away from sources of fire and heat.
- m) Check the condition of the door and observation window regularly for any scratches or other damage. If any damage is found on the door or window, creating negative pressure in the unit is prohibited. This also applies to the other functional elements of the unit.
- n) Use of the unit in rooms with flammable or explosive gas mixtures is prohibited. The unit is not designed as explosion-proof. When switching the unit off / on, ignition / explosion may occur due to an electric arc generated in the main switch of the unit.
- o) Pay particular attention to the physical and chemical properties (e.g. flash point, etc.) of the contents of the oven. Paints, enamels or other similar materials, whose reagents may form an explosive mixture with air, must not be dried, evaporated or fired. There must be no situation in which flammable

gas or a flammable gas-air mixture is formed in the chamber of the vacuum oven or in its immediate vicinity.

- p) When starting the vacuum oven for the first time, ensure constant supervision of the operation of the unit until stable operating conditions are achieved.
- q) Large amounts of dust or corrosive vapors inside the drying chamber or around the unit can produce deposits inside the unit and lead to short circuits or damage to the electronics. Therefore, it is important to take all appropriate precautions against excessive accumulation of dust or corrosive vapors inside and outside the unit.
- r) In order to ensure proper air circulation in the drying chamber, sufficient clearance must be provided between the load and the walls inside the chamber. Do not place any load on the bottom, side walls, or ceiling of the chamber.
- s) If problems such as smoke, or a strange odor emanating from the unit occur, immediately turn off the unit and disconnect power completely.
- t) Do not touch the door or glass during or immediately after operation. High temperature can cause severe burns.


CAUTION! Although the product has been designed to be safe, with adequate safeguards, and despite the additional safety features provided to the user, there is still a slight risk of accident or injury when handling the unit. You are advised to use caution and common sense when using this product.

3. Rules of use

The unit is used for the thorough drying of materials with complex compositions and those that oxidize/decompose easily. It can also be used for sterilization.


The user is responsible for any damage resulting from misuse.

3.1. Description

1. Door
2. Drying chamber
3. Door latch
4. Control panel with display
5. Vacuum gauge
6. Air Inlet Valve ("AIR INLET; BALANCE HOLE")
7. Air outlet valve ("AIR OUTLET")

CONTROL PANEL WITH DISPLAY [4].

- 8. "PV" temperature value display.
- 9. "SV" time display
- 10. "RUN/AT" indicator light
- 11. "OUT" indicator light
- 12. "MK" indicator light
- 13. "ALM" indicator light
- 14. SET button
- 15. "SHIFT/AT" function button
- 16. Button decreasing the set value
- 17. Button increasing the set value

3.2. Preparation for operation

The ambient temperature must not exceed 40°C and ambient humidity should not exceed 50%. Use the unit in properly ventilated spaces. Keep the unit away from any hot surfaces. Always operate the unit on a level, stable, clean, fireproof and dry surface that can support the weight of the unit and out of the reach of children and persons with limited mental, sensory or intellectual functions. It is important that the unit is in an exact horizontal position. Do not place the unit near open fire sources, electric stoves, coal stoves or other heat sources, or in direct sunlight. Maintain a minimum clearance of 10 cm from any wall of the unit. The work area should provide immediate access to the air valves and the mains plug. Ensure that the power supply to the unit corresponds to that specified on the identification plate! Wash and disinfect the drying chamber and shelves before first use.

3.3. Working with the unit

1. Place the material to be dried on the shelf. Keeping at least 1/3 of the shelf area empty for optimal conditions is recommended.
2. Close the air inlet valve [6] and the door. Make sure they are tight!
3. Connect the vacuum pump hose to the air outlet valve [7].
4. Unscrew the valve [7].
5. Switch on the vacuum pump.
6. Observe the vacuum gauge. A pointer going counterclockwise indicates venting of the chamber.
7. When the vacuum gauge reads below zero (vacuum has been built up in the chamber), turn off the valve [7] and switch off and disconnect the vacuum pump. Pay special attention to the sequence of actions leading to the pump disconnection: a change may cause the pump oil to leak into the drying chamber.
8. Connect the oven to the power source. The control panel will be illuminated.
9. Set the desired temperature:
 - (a) Press the setting button [14].
 - b) Use the arrows [16] and [17] to set the desired temperature value.
 - c) Confirm by pressing the setting button [14] again.
10. Set the desired operating time:
 - (a) Press the setting button [14] twice.
 - b) Use the arrows [16] and [17] to set the desired time value.
 - c) Confirm by pressing the setting button [14].
 - d) The unit starts timing when the set temperature value is reached.

During both time and temperature setting, holding down one of the arrows will rapidly increase/decrease the set value.

After one minute of inactivity, the unit will return to the basic operating mode.

CAUTION! If only the temperature value is set (time value = 0), the unit will operate in continuous mode.

3.3.1. Alarms:

If the current temperature value reaches the acceptable limit, the unit will emit a sound signal and the "ALM" signal lamp will light up. If you change the temperature setting, the lamp will remain lit but the unit will not beep.

The beep can be silenced at any time by pressing any of the buttons on the control panel.

3.3.2. AT/PID function

The AT/PID function is designed to eliminate the possibility of control errors based on parameter group settings (proportional, integral, differential).

1. Hold down the [15] key for 6 seconds. The "RUN/AT" indicator light starts flashing.
2. When finished, the "RUN/AT" indicator light will turn off and the unit will return to normal operation mode.
3. To abort the process, hold down the [15] button again.

3.3.3. Adjusting the internal parameters of the unit:

Hold down the "SET" button for 3 seconds. The display will show "Lc" indicating that a password must be entered. There are 3 passwords that launch 3 different setup menus. Confirm by pressing "SET" button.

PASSWORD: 3

Parameter indication	Parameter name	Parameter function description	(Range) Factory value
Lc-	Password	When "Lc=3", you	0

		can view and modify the parameter value.	
AL-	Over temperature Deviation alarm	When "Temperature measurement value>Temperature setting value+AL", the alarm light is on, the buzzer sounds (see V.3), and the heating output is disconnected.	(0~100°C) 5.0
T-	Control cycle	Heating control period.	(1~60s)
P1-	Proportional zone at low temperature	Time proportional action adjustment.	(1.0~Range value) 35.0
I1-	Integration time in low temperature area	Integral adjustment.	(1~1000s) 200
d1-	Differential time at low temperature	Differential action adjustment.	(0~1000s) 200
P2-	Proportional band	Time proportional action adjustment.	(1.0~Range value) 35.0
I2-	Integration time	Integral adjustment.	(1~1000s) 200
d2-	Differential time	Differential action adjustment.	(0~1000s) 200
dc-	Low temperature zone inflection point	When the set temperature is $\leq dc$, it belongs to the low temperature zone	(0~Range value)80.0
Pb-	Zero adjustment	Correct the error produced by the sensor (low temperature) measurement. Pb=actual temperature value-meter measurement value	(-50~50°C) 0
PK-	Full scale adjustment	Correct the error caused by the sensor (high temperature) measurement. PK=1000* (actual temperature value-instrument	(-999~999) 0

		measured value) / instrument measured value	
Et-	Timing function	When ET=0, there is no timing function; when it is 1, it starts timing when power on, when it is 2, it starts timing when it is set	(0~2) 2

Note 1: If the relay output is selected, the heating control period should be 20 seconds, the current output should be 1 second, and other models should be 5 seconds.

Note 2: If the FCH-600X and FCH-610X series, the factory setting of the timer function is 2, and the other models are 0.

PASSWORD: 9

Parameter indication	parameter name	Parameter function description	(Range) Factory value
Lc-	password	When "Lc=9", you can view and modify the parameter value.	0
Co-	Turn off heating output deviation	When "temperature measurement value \geq temperature set value + Co", turn off the heating output.	(0.0~50.0°C) 5.0
Hn-	Constant temperature timing	0: minute timing; 1: hour timing	(0~1) 0
En-	End of operation	En=0 Turn off the output at the end of operation; En=1 Continue to keep constant temperature after running;	(0~1) 0
Lt-	Maximum power output	The maximum power percentage of heating output;	(0~100)100
oP-	Gating function	0: Turn off the door open judgment function; 1: Turn on the door open judgment function.	(0~1) 1
rH-	Range value	Set according to the temperature	FCD:(0~400.0°C)250.0

		measurement range.	
ad-	mailing address	The communication address of this machine.	(1~32) 1

English name and parameter indication symbol comparison table

Parameter indicator	SP	St	Lc	AL	T	P	■	d
English name	SP	St	Lc	AL	T	P	■	d
Parameter indicator	Pb	Pk	Co	Hn	oP	rH	En	Lt
English name	Pb	Pk	Co	Hn	oP	rH	En	Lt

PASSWORD: 23 (invalid for this model)

3.4. Cleaning and maintenance

- a) Cleaning and maintenance
Pull the mains plug and let the unit cool down completely before cleaning, adjusting or replacing accessories and when the unit is not in use.
- b) Use only non-corrosive cleaning agents for cleaning the surfaces.
- c) After each cleaning, all the parts should be dried well before the unit is used again.
- d) Store the unit in a dry and cool place protected from moisture and direct sunlight.
- e) Do not spray the unit with a stream of water or immerse it in water.
- f) Perform regular inspections of the unit checking technical fitness and any damages.
- g) Use a soft, damp cloth for cleaning.
- h) Do not use sharp and/or metal objects (e.g. a wire brush or metal spatula) for cleaning, as these may damage the surface of the material from which the unit is made.
- i) Do not clean the unit with acidic substances, medical products, diluent, fuel, oil or other chemicals as this may cause damage to the unit.
- j) Wipe down the oven shelves after each use to get rid of debris and ventilate the compartment by leaving the door open after use.

DISPOSAL OF USED UNITS.

At the end of its useful life, this product should not be disposed of with normal household waste but should be taken to a collection point for the recycling of

electrical and electronic equipment. This is indicated by the symbol on the product, operating instructions or packaging. The materials used in this unit are recyclable according to their marking. You will be making an important contribution to protecting our environment by reusing, recycling or otherwise disposing of used units.

Your local administration will provide you with information about the appropriate disposal point for used units.

3.5. Troubleshooting

Problem	Possible cause	Action
No power supply	Incorrectly connected or damaged power cord	Check if the power cord is damaged - replace if necessary. If the power cord is not damaged, check the connection to the power source - remove and re-insert the power cord plug into the electrical outlet.
	The fuse is damaged	Replace
The temperature in the chamber does not go up	The set temperature is too low	Adjust the temperature value
	Heater/temperature controller/temperature sensor is defective	Check and replace the heater/temperature controller/temperature sensor
	Temperature sensor loose	Tighten the nut on the sensor
No vacuum in the chamber	Air inlet valve is loosened	Loosen the nut on the temperature sensor
	The vacuum pump is incorrectly connected or the hose is defective	Check the condition of the pump hose and its connection
	Vacuum pump defective (e.g. short circuit)	Connect another functioning vacuum pump (or repair the existing one)
	Chamber door not closed	Close chamber door
	Chamber door leakage	Replace the seal on the chamber door
	Vacuum meter defective	Replace the vacuum meter

The message "-----" appears in the display	Temperature sensor is defective	Replace the temperature sensor
Vacuum meter does not indicate "0" or values out of range	Vacuum pump short-circuit/defect	Repair/connect another functioning vacuum pump
	Vacuum meter defective	Replace the vacuum gauge if the vacuum pump connected is operational

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