USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value		
Product name	Floor scale		
Model	SBS-WCS-300	SBS-PS-300	SBS-PS-150
Power supply	Input: 100-240V AC/50-60Hz 0.5A Output: 12V 1A		
Battery type	Li-ion 1300 mAh/7.4V		
Protection class	II		
Protection rating IP	IP54		
Load [kg]	300	300	150
Maximum load [kg]	300.9	300.9	150.45
Minimum load [kg]	2	2	1
Resolution [g]	100	100	50
Overload	+9d		
Ambient temperature [°C] during operation	0-40		
Ambient humidity [%] during operation	25-95		
Zero Range	0-4% to full capacity		
Tare range	4% to full capacity		
Maximum safety load	120% FS		
Ultimate Load	150% FS		
Units	Kg; lb		
Automatic shut-down [s]	Can be set to OFF, or ON after 120s / 180s / 300s of inactivity		
Dimensions [mm]	1110 x 1000 x 1080	1000 x 500 x 70	1000 x 500 x 70
Weight [kg]	28.5	16.8	14.5

Parameter description	Parameter value
Product name	Floor scale
Model	SBS-PS-1800
Unit dimensions [length x Width x height + leg height (if applicable); mm]	965 x 510 x 70
Dimensions of the weighing platform [mm]	965 x 510
Weight [kg]	11
Minimum [g] / maximum load [kg]	1000 / 180
Measurement accuracy [kg]	0.05
Battery type	4 x AAA
Power supply parameters	AC/DC 9V 300mA
Available measurement units	kg; lb
Operating ambient temperature [°C]	0-40
Ambient humidity [%] during operation	25-95
Zero Range	0-4% to full capacity

Tare range	4% to full capacity
Display type	LCD, 5 digits
Length of cable connecting the display to the weighing platform [m]	2.5
Sleep mode	Automatic (delay time adjustable)

Parameter description	Parameter value	
Product name	Postal scale	
Model	SBS-PT-60SNL	SBS-PT-150SNL
Unit dimensions [length x width x height + leg height; mm]	355 x 405 x 55 + 5	
Length of cable connecting the display to the weighing platform (min-max; m)	1.1-2.2	
Weight [kg]	4.5	
Battery type	6 x AA	
Power supply	DC 9V	
Display type	LCD	
Operating temperature [°C]	1-40	
Minimum load [kg]	0.4	0.4
Maximum load [kg]	60	150
Maximum overload value [kg]	90	225
Accuracy [kg]	0.02	0.05
Measurement unit	kg; lb	
Automatic sleep mode	Can be set to OFF, or ON after 120s / 180s / 300s of inactivity	
Zeroing range [kg]	2.4	6

Parameter description	Parameter value	
Product name	Wall scale	
Model	SBS-PWS-150	SBS-PWS-300
Battery / power supply	Lithium battery (1300mAh) / DC (12V/1000mA)	
Maximum load [kg]	150	300
Minimum load [kg]	1	2
Graduation [g]	50	100
Dimensions of the weighing platform [mm]	550 x 450	
Weight [kg]	17.3	

1. General Description

The manual is intended to assist in safe and reliable use. The product is designed and manufactured strictly according to technical specifications using the latest technology and components and maintaining the highest quality standards.

**CAREFULLY READ AND UNDERSTAND THIS MANUAL
BEFORE STARTING THE WORK.**

To ensure the long and reliable operation of the device, make sure to operate and maintain it properly in accordance with the guidelines in this instruction manual. The technical data and specifications in this manual are up-to-date. The manufacturer reserves the right to make changes in order to improve the quality. Taking the technical progress and the possibility of reducing noise into account, the unit is designed and built in such a way that risks resulting from noise emissions are reduced to the lowest possible level.

Legend

	The product satisfies the relevant safety standards.
	Read instructions before use.
	The product must be recycled.
	WARNING! or CAUTION! or REMEMBER! Applicable to the given situation. (general warning sign)
	ATTENTION! Electric shock warning!
	Class II protection device with double insulation.
	For indoor use only.

PLEASE NOTE! Drawings in this manual are for illustration purposes only and in some details may differ from the actual product.

2. Usage safety

CAUTION! Read all safety warnings and instructions. Failure to follow the warnings and instructions may result in electric shock, fire, and/or serious injury or death.

The term “device” or “product” in the warnings and the description of the instructions refers to

Floor scale / Postal scale / Wall scale

2.1 Electrical safety

- a) The plug must fit the socket. Do not modify the plug in any way. Using original plugs and matching sockets reduces the risk of electric shock.
- b) Avoid touching earthed elements such as pipes, heaters, boilers and refrigerators. There is an increased risk of electric shock if the earthed device is exposed to rain, comes into direct contact with a wet surface or is operating in a damp environment. Water getting into the device increases the risk of damage to the device and of electric shock.
- c) Do not touch the device with wet or damp hands.
- d) Use the cable only for its designated use. Never use it to carry the device or to pull the plug out of a socket. Keep the cable away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cables increase the risk of electric shock.
- e) If using the device in a damp environment cannot be avoided, a residual current device (RCD) should be applied. The use of an RCD reduces the risk of electric shock.
- f) Do not use the device if the power cord is damaged or shows obvious signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer’s service centre.

- g) To avoid electric shock, do not immerse the cord, plug or device in water or other liquids. Do not use the device on wet surfaces.
- h) ATTENTION! DANGER TO LIFE! While cleaning, never immerse the device in water or other liquids.
- i) Do not use the device in very humid environments or in the direct vicinity of water tanks.
- j) Replace all batteries at the same time – Do not replace only a portion of the batteries as this may cause a malfunction.
- k) If the scale is not to be used for a long period of time, remove all batteries from the battery compartment to avoid leakage, which may cause damage to this scale.

2.2 Safety in the workplace

- a) If you are unsure about whether the product is operating correctly or if you find damage, please contact the manufacturer's service centre.
- b) Only the manufacturer's service centre may make repairs to the product. Do not attempt to make repairs yourself!
- c) Regularly inspect the condition of the safety labels. If the labels are illegible, they must be replaced.
- d) Please keep this manual available for future reference. If this device is passed on to a third party, the manual must be passed on with it.
- e) Keep packaging elements and small assembly parts in a place not available to children.
- f) Keep the device away from children and animals.
- g) If this device is used together with another equipment, the remaining instructions for use shall also be followed.

Remember! When using the device, protect children and other bystanders.

2.3 Personal safety

- a) The device is not designed to be handled by persons (including children) with limited mental and sensory functions or persons lacking relevant experience and/ or knowledge unless they are supervised by a person responsible for their safety or they have received instruction on how to operate the device.
- b) The device is not a toy. Children must be supervised to ensure that they do not play with the device.

2.4 Safe use of the product

- a) Do not use the device if the ON/OFF switch does not function properly (does not switch the device on and off). When not in use, store in a safe place, away from children and people not familiar with the device who have not read the user manual. The device may pose a hazard in the hands of inexperienced users.
- b) Keep the device in perfect technical condition. Device repair or maintenance should be carried out by qualified persons, only using original spare parts. This will ensure safe use.
- c) To ensure the operational integrity of the device, do not remove factory-fitted guards and do not loosen any screws
- d) Do not move, adjust or rotate the device in the course of work.
- e) Clean the device regularly to prevent stubborn grime from accumulating

- f) The device is not a toy. Cleaning and maintenance may not be carried out by children without supervision by an adult person.
- g) It is forbidden to interfere with the structure of the device in order to change its parameters or construction.
- h) Keep the device away from sources of fire and heat.
- i) Do not overload the device.
- j) Keep unused batteries out of reach of children.
- k) Before making a measurement, make sure that the device is placed on a stable surface and is not exposed to vibrations and sudden temperature changes, which may adversely affect the measurement results.
- l) Do NOT leave the load on the scale for long. This will decrease the scale's accuracy and shorten the load cell's life.
- m) Do NOT use the scale in thunder or rain.
- n) Avoid extreme temperatures. Do not place the unit in direct sunlight or near the air conditioning.
- o) Avoid places where moisture can lead to condensation. Avoid direct contact with water. Do not place the unit close to windows, doors or fans that may cause unstable results due to air currents.
- p) The scale should not be subjected to any kind of shakes, mechanical pressure or falling.
- q) Do not use cell phones or other radio frequency devices in the vicinity of the scales in use. It can cause the measurements to be inaccurate.
- r) It is recommended to place the weight in the centre of the scales to ensure the most accurate measurement.

ATTENTION! Despite the safe design of the device and its protective features, and despite the use of additional elements protecting the operator, there is still a slight risk of accident or injury when using the device. Stay alert and use common sense when using the device.


3. Use guidelines

The device is designed to measure the weight of elements placed on the platform. Do not exceed the maximum load of the scale - it may damage it!

The user is responsible for any damage resulting from misuse.

3.1 Device description

3.1.1. Product overview

1 - Control panel with display

2 - Weighing platform

3 - Mounting plate

3.1.2. Control panel with display


1 - Display

2 - Control buttons

3 - RS232 socket

4 - ON/OFF switch

3.1.3. Description of control buttons

- **“ZERO” button:** Resets the display to zero when the load is below 4% of full capacity, or when a negative reading is displayed.
- **“UNITS” button:** Switches the weight unit between “kg” and “lb”.
- **“TARE” button:** Clears the tare weight and resets the display to zero when the load is between 4% and 100% of full capacity.
- **“HOLD” button:** Holds (locks) or releases (unlocks) the displayed weight.

3.2 Preparing for use

3.2.1. Appliance location

The ambient temperature must not exceed 40°C and ambient humidity should not exceed 85%. Keep the unit away from any hot surfaces. Always operate the unit on a level, stable, clean, fireproof and dry surface and out of the reach of children and persons of impaired mental, sensory and intellectual functions. Maintain a minimum clearance of 10 cm from any wall of the unit. Place the unit in such a way that the mains plug can be reached at any time. Ensure that the power supply to the unit corresponds to that specified on the identification plate!

Insert the batteries into the compartment according to polarity. The chamber is located at the back of the scale.

3.2.2. Calibration

Notice: Calibration has been done from the factory. Do not recalibrate the appliance unless it is moved to a different altitude or gravity, or if it shows inaccurate readings.

Press and hold the **“ZERO”** and **“TARE”** buttons at the same time, and then turn on the power using the **“ON/OFF”** button. The display will show **“CAL”** followed by the AD value*. Press the **“TARE”** button, and the display will show the pre-set calibration weight value ($\geq 2/3$ of the full capacity). Put the calibration weight on the platform, then press the **“TARE”** button. The display will first show **“-----”**, and then the AD value. The calibration is now complete.

* AD value (Analog-to-Digital Conversion Value) is the raw digital signal from the load cell. It is mainly used during calibration and troubleshooting, not during normal weighing.

3.3 Installation of the unit

For models SBS-PWS-150 & SBS-PWS-300, the mounting plate of the scale should be attached to the wall with the use of appropriate mounting pins/anchors, depending on the material of the wall to which the scale will be attached. Remember that the wall (with mounting pins/anchors) needs to withstand the maximum weight load depending on the model (150 kg for model SBS-PWS-150 and 300 kg for model SBS-PWS-300).


Attach the scale using all mounting holes on the mounting plate.

3.4 Working with the unit

3.4.1. Before weighing


Whenever possible, please allow the scale to warm up for several seconds after first turning the power on, so that the scale will function properly and accurately.

3.4.2. Error messages

	Low battery.
Err-0	Overload.
Err-2	Exceeded zero track range when turning on the power.
Err-5	The scale is unstable when turning on the power.
Err-C	The calibration code entered is insufficient, or an operation error has occurred.

3.4.3. Setup menu (1st level)

1) Settings instructions (Diagram version)

2) Settings instructions (Text version)

- Turn on the unit while holding down the "ZERO" button. After turning on the unit and still holding the "ZERO" button, press the "UNITS" button three times. The unit will go to the list of available functions to adjust.
- To adjust the displayed parameter, press the "UNITS" button. The parameters are adjusted using the "TARE" button.


- To return to the main settings menu or skip a parameter and move to the next one, press the **“ZERO”** button.

3) Available settings:

Displayed messages	Corresponding function	Function adjustment options
CAPu	Maximum load unit	<ul style="list-style-type: none"> • “kg” • “lb”
CAP	Maximum load	<p>When the “UNITS” button is pressed, “0000” will be displayed.</p> <ul style="list-style-type: none"> • Press the “HOLD” button to mark the desired digit in the displayed number sequence. • Press “TARE” to increment the marked digit.
CALU	Calibration unit	<ul style="list-style-type: none"> • “kg” • “lb”
CAL	Calibration value	<p>When the “UNITS” button is pressed, “0000” will be displayed.</p> <ul style="list-style-type: none"> • Press the “HOLD” button to mark the desired digit in the displayed number sequence. • Press “TARE” to increment the marked digit.
DIVS	Division	<ul style="list-style-type: none"> • 0.5 • 1 • 2
ADC	Analog-to-Digital Converter number selection (in bits)	<ul style="list-style-type: none"> • 18 • 21

3.4.4. Setup menu (2nd level)

1) Settings instructions (Diagram version)

2) Settings instructions (Text version)

- Turn on the unit while holding down the **“ZERO”** button. When the unit is started, let go of the **“ZERO”** button.
- Press the **“TARE”** button three times. The unit will go to the list of available functions to adjust.
- To adjust the displayed parameter, press the **“UNITS”** button. The parameters are adjusted using the **“TARE”** button.
- To return to the main settings menu or skip a parameter and move to the next one, press the **“ZERO”** button.


3) Available settings:

Displayed messages	Corresponding function	Function adjustment options
A-OFF	Automatic shutdown (adjusting the time after which the function will be activated)	Automatic shutdown will occur after: <ul style="list-style-type: none"> • 120 seconds • 180 seconds • 300 seconds • “OFF” (disables the function)
b.L	Backlight	<ul style="list-style-type: none"> • “ON” – on • “OFF” – off • “STB” – on if there is an animal/object on the platform; turns off automatically 12 seconds after the reading is stabilized. • “AUTO” – on if there is an animal/object on the platform; off if the reading is 0.
UNIT	Weighing unit (selecting the weighing unit in this menu is responsible for its appearance when the “TARE” button is pressed during measurement)	<ul style="list-style-type: none"> • “ON kg” or “OFF kg” – weighing in kilograms on/off • “ON lb” or “OFF lb” – weighing in pounds on/off • “ON oz” or “OFF oz” – weighing in ounces on/off
Por.u	Default weighing unit	Selected weighing unit shown after startup: <ul style="list-style-type: none"> • “KG” – kilograms • “LB” – pounds • “OZ” – ounces • “KG.LB.OZ” – selection to be made after each startup
G.SW	Use of the gravitational force on the weighed object	<ul style="list-style-type: none"> • “ON” – on • “OFF” – off
Bmi.sw	Use of body mass index (BMI)	<ul style="list-style-type: none"> • “ON” – on • “OFF” – off
bAud	Baud rate	<ul style="list-style-type: none"> • 1200 • 2400 • 4800 • 9600
PAr	Parity bit control	<ul style="list-style-type: none"> • “nonE” – none • “odd” – odd • “EwEn” – even

Mod	Transmission mode	<ul style="list-style-type: none"> • “kEy” – sending data after receiving a command from the computer • “Stb” – data transfer after the reading stabilizes • “coN” – continuous data transfer • “OFF” – disables the function
Ani.sw	Animal weighing mode	<ul style="list-style-type: none"> • “ON” – enabled • “OFF” – disabled
FiLd.L	Filter Level	<ul style="list-style-type: none"> • “STr” – strong • “WEEK” – weak • “nor” – normal
Lock.r	Locking range	<ul style="list-style-type: none"> • 10 – 10d • 30 – 30d • 50 – 50d • “OFF” – disables the function
Pzr	Zero range when powered on - Percentage of weight on the scale, which will be equal to the value “0” after activation.	<ul style="list-style-type: none"> • 3 • 4 • 10 • 20
Kzr	Key zero range - Percentage of weight on the scale that will enable the zeroing function	<ul style="list-style-type: none"> • 3 • 4 • 5 • 10

3.4.5. RS232 transmission

1) RS232 interface: DB9

2) RS232 transmission

- Mode: Simplex asynchronous serial
- Baud rate: 9600
- Data bit: 8
- Parity bit: NONE
- Stop bit: 1
- Data format: ASCII

- 1 M= overload mode**
- 2 W= normal weighing mode**

: =3 weight information

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
W/M	:	±																CR	LF

3 Weight information 4 Unit formation

5 Example: 50kg

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
W	:	±							5	0	.	0	k	g				CR	LF

- 1- M=overload mode
- 2- W=normal weighing mode
- 3- Weight information
- 4- Unit formation
- 5- Example: 50kg

3.5 Cleaning and maintenance

3.5.1. General instructions

- a) Disconnect the unit from charging before any cleaning.
- b) Use only non-corrosive cleaning agents for cleaning the surfaces.
- c) Store the unit in a dry and cool place protected from moisture and direct sunlight.
- d) Do not spray the unit with a stream of water or immerse it in water.
- e) Perform regular inspections of the unit, checking technical fitness and any damages.
- f) Do not leave the batteries in the unit when it will not be used for a long time.
- g) Use a soft, damp cloth for cleaning.
- h) Do not use sharp and/or metal objects (e.g., a wire brush or metal spatula) for cleaning, as these may damage the surface of the material from which the unit is made.
- i) Do not clean the unit with acidic substances, diluent, fuel, oil or other chemicals as this may cause damage to the unit.

3.5.2. Safe removal of batteries and rechargeable batteries

Batteries are installed in the devices.

Remove used batteries from the device using the same procedure by which you installed them.

Recycle batteries with the appropriate organisation or company.

3.5.3. Disposing of used devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

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