USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value
Product name	Lower Temp Liquid Heating & Cooling Reaction Bath
Model	SBS-LCC-7000H
Voltage [V] / Frequency [Hz]	230~ / 50
Power [W]	2000
Temperature range [°C]	-20 to 100
Tank volume [L]	5
Pump flow [L/min]	20
Refrigerant	R290
Dimensions (Width x Length x Height) [mm]	385x505x1010
Weight [kg]	42.4
Protection class	I

Purpose

The product is used to provide precise temperature control for various scientific, industrial, and experimental processes. This equipment is designed to either heat or cool a circulating liquid within a controlled range, typically at lower temperatures, enabling it to maintain stable thermal conditions during chemical reactions, material testing, or other laboratory procedures.

Operation


WARNING

- Ensure that frozen liquid or water has been added according to the reaction temperature before starting the power.
- The power supply must meet the equipment's specifications and be properly grounded.
- Do not use the pump if the water in the container is too low, dry, or frozen.
- Maintain a minimum spacer distance of 700 mm and ensure that the vents around the machine are not blocked.
- Turn off the power before opening the shell. Only qualified personnel should perform maintenance, as the closed refrigeration system contains high-pressure gas.

Panel description

A: PV
B: SV

- A- PV: Real time temperature
- B- SV: Setting temperature
- 1- Window display
- 2- Temperature setting
- 3- Reduce / Stirring
- 4- Backlight / Plus
- 5- Cycle
- 6- Cooling
- 7- Heating
- 8- Power
- 9- Stirring speed
- 10- Swift / Auto-adjust

Panel Function

- 1) "Setting" Button : Used to set or view the temperature and related parameters.
- 2) "Plus/Backlight" Button : In the non-set state, pressing this button turns the backlight on or off. In the set state, pressing this button turns the display power supply on or off (this does not cut off the instrument's own power).
- 3) "Cooling" Button : Opens or closes the refrigeration function, indicating when refrigeration is allowed.
- 4) "Heating" Button : Opens or closes the heating function, indicating when heating is allowed.
- 5) "Cycle" Button : Opens or closes the loop output.

Adjustment Method

1. Press the **SET** key, and the **SV** window begins flashing.
2. Use the **▲▲▼** buttons simultaneously to set the desired temperature (from -20°C to 100°C), then press the **SET** button to confirm and end the process.

3. If you want to change other parameters, press the **SET** key for 5 seconds to enter the menu and adjust the parameters as required.
4. Once done, press the **SET** button for 5 seconds to exit the parameter setting mode.
5. The instrument will either automatically exit after 30 seconds and return to the working state, or if not saved, the modified data will be lost.

How to set the temperature

1. Power On:

Press the power button for 4 seconds. The controller will display "In-P" and show the range value. After 3 seconds, the window will enter the normal state.

2. Temperature Setting:

Press the "Settings" button to enter the temperature setting mode. The top row of the display will show the prompt "SP," and the bottom row will show the set temperature. Use the "Reduce" and "Plus" buttons to adjust the temperature value. After setting the desired temperature, press the "Settings" button again to save the value.

3. Alarm Conditions:

- Temperature Sensor Failure or Out of Range: The buzzer will sound, and the PV window will flash "-E0-".
- Overtemperature Alarm: The buzzer will sound, the "ALM" alarm lamp will light up, and the PV window will flash "-E1-".
- Low Temperature Alarm: The buzzer will sound, the "ALM" alarm lamp will light up, and the PV window will flash "-E2-".
- Silencing the Buzzer: Press the "Reduce" button to silence the buzzer when it sounds.

4. Automatic Exit:

If no button is pressed within 1 minute while in the setting mode, the controller will automatically return to the normal display state.

5. Temperature Control:

Two-position control is applied to maintain the temperature within a specific range. You can adjust the range by changing the **P** value. For example, if $SV = 100^{\circ}\text{C}$ and $P = 30^{\circ}\text{C}$, the actual temperature range will be $100 \pm \frac{30}{2}^{\circ}\text{C}$.

Process for Usage

1. Open the radiator cap and pour the appropriate medium (e.g., water or other liquids) as per your thermocouple temperature requirements. Ensure the liquid level is less than 2 cm from the bath edge.
2. Use insulated hoses to connect the machine's outlets and inlets accordingly.

3. Plug the machine into a properly grounded power outlet.
4. After turning on the machine's power, the display will show the ambient temperature. Turn on the cooling switch. The compressor will start working after a 3-minute delay.
5. When the cooling system is active, open the internal recycling valve if internal recycling is required.
6. Once the machine's temperature matches the set temperature, open the valve and switch on the cooling equipment, then proceed with your experiment.
7. After completing the experiment, turn off the LCOP button, cooling button, and finally the power button in sequence.
8. If the machine will not be used for an extended period, release the coolant liquid and clean the inside of the machine to prevent any buildup or damage.

Appendix

1. Operating Temperature:

- The temperature of the circulating liquid must not be lower than the lowest operating temperature of the machine's load.
- If water is used as the circulating fluid, the inlet liquid temperature must not exceed 40°C.
- If aqueous ethanol is used, the inlet liquid temperature must not exceed 10°C.

2. Low Antifreeze Ratios:

- Ethylene glycol solution: Glycol/Water = 55/45 for temperatures as low as -40°C.
- Glycerol solution: Glycerol/Water = 70/30 for temperatures as low as -30°C.

3. Antifreeze Density Requirements:

When using antifreeze, ensure the following density requirements are met:

- Ethylene glycol solution: P20 = 1079 kg/m³ at -40°C.
- Glycerol solution: P20 = 1183.6 kg/m³ at -30°C.
- If the viscosity of the ethylene glycol solution increases or the experiment is conducted below -40°C, use 95% pure alcohol as the antifreeze.

4. Minimum Temperature Requirement:

If using a non-frozen liquid instead of water, ensure that the minimum temperature is above 5°C.

Troubleshooting

Fault Symptoms	Failure Cause	Trouble Shooting
----------------	---------------	------------------

Machine can't work after turning the power	No power	Check power source
	Leakage protector activated	Reset the leakage protector
	Fuse burn -out	Replace the fuse
	Power cord disconnected	Tighten the power cord
Panel does not display properly	Wire loose	Test and secure the wire connections
	Platinum Short circuit or open circuit	Change platinum
	Panel fault	Test and repair the panel
Panel digital display flashes	Power supply interference	Remove interference, use a stable voltage
	Panel Fault	Test and repair the panel
Shell charged	Ground wire poor contact	Ensure proper connection of the ground wire
	Other wire disconnected	Check and reconnect the wire
	Heat insulation performance unstable	Replace the heater
No cooling	Cooling switch not turned on	Turn on the cooling switch
	Condenser not working	Replace the condenser
	Compressor does not start	Check the wiring and compressor
	Freon leak	Add Freon and try again
Cooling, but slowly	Freon Leak	Add Freon and try again
	Condenser not working	Replace the condenser
	Pipe blockage	Clear the pipeline

Disposing of Used Devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

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