USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

OVERHEAD STIRRER

DE	Produktname	Laborrührer
EN	Product name	Overhead stirrer
PL	Nazwa produktu	Mieszadło laboratoryjne
CZ	Název výrobku	Horní míchačka
FR	Nom du produit	Agitateur en tête
IT	Nome del prodotto	Agitatore sopraelevato
ES	Nombre del producto	Agitador de cabeza
HU	Termék neve	Fej feletti keverő
DA	Produktnavn	Overhead omrører
FI	Tuotteen nimi	Yläsekoitin
NL	Productnaam	Bovenroerder
NO	Produktnavn	Overhead omrører
SE	Produktnamn	Överliggande omrörare
PT	Nome do produto	Agitador de sobrecabeça
SK	Názov produktu	Horná miešačka
BG	Име на продукта	Горна бъркалка
EL	Όνομα προϊόντος	Κορυφαίος αναδευτήρας
HR	Naziv proizvoda	Gornje miješalo
LT	Produktu pavadinimas	Viršutinis maišytuvas
RO	Numele produsului	Agitator superior
SL	Ime izdelka	Zgornje mešalo
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-MR-2540
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated for your convenience using machine translation. Reasonable efforts have been made to provide an accurate translation; however, no automated translation is perfect nor is it intended to replace human translators. The official User Manual is the English version. Any discrepancies or differences created in the translation are not binding and have no legal effect for compliance or enforcement purposes. If any questions arise related to the accuracy of the information contained in the User Manual, please refer to the English version of those contents which is the official version.

FOR YOUR SAFETY AND CORRECT OPERATION, PLEASE CAREFULLY READ THIS INSTRUCTION BEFORE USING IT.

NOTE: All of the information reported herein is based on data available at the moment of printing. The factory reserves the right to modify its own products at any moment without notice and incur any sanction. So, it is suggested to always verify possible updates.


1. TECHNICAL DATA

Parameter description	Parameter value
Product name	Overhead Stirrer
Model	SBS-MR-2540
Voltage/frequency	230 V / 50 Hz
Power	80 W
Platform size cm	20 x 30
Max rotation speed rev/min	100 – 2500 rpm
IP protection class	IP42



2. PRODUCT OVERVIEW

177 cm

3. General product information


Overhead stirrer. Users must read this manual carefully, follow the instructions and procedures, and read all warnings when using this instrument.

Receiving Inspection

Carefully unpack the device and check for any damage caused during transport. If this happens, contact the manufacturer for technical support.


Important Notice!!!

When installing, don't forget to install the metal gasket when tightening the two screws on the base plate and the back of the main body.

A Metal gasket

Control Panel Instructions

1. “Run/Stop Key”: Press the key, you can run or stop the controller.
2. “Code Switch”: Setting the speed value and inter parameters.

4. BUTTONS AND FUNCTIONS

On power-up, the controller display **FS2P** on the upper display window and **S1.1** on the lower display window. After three seconds, the controller will be into normal display status.

Speed set and view (Note: Refer to the internal parameters table-4: nod parameter)

1) When speed Set Mode is **0** :

1> When the timing function is disabled (internal parameters table-3:

ET=0)

Press the **Code Switch** to enter setting mode. The controller displays - **SP-** on the upper display window and speed set value on the lower display window. Rotate the **Code Switch** to change the set value that you want. Press **Code Switch** again, the controller will exit setting mode and set value will be auto saved.javascript:void(0);

2> When the timing function is allowed

- ① When the motor is allowed forward or reverse running

(Internal parameters table-2: **Fr=0 or Fr=1**)

Press the **Code Switch** to enter setting mode. The controller displays **-SP-** on the upper display window and speed set value on the lower display window. Rotate the **Code Switch** to change the set value that you want. Press **Code Switch** again, the controller displays **-ST-** on the upper display window and time set value on the lower display window. Rotate the **Code Switch** to change the set value that you want. Press **Code Switch** again, the controller will exit setting mode and set value will be auto saved.

javascript:void(0);

- ② When the motor is allowed forward and reverse running

(Internal parameters table-2: **Fr=2**)

Press the **Code Switch** to enter setting mode. The controller displays **-SP-** on the upper display window and speed set value on the lower display window. Rotate the **Code Switch** to change the set value that you want.

Press **Code Switch** again, the controller displays **-ST-** on the upper display window and total time set value on the lower display window.

Press **Code Switch** again, the controller displays **-FT-** on the upper display window and forward time set value on the lower display window.

Press **Code Switch** again, the controller displays **-PT-** on the upper display window and stop time set value on the lower display window.

Press **Code Switch** again, the controller displays **-rT-** on the upper display window and reverse time set value on the lower display window.

When the total time is **0**, it indicates there is no timing function, the controller continuous running. When the total time is not **0**, after the timer stops running, the lower display area displays **End**. Press **Run/Stop** key, you can restart the controller running.

When speed Set Mode is **1** :

Rotate the **Code Switch** to enter setting mode. The controller displays **-SP-** on the upper display window and speed set value on the lower display window. Rotate the **Code Switch** to change the set value that you want. After 1 second, the controller will exit setting mode, and speed set value will be auto saved. On stop state, Press **Run/Stop** key, **Run** identifier lights, the controller will start running. On running state, Press **Run/Stop** key, **Stop** identifier lights, the controller will stop running.

When speed Set Mode is **2** :

Rotate the **Code Switch** to enter setting mode. The controller displays **-SP-** on the upper display window and speed set value on the lower display window. Rotate the **Code Switch** to change the set value that you want. After 1 second, the controller will exit setting mode and set value will be auto saved.

When the controller happen stall, Hall logical error, under bus voltage, over bus voltage, communication failure, the controller will automatically stop running and displays the Fault code.

Fault Code Description :

- 1) The upper window displays **Er-1 : Power module Fault**
- 2) The upper window displays **Er-2 : Stall Fault**
- 3) The upper window displays **Er-3 : Hall Fault**
- 4) The upper window displays **Er-4 : Under bus voltage Fault**
- 5) The upper window displays **Er-5 : Over bus voltage Fault**
- 6) The upper window displays **Er-6 : Communication failure Fault**

Note: All the internal parameters have been adjusted when factory test. Forbidden to modify them except Sensor Correction parameter.

Press **Code Switch** for three seconds, the controller display **Lc** on the upper display window, password value is displayed on the lower display window. Rotate the **Code Switch** to change the password value. Press **Code Switch** again, if the password value is correct, controller will automatically enter the internal speed parameters state. Press **Code Switch** for three seconds, you can exit the internal speed parameters state, and each parameter will be auto saved.

Internal parameters table-1

Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=3 , enter the parameters list	0
Pd-	Proportional	Adjustment of proportional function	(1~100) 10
Id-	Integration	Adjustment of integration function	(1~100) 5
InT	Acceleration time	The time of from 0 to maximum speed	(1~60) 10
dET	Deceleration time	The time of from maximum speed to 0	(1~60) 10
SdL	Minimum set value	The speed minimum set value	(20~6000) 20
SdH	Maximum set value	The speed maximum set value	(20~6000) 2000

Internal parameters table-2

Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=9, enter the parameters list	0
EAr	Gear ratio	Large diameter gear/pinion diameter	(1.0~10.0) 1.0
PoL	Motor poles	The motor poles	(1~32) 4
dIF	Motor reference direction	0 : clockwise rotation 1 : counterclockwise rotation	(0~1) 0
FdS	Speed feedback	Speed feedback value	(0.1~10.0) 1.0
FdC	Current feedback	Speed feedback value	(0.1~10.0) 1.0
FrE	Carrier Frequency	PWM carrier Frequency	(5~15) 15
Po-	Motor Power	Motor Power	(1~750) 40
CL-	Over current multiples	Allow rated current multiples	(1.0~10.0) 5.0
Fr-	Motor rotation direction	0 : Motor runs only clockwise 1 : Motor runs only counterclockwise 2 : Motor runs only forward or reverse	(0~2) 0
db-	Speed dead zone	The dead zone of speed	(0~100) 5

Internal parameters table-3

Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=101, enter the parameters list	0
Pb-	Zero correction	Update the measurement error Pb= actual value – measured value	(-99.9~99.9) 0.0
PK-	Full correction	Update the measurement error PK=1000× (actual value—measured value) /measured value	(-999~999) 0
SPL	Minimum set value	The temperature minimum set value	(-150.0~0.0) -150.0 °C
SPH	Maximum set value	The temperature maximum set value	(0.0~300.0) 300.0 °C
dot	Decimal display	0 : Disable decimal display 1 : Enable decimal display	(0~1) 1
TSn	Temperature sensor selection	0 : Disable temperature sensor 1 : Enable temperature sensor	(0~1) 0

Internal parameters table-4

Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=109, enter the parameters list	0
ET-	Timing selection	0 : Disable timing 1 : Enable timing	(0~1) 0
Hn-	Timing Mode	0 : Minute timing 1 : Hour timing	(0~1) 0
rT-	Time Correction	Correction value = 【Running time second) -Actual time (second) 】 *10 / Actual time (minute)	(-999~999) 0

Internal parameters table-5


Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=103, enter the parameters list	0
nod	Speed setting mode	0 : Press the Code Switch to set speed, Press Run/Stop key to control motor run or stop 1 : Rotate the Code Switch to set speed, Press Run/Stop key to control motor run or stop 2 : Rotate the Code Switch to set speed, according to set speed to control motor run or stop	(0~2) 0

Internal parameters table-6

Parameter indicator	Name	Instruction of the function	Setting range
Lc-	password	Lc=110, enter the parameters list	0
rEST	Restore factory parameter	Restore parameters to default value	(0~1) 0


Wiring

Driver Board Wiring

- A Thick Wire
- B Thin Wire

Display Board Wiring

A PT100
B Sensor

5. Faults

Instruments can't be power ON

1. Check whether the power cable is plugged
 2. Check whether the fuse is broken or loose
- Fault in power on self-test
3. Switch OFF the unit, then switch ON and reset the instruments to factory default setting.
 - Temperature cannot reach set point
 4. Check whether the safety temperature value is set too low
 - Stir speed cannot reach set point
 5. Excessive medium viscosity may cause abnormal speed reduction of the motor
 - Heating cannot be started after setting the temperature, or stirring cannot be started when adjusting the control knob.
 6. Check the control panel has damages which may have arisen during transport.

If these faults are not resolved, please set the instruments to factory

default setting, or take the unit to your technical service center, or contact the manufacturer.

6. Maintenance and Cleaning

- Proper maintenance can keep instruments working in a good state and lengthen their lifetime.
- Be careful not spray the cleanser into the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use cleanser that we advised as below:

<u>Dyes</u>	<u>Isopropyl alcohol</u>
<u>Construction materials</u>	<u>Water containing tenside/ isopropyl alcohol</u>
<u>Cosmetics</u>	<u>water containing tenside/ isopropyl alcohol</u>
<u>Foodstuffs</u>	<u>Water containing tenside</u>
<u>Fuels</u>	<u>Water containing tenside</u>

- Wear the proper protective gloves during cleaning of the instrument.
- Before using other methods for cleaning or decontamination, the user must contact the manufacturer to ascertain that this method does not destroy the instrument.
- The enamel makes the hotplate easier to care for and more resistant to acids and bases. Because of it, however, the heating plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.
- The instrument must be cleaned and put into the initial packaging carton before being sent to service for repair, avoiding the contamination of hazardous.
- Use the instrument in a dry clean room and temperature stable environment.

7. Safety

	Connect the device to an earthed power supply to ensure safety of machine and experiment; connect the power as the machine required.
	This equipment is forbidden to use in inflammable and explosive, poisonous and strong corrosive experiments.
	Make sure horizontal installation.
	Non-professionals are not allowed to disassemble and repair this machine.
	Pay attention to the set temperature while dealing with the inflammable matters.
	Make sure you dry the resin container, if the temperature is setting too high by accident, the container would be dissolved and then fall on the heater to cause fire.
	Overfilled sample will lead to overheating of working room under parts, which will dissolve the inflammable material and cause fire.
	While the machine is working, don't touch the top, window and exhaust port of the device to protect from high-temperature burns.
	Read the instruction book before operation.

- When working, wear the personal guard to avoid the risk from:
 - Splashing and evaporation of liquids
 - Release of toxic or combustible gases.
- Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface, do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- Gradually increase the speed, reduce the speed if:
 - The stirring bar breakaway because of too high speed
 - The instrument is not running smoothly, or container moves on the base plate.
- Temperature must always be set to at least 25°C lower than the fire points of the media used.
- Beware of hazards due to:
 - Flammable material or media with a low boiling temperature
 - Overfilling of media
 - Unsafe container
- Process pathogenic materials only in closed vessels.
- Check the instrument and accessories beforehand for damage each time you use them.

Do not use damaged components. Safe operation is only guaranteed with the accessories described in the “Accessories” chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.

- The instrument can only be disconnected from the main power supply by pulling out the mains plug or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the mains power supply cable does not touch the plate. Do not cover the device.

Storage and transportation

- Keep it in dry and clean room with good ventilation and no corrosive gas
- prevent it from getting wet by the rain and avoid violent collision in transportation.

Working condition

Ambient temperature: 5~40°C.

Ambient humidity: ≤90%.

Voltage: 220V±10%, 50/60Hz or 110V±10%, 50/60Hz

**Note:**

If there is any apparent damage to the system,
Please do not plug it into the power line.

Trial Runs

- Make sure the required operating voltage and power supply voltage match.
- Ensure the socket must be earthed reliably.
- Ensure the power be off
- Plug in the power cable, ensure the power is on and begin initializing.
- Add the medium into the vessel with a stirring bar if with the magnetic stirrer function.
- Put the vessel on the work plate.
- Set the rated stirring speed and start stirring.
- Observe LCD display.
- Stop the stirring functions.

If these operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact the manufacturer for technical support.

PART LIST

No.	Description	Qty	Remark
1.	Main unit	1	
2.	Power cable	1	
3.	Rack with rods	1	

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