## USER MANUAL

Bedienungsanleitung  
Instrukcja obsługi  
Návod k použití  
Manuel d'utilisation  
Istruzioni per l'uso  
Manual de instrucciones  
Használati útmutató  
Brugsanvisning  
Käyttöohje  
Gebruiksaanwijzing  
Bruksanvisning  
Instruções de utilização  
Používateľská príručka  
Ръководство за употреба  
Οδηγίες χρήσης  
Upute za uporabu  
Naudojimo instrukcija  
Manual de utilizare  
Navodila za uporabo

---

---

THERMOSTATIC WATER BATHS

<b>DE</b>	Produktname	Thermostatisches Wasserbad
<b>EN</b>	Product name	Thermostatic water bath
<b>PL</b>	Nazwa produktu	Termostatyczna łaźnia wodna
<b>CZ</b>	Název výrobku	Termostatická vodní lázeň
<b>FR</b>	Nom du produit	Bain-marie thermostatique
<b>IT</b>	Nome del prodotto	Bagno d'acqua termostatico
<b>ES</b>	Nombre del producto	Baño de agua termostático
<b>HU</b>	Termék neve	Termosztatikus vízfürdő
<b>DA</b>	Produktnavn	Termostatisk vandbad
<b>FI</b>	Tuotteen nimi	Termostaattinen vesikylpä
<b>NL</b>	Productnaam	Thermostatisch waterbad
<b>NO</b>	Produktnavn	Termostatisk vannbad
<b>SE</b>	Produktnamn	Termostatisk vattenbad
<b>PT</b>	Nome do produto	Banho-maria termostático
<b>SK</b>	Názov produktu	Termostatická vodná lázeň
<b>BG</b>	Име на продукта	Термостатична водна баня
<b>EL</b>	Όνομα προϊόντος	Θερμοστατικό λουτρό νερού
<b>HR</b>	Naziv proizvoda	Termostatska vodena kupelj
<b>LT</b>	Produktu pavadinimas	Termostatinė vandens vonia
<b>RO</b>	Numele produsului	Baie de apă termostatică
<b>SL</b>	Ime izdelka	Termostatska vodna kopel
<b>DE</b> Modell   <b>EN</b> Product model   <b>PL</b> Model produktu   <b>CZ</b> Model výrobku   <b>FR</b> Modèle   <b>IT</b> Modello   <b>ES</b> Modelo   <b>HU</b> Modell   <b>DA</b> Model   <b>FI</b> Tuotteen malli   <b>NL</b> Productmodel   <b>NO</b> Produktmodell   <b>SE</b> Produktmodell   <b>PT</b> Modelo do produto   <b>SK</b> Model   <b>BG</b> Модел на продукт   <b>EL</b> Μοντέλο προϊόντος   <b>HR</b> Model proizvoda   <b>LT</b> : Gaminio modelis   <b>RO</b> : Model de produs   <b>SL</b> : Model izdelka		SBS-TWB-2000 SBS-TWB-2002 SBS-TWB-2004 SBS-TWB-2006
<b>DE</b> Hersteller   <b>EN</b> Manufacturer   <b>PL</b> Producent   <b>CZ</b> Výrobce   <b>FR</b> Fabricant   <b>IT</b> Produttore   <b>ES</b> Fabricante   <b>HU</b> Termelő   <b>DA</b> Producent   <b>FI</b> Valmistaja   <b>NL</b> Producent   <b>NO</b> Produsent   <b>SE</b> Tillverkare   <b>PT</b> Fabricante   <b>SK</b> Výrobca   <b>BG</b> Производител   <b>EL</b> Κατασκευαστής   <b>HR</b> Proizvođač   <b>LT</b> Gamintojas   <b>RO</b> Producător   <b>SL</b> Proizvajalec		expondo Polska sp. z o.o. sp. k.
<b>DE</b> Anschrift des Herstellers   <b>EN</b> Manufacturer Address   <b>PL</b> Adres producenta   <b>CZ</b> Adresa výrobce   <b>FR</b> Adresse du fabricant   <b>IT</b> Indirizzo del produttore   <b>ES</b> Dirección del fabricante   <b>HU</b> A gyártó címe   <b>DA</b> Producentens adresse   <b>FI</b> Valmistajan osoite   <b>NL</b> Adres producent   <b>NO</b> Produsentens adresse   <b>SE</b> Tillverkarens adress   <b>PT</b> Endereço do fabricante   <b>SK</b> Adresa výrobcu   <b>BG</b> Адрес на производителя   <b>EL</b> : Διεύθυνση κατασκευαστή   <b>HR</b> Adresa proizvođača   <b>LT</b> Gamintojo adresas   <b>RO</b> Adresa producătorului   <b>SL</b> Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra   Poland, EU

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*This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via [info@expondo.com](mailto:info@expondo.com).*

## Technical data

Parameter description	Parameter value			
Product name	Thermostatic water bath			
Model	SBS-TWB-2000	SBS-TWB-2002	SBS-TWB-2004	SBS-TWB-2006
Rated voltage [V~] / Frequency [Hz]	220-240 / 50-60			
Rated power [W]	400	500	1000	1620
Tem-range [°C]	RT+5~100			
Number of chambers	1	2	4	6
Dimensions [Width x Depth x Height; mm]	170 x 154 x 210	318 x 168 x 210	350 x 318 x 210	524 x 322 x 210
Weight [kg]	3.3	4.5	6	7.5

SBS-TWB-2000	SBS-TWB-2002
	
SBS-TWB-2004	SBS-TWB-2006
	

## I. Summary

The water bath is suitable for applications such as distillation, concentration, drying, and thermostatic heating. It is widely used in medical units, universities, colleges, scientific research facilities, and laboratories. Additionally, it is applicable in industrial and mining enterprises, including chemical, printing, dyeing, and pharmaceutical industries.

## II. Structure features

1. **Construction:** The enclosure of the water bath is formed and machined using high-quality steel plates. The surface is treated with a static electric spray process, enhancing durability and sturdiness. The inner container is reinforced with stainless steel stretching for improved durability.
2. **Corrosion resistance:** Both the liner and upper cover are made from high-quality stainless steel plates, providing strong resistance to corrosion.
3. **Efficient heating:** A U-shaped heating pipe is used for direct water heating, allowing for a quick temperature rise with minimal thermal loss.

- 
4. **Temperature control:** Equipped with a single-row digital display or intelligent temperature controller, the device allows for simple operation and delivers a stable and favorable heating effect.

### III. Working conditions

1. **Temperature range:** 5°C to 40°C
2. **Relative humidity:** Less than 85% RH
3. **Power supply:** Voltage 220–240V, frequency 50–60 Hz
4. **Environmental conditions:** The equipment should be kept free from strong vibrations and corrosive gases.

### IV. Precautions

1. **Water level:** Before use, add water to a depth of 50 mm, ensuring it covers the clapboard. Then connect the unit to a power supply and begin heating. Operating with insufficient water is prohibited.
2. **Safety:** Avoid touching the heating pipe with your hands during operation to prevent scalding.
3. **Post-use maintenance:** After each use, discharge the water, dry the unit, and keep it clean to extend its service life.

## SBS-TWB-2000

1	PV: Measure value
2	SV: Set Value
3	LED indicators
4	SET key
5	Shift key
6	Down key
7	Up key

1	PV: Measure value
2	SV: Set Value
3	LED indicators
4	SET key
5	Shift key
6	Down key
7	Up key
8	Power switch
9	Drain switch

## V Instructions for use

1. **Positioning:** Place the instrument on a flat, horizontal surface.
2. **Filling the water tank:** Open the cover and add pure or distilled water to the water tank. Ensure the water level is above the heating pipe and temperature sensor.
3. **Power connection:** Connect to an appropriate power source, then turn on the switch to activate the power supply.
4. **Temperature display:** The upper row of the display shows the current (test) temperature, while the lower row shows the set temperature.
5. **Adjusting temperature:** Briefly press the "Set" key to enter the settings mode. Use the Shift ←, Up ↑ and Down ↓, keys to adjust the temperature.
6. **Setting the timer:** Press the Upper ↑ key to set the time. When "ST" appears in the upper row and "OFF" in the lower row, press the Upper ↑ key again to display the time in the lower row. Use the Shift ← key to adjust the time in minutes. There are two timing modes available:

- 
- **Constant temperature timing:** Starts after the temperature stabilizes.
  - **Set-time timing:** Starts after setting the time.

Once the set time is reached, heating will stop. To restart the operation, turn off and on the power switch.

7. **Automatic tuning function:** If the test temperature fluctuates, activate the self-tuning function to adjust it automatically. Press the Plus ↑ and Minus ↓ keys until the indicator light turns on.

## VI Troubleshooting

Failure	Cause	Handling method
No power supply	<ol style="list-style-type: none"> <li>1. Bad contact between plug and socket</li> <li>2. The fuse is burnt.</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace the plug or socket tube.</li> <li>2. Replace the fuse with same specification.</li> </ol>
No temp. rise	<ol style="list-style-type: none"> <li>1. The temp. Controller is broken</li> <li>2. The sensor is broken</li> <li>3. The set temperature is lower than water temperature</li> <li>4. The heating pipe is burnt</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace the instrument</li> <li>2. Replace the sensor</li> <li>3. Reset the temperature</li> <li>4. Replace the heating pipe</li> </ol>
The big difference between display temp. and actual temp.	<ol style="list-style-type: none"> <li>1. The temp. controller</li> <li>2. The temp. sensor is broken.</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace the temp. controller</li> <li>2. Replace the temp. sensor.</li> </ol>

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SBS-TWB-2000	SBS-TWB-2002
 A white rectangular water bath with a digital display showing 22.1 and a green power button on the left side.	 A white rectangular water bath with a digital display showing 22.5 and a green power button on the left side.
SBS-TWB-2004	SBS-TWB-2006
 A white rectangular water bath with a digital display showing 22.5 and a green power button on the left side.	 A white rectangular water bath with a digital display showing 22.7 and a green power button on the left side.

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