USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

ELECTROMOTION PIPETTE CONTROLLER

DE	Produktname	ELEKTROMOTORISCHER PIPETTENCONTROLLER
EN	Product name	ELECTROMOTION PIPETTE CONTROLLER
PL	Nazwa produktu	PIPETOR AKUMULATOROWY
CZ	Název výrobku	ELEKTROMOTORICKÝ PIPETOVÝ KONTROLER
FR	Nom du produit	CONTROLEUR DE PIPETTE ELECTROMECHANIQUE
IT	Nome del prodotto	CONTROLLER PER PIPETTE ELETTROMECCANICO
ES	Nombre del producto	CONTROLADOR DE PIPETA ELECTROMECAÁNICO
HU	Termék neve	ELEKTROMOS PIPETTAVEZÉRLŐ
DA	Produktnavn	ELEKTROMOTORISK PIPETTEKONTROL
FI	Tuotteen nimi	FSÄHKÖMOOTTORINEN PIPETTIOHJAIN
NL	Productnaam	ELEKTROMOTORISCHE PIPETCONTROLLER
NO	Produktnavn	ELEKTROMOTORISK PIPETTEKONTROLLER
SE	Produktnamn	ELEKTROMOTORISK PIPETTKONTROLL
PT	Nome do produto	CONTROLADOR DE PIPETA ELETROMECAÁNICO
SK	Názov produktu	ELEKTROMOTORICKÝ PIPETOVÝ OVLÁDAČ
BG	Име на продукта	ЕЛЕКТРОМОТОРЕН ПИПЕТЕН КОНТРОЛЕР
EL	Όνομα προϊόντος	ΗΛΕΚΤΡΟΚΙΝΗΤΟΣ ΕΛΕΓΚΤΗΣ ΠΙΠΕΤΑΣ
HR	Naziv proizvoda	ELEKTROMOTORNI KONTROLER PIPETE
LT	Produkto pavadinimas	ELEKTROMOTORINIS PIPETĖS VALDIKLIS
RO	Numele produsului	CONTROLER PIPETĂ ELECTROMECHANIC
SL	Ime izdelka	ELEKTROMOTORNI KRMILNIK PIPETE
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LAB-110
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated for your convenience using machine translation. Reasonable efforts have been made to provide an accurate translation; however, no automated translation is perfect nor is it intended to replace human translators. The official User Manual is the English version. Any discrepancies or differences created in the translation are not binding and have no legal effect for compliance or enforcement purposes. If any questions arise related to the accuracy of the information contained in the User Manual, please refer to the English version of those contents which is the official version.


Technical data

Parameter description	Parameter value
Product name	Electromotion pipette controller
Model	SBS-LAB-110
Input	230 V~/50 Hz
Battery	Lithium 2.4V DC / 700 mAh
Battery life (fully charged)	~8 h
Capacity	0.1 – 100 ml
Display	LCD
Suction speed	25ml < 5s
Dimensions (LxWxH)	46 x 44.5 x 30.2 cm
Weight	0.25 kg

1. Product overview

1. Suction tube adapter
2. Suction button
3. Drainage button
4. LCD Display
5. Speed control buttons
6. Charging port

2. Using the device

Application

This device is used for moving liquid from the original container to another container within a certain volume range.

Device operating environment

- For indoor use only
- Altitude: 2000 m
- The working temperature range of the device is +5°C~+40°C.
- The applicable temperature range of the device is 80%.
- Ensure that the area around the device is free from airflow or vibrations that could impact its performance.

Usage

Gently push the suction tube into the silicone adapter, and carefully set it in place. For proper insertion of suction tube, separately hold the mouth of the adapter tube and the upper end of the suction tube.

By using the speed control button below the LCD screen, you can set the desired suction and discharge speeds. The suction and discharge speeds have 6 steps and can be adjusted with the speed control buttons marked with up and down arrows. The LCD display shows the selected speed. In addition, gravity can also be used for draining liquids.

Select the appropriate speed, press the suction button, and draw the liquid slightly over the desired volume. Remove the suction tube from the liquid in the container. Then press the liquid feeding button to discharge part of the liquid, so that the bottom of the liquid in the suction tube reaches the desired level.

Two different drainage modes:

Gravity drainage: touch the drainage button and press it slightly. Like this, the drainage function will act under its own weight, and the liquid will be discharged naturally.

Electric drainage: pressing the discharge button further will trigger the micro pump to speed up the discharge speed.

(Note: in either case, the discharge speed can be sensitively adjusted by finger pressure.)

Safety Precautions

Charge the device for 10 hours before its first use to ensure a sufficient power supply. For future charges, the duration is 2-3 hours. Before plugging the AC charger into the power outlet, ensure that the input voltage and connector type are compatible with local standards.

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