USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

MAGNETIC HOTPLATE STIRRER

DE	Produktname	Magnetrührer
EN	Product name	Magnetic hotplate stirrer
PL	Nazwa produktu	Mieszadło magnetyczne z podgrzewaniem
CZ	Název výrobku	Magnetické míchačky s ohřevem
FR	Nom du produit	Agitateur magnétique chauffant
IT	Nome del prodotto	Agitatore magnetico riscaldante
ES	Nombre del producto	Agitador magnético con placa calefactora
HU	Termék neve	Mágneses fűtőlap keverő
DA	Produktnavn	Magnetisk varmeplade omrører
FI	Tuotteen nimi	Magneettinen kuumalevyn sekoitin
NL	Productnaam	Magnetische verwarmingsplaatmenger
NO	Produktnavn	Magnetisk varmeplateomrører
SE	Produktnamn	Magnetisk varmplatta omrörare
PT	Nome do produto	Agitador magnético com placa de aquecimento
SK	Názov produktu	Magnetická miešačka s ohrevom
BG	Име на продукта	Магнитна бъркалка с нагревателна плоча
EL	Όνομα προϊόντος	Μαγνητικός αναδευτήρας με θερμαινόμενη πλάκα
HR	Naziv proizvoda	Magnetna mešalna plošča s segrevanjem
LT	Produktu pavadinimas	Magnetinė kaitinimo plokštelė su maišykle
RO	Numele produsului	Agitator magnetic cu plită de încălzire
SL	Ime izdelka	Magnetni mešalnik s grelno ploščo
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LAB-112
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

1. Technical data

Parameter description	Parameter value
Product name	Magnetic Hotplate Stirrer
Model	SBS-LAB-112
Rated voltage [V~] / Frequency [Hz]	230 / 50
Rated power [W]	400
Rotation Speed Range [rpm]	50 - 1500
Dimensions [Width x Depth x Height; mm]	240 x 150 x 400
Weight [kg]	1.5

2. General description

The user manual is designed to assist in the safe and trouble-free use of the device. The product is designed and manufactured in accordance with strict technical guidelines, using state-of-the-art technologies and components. Additionally, it is produced in compliance with the most stringent quality standards.

DO NOT USE THE DEVICE UNLESS YOU HAVE THOROUGHLY READ AND UNDERSTOOD THIS USER MANUAL.

To increase the product life of the device and to ensure trouble-free operation, use it in accordance with this user manual and regularly perform maintenance tasks. The technical data and specifications in this user manual are up to date. The manufacturer reserves the right to make changes associated with quality improvement. The device is designed to reduce noise emission risks to a minimum, taking into account technological progress and noise reduction opportunities.

2.1 Legend

Icon	Description
	The product satisfies the relevant safety standards.
	Read instructions before use.
	The product must be recycled.
	WARNING! or CAUTION! or REMEMBER! Applicable to the given situation. (general warning sign)
	Wear protective goggles.
	Wear protective gloves.
	Wear protective clothing.
	ATTENTION! Electric shock warning!
	ATTENTION! Fire hazard - flammable materials!
	WARNING! Toxic substances, danger of poisoning!
	ATTENTION! Hot surface, risk of burns!
	Only use indoors.

PLEASE NOTE! DRAWINGS IN THIS MANUAL ARE FOR ILLUSTRATION PURPOSES ONLY AND IN SOME DETAILS MAY DIFFER FROM THE ACTUAL PRODUCT.

3. Usage safety

ATTENTION! READ ALL SAFETY WARNINGS AND ALL INSTRUCTIONS. FAILURE TO FOLLOW THE WARNINGS AND INSTRUCTIONS MAY RESULT IN ELECTRIC SHOCK, FIRE AND/OR SERIOUS INJURY OR EVEN DEATH.

The terms "device" or "product" are used in the warnings and instructions to refer to:

Magnetic hotplate stirrer

3.1 Electrical safety

- a) The plug must fit the socket. Do not modify the plug in any way. Using original plugs and matching sockets reduces the risk of electric shock.
- b) Avoid touching earthed elements such as pipes, heaters, boilers and refrigerators. There is an increased risk of electric shock if the earthed device is exposed to rain, comes into direct contact with a wet surface or is operating in a damp environment. Water getting into the device increases the risk of damage to the device and of electric shock.
- c) Do not touch the device with wet or damp hands.
- d) Use the cable only for its designated use. Never use it to carry the device or to pull the plug out of a socket. Keep the cable away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cables increase the risk of electric shock.
- e) If using the device in a damp environment cannot be avoided, a residual current device (RCD) should be applied. The use of an RCD reduces the risk of electric shock.
- f) Do not use the device if the power cord is damaged or shows obvious signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service centre.
- g) To avoid electric shock, do not immerse the cord, plug or device in water or other liquids. Do not use the device on wet surfaces.

3.2 Safety in the workplace

- a) Make sure the workplace is clean and well lit. A messy or poorly lit workplace may lead to accidents. Try to think ahead, observe what is going on and use common sense when working with the device.
- b) Do not use the device in a potentially explosive environment, for example in the presence of flammable liquids, gases or dust. The device generates sparks which may ignite dust or fumes.
- c) If you discover damage or irregular operation, immediately switch the device off and report it to a supervisor without delay.
- d) If you are unsure about whether the product is operating correctly or if you find damage, please contact the manufacturer's service centre.
- e) Only the manufacturer's service centre may make repairs to the product. Do not attempt to make repairs yourself!
- f) In case of fire, use a powder or carbon dioxide (CO₂) fire extinguisher (one intended for use on live electrical devices) to put it out.
- g) Children or unauthorised persons are forbidden to enter a work station. A distraction may result in loss of control over the device.
- h) Use the device in a well-ventilated space.
- i) Regularly inspect the condition of the safety labels. If the labels are illegible, they must be replaced.

- j) Please keep this manual available for future reference. If this device is passed on to a third party, the manual must be passed on with it.
- k) Keep packaging elements and small assembly parts in a place not available to children.
- l) Keep the device away from children and animals.
- m) If this device is used together with another equipment, the remaining instructions for use shall also be followed.

REMEMBER! WHEN USING THE DEVICE, PROTECT CHILDREN AND OTHER BYSTANDERS.

3.3 Personal safety

- a) Do not use the device when tired, ill or under the influence of alcohol, narcotics or medication which can significantly impair the ability to operate the device.
- b) The machine may be operated by physically fit persons who are able to handle the machine, are properly trained, who have reviewed this operating manual and have received training in occupational health and safety.
- c) The machine is not designed to be handled by persons (including children) with limited mental and sensory functions or persons lacking relevant experience and/or knowledge unless they are supervised by a person responsible for their safety or they have received instruction on how to operate the machine.
- d) When working with the device, use common sense and stay alert. Temporary loss of concentration while using the device may lead to serious injuries.
- e) Use personal protective equipment as required for working with the device, specified in section 1 "Legend". The use of correct and approved personal protective equipment reduces the risk of injury.
- f) To prevent the device from accidentally switching on, make sure the switch is on the OFF position before connecting to a power source.
- g) Do not wear loose clothing or jewellery. Keep hair, clothes and gloves away from moving parts. Loose clothing, jewellery or long hair may get caught in moving parts.
- h) Do not put your hands or other items inside the device while it is in use!

3.4 Safe device use

- a) Do not use the device if the "ON/OFF" switch does not function properly (does not switch the device on and off). Devices which cannot be switched on and off using the "ON/OFF" switch are hazardous, should not be operated and must be repaired.
- b) Disconnect the device from the power supply before commencement of adjustment, cleaning and maintenance. Such a preventive measure reduces the risk of accidental activation.
- c) Keep the device out of the reach of children.

- d) Device repair or maintenance should be carried out by qualified persons, only using original spare parts. This will ensure safe use.
- e) To ensure the operational integrity of the device, do not remove factory-fitted guards and do not loosen any screws.
- f) Do not touch articulated parts or accessories unless the device has been disconnected from the power source.
- g) Do not move, adjust or rotate the device in the course of work.
- h) Do not leave this appliance unattended while it is in use.
- i) Clean the device regularly to prevent stubborn grime from accumulating.
- j) Do not run the device when empty.
- k) It is forbidden to interfere with the structure of the device in order to change its parameters or construction.
- l) Keep the device away from sources of fire and heat.
- m) Do not overload the device.
- n) NOTE! When the heating function is being used, the workplate warms up to a high temperature. Do not touch the hot parts of the device as they may cause burns.
- o) Heated liquid can cause burns. It is advisable to use gloves when handling a vessel filled with hot liquid. Take special care not to spill the liquid.
- p) When working with flammable materials, pay special attention that the set heating temperature is not too high. The set heating temperature must be at least 25°C lower than the ignition temperature of the stirred substance.
- q) Make sure that the container can withstand the set heating temperature. A melted container may cause a fire. Never heat an empty container. Control the amount of liquid in containers. The liquid may evaporate.
- r) Spilling the flammable contents of the container onto the heating elements of the device may cause a fire.
- s) Special attention must be paid to the hazards arising from:
 - Flammable materials and substances with low boiling temperature
 - Overfilled containers
 - Containers made of brittle material or of low melting temperature
- t) Pathogenic materials must be stirred only in closed containers.
- u) Do not cover the device while it is in operation.
- v) Make sure that the power cord does not touch any heated parts of the device.
- w) Keep the device away from strong magnetic fields.
- x) Do not use metal vessels or vessels with a thick bottom. This will make work with the device more difficult.

- y) When mixing viscous liquids reduce the rotational speed and mix a small volume of liquid only.
- z) Follow biosafety regulations and check the resistance of the vessel when working with chemical, corrosive, infectious, pathogenic and biological materials.


ATTENTION! DESPITE THE SAFE DESIGN OF THE DEVICE AND ITS PROTECTIVE FEATURES, AND DESPITE THE USE OF ADDITIONAL ELEMENTS PROTECTING THE OPERATOR, THERE IS STILL A SLIGHT RISK OF ACCIDENT OR INJURY WHEN USING THE DEVICE. STAY ALERT AND USE COMMON SENSE WHEN USING THE DEVICE.

4. Use guidelines

The device is designed for mixing and heating solutions in laboratory settings, making it ideal for chemical, biological, and pharmaceutical applications. It combines a magnetic stirrer with a hotplate, allowing users to heat and stir liquids simultaneously, ensuring even mixing and temperature distribution. The device is particularly useful for preparing chemical reactions, dissolving solids, and maintaining consistent temperatures in experiments.

The user is liable for any damage resulting from unintended use of the device.

4.1 Product overview

1) Bracket

2) Cross clamp

- 3) Temperature probe
- 4) Heating plate
- 5) Operation panel

4.2 Preparing for use

4.2.1 Appliance location

- Indoor use only
- Altitude: $\leq 2000\text{m}$
- The working temperature range of the instrument is $+5^{\circ}\text{C}\sim+40^{\circ}\text{C}$
- The applicable temperature range of the instrument is $<80\%$
- There is no vibration or airflow that affects performance around

4.3 Device use

4.3.1 Safety precautions

- It is strictly forbidden to unplug and plug the power connector and toggle the power button when the hands are wet with liquid!
- It is strictly forbidden to unplug the power plug when the instrument is live!
- It is strictly forbidden to maintain and clean the instrument when it is powered on!
- It is strictly forbidden to install the instrument on an uneven, swaying and vibrating work surface!

4.3.2 Operation instructions

- 1) Turn on the power switch.
- 2) Press the **Menu** key to enter the temperature setting mode. Use the Up / Down keys to adjust the temperature. If the temperature is set to "--:--", heating is disabled.
- 3) Press the **Menu** key again to select the stirring speed setting mode. Use the Up / Down keys to adjust the speed. If the speed is set to "00:00", stirring is disabled.
- 4) Press the **Menu** key again to enter the timing setting mode. Use the Up / Down keys to adjust the time. If the time is set to "--:--", the device operates in continuous mode (long-term action).

Note: In continuous mode, the time area on the display will be hidden. If a specific time is set, the device operates in timing mode (the countdown time will be displayed).

- 5) Press the **Start / Stop** key to begin or stop the operation.

4.4 Cleaning and maintenance

4.4.1 General instructions

- a) Unplug the mains plug and allow the device to cool completely before each cleaning, adjustment or replacement of accessories, or if the device is not being used.
 - Wait for the rotating elements to stop.
- b) Always unplug the device before cleaning or putting it away.
- c) Use only non-corrosive cleaners to clean the surface.
- d) After cleaning the device, all parts should be dried completely before using it again.
- e) Store the unit in a dry, cool place, free from moisture and direct exposure to sunlight.
- f) Do not spray the device with a water jet or submerge it in water.
- g) Do not allow water to get inside the device through vents in the housing of the device.
- h) The device must be regularly inspected to check its technical efficiency and spot any damage.
- i) Use a soft, damp cloth for cleaning.
- j) Do not use sharp and/or metal objects for cleaning (e.g., a wire brush or a metal spatula) because they may damage the surface material of the appliance.
- k) Do not clean the device with an acidic substance, agents of medical purposes, thinners, fuel, oils or other chemical substances because it may damage the device.

4.4.2 Disposing of used devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

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