## USER MANUAL

Bedienungsanleitung  
Instrukcja obsługi  
Návod k použití  
Manuel d'utilisation  
Istruzioni per l'uso  
Manual de instrucciones  
Használati útmutató  
Bruksanvisning  
Käyttöohje  
Gebruiksaanwijzing  
Bruksanvisning  
Instruções de utilização  
Používateľská príručka  
Ръководство за употреба  
Οδηγίες χρήσης  
Upute za uporabu  
Naudojimo instrukcija  
Manual de utilizare  
Navodila za uporabo

---

---

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<b>DE</b>	Produktname	Magnetrührer
<b>EN</b>	Product name	MINI MAGNETIC STIRRER
<b>PL</b>	Nazwa produktu	Małe mieszadło magnetyczne
<b>CZ</b>	Název výrobku	Mini magnetické míchačky
<b>FR</b>	Nom du produit	Mini agitateur magnétique
<b>IT</b>	Nome del prodotto	Mini agitatore magnetico
<b>ES</b>	Nombre del producto	Mini agitador magnético
<b>HU</b>	Termék neve	Mini mágneses keverő
<b>DA</b>	Produktnavn	Mini magnetisk omrører
<b>FI</b>	Tuotteen nimi	Mini magneettinen sekoitin
<b>NL</b>	Productnaam	Mini magnetische roerder
<b>NO</b>	Produktnavn	Mini magnetisk omrører
<b>SE</b>	Produktnamn	Mini magnetomrörare
<b>PT</b>	Nome do produto	Mini agitador magnético
<b>SK</b>	Názov produktu	Mini magnetická miešačka
<b>BG</b>	Име на продукта	Мини магнитна бъркалка
<b>EL</b>	Όνομα προϊόντος	Μίνι μαγνητικός αναδευτήρας
<b>HR</b>	Naziv proizvoda	Mini magnetno mešalo
<b>LT</b>	Produktu pavadinimas	Mini magnetinė maišyklė
<b>RO</b>	Numele produsului	Mini agitator magnetic
<b>SL</b>	Ime izdelka	Mini magnetni mešalnik
<b>DE</b> Modell   <b>EN</b> Product model   <b>PL</b> Model produktu   <b>CZ</b> Model výrobku   <b>FR</b> Modèle   <b>IT</b> Modello   <b>ES</b> Modelo   <b>HU</b> Modell   <b>DA</b> Model   <b>FI</b> Tuotteen malli   <b>NL</b> Productmodel   <b>NO</b> Produktmodell   <b>SE</b> Produktmodell   <b>PT</b> Modelo do produto   <b>SK</b> Model   <b>BG</b> Модел на продукт   <b>EL</b> Μοντέλο προϊόντος   <b>HR</b> Model proizvoda   <b>LT</b> : Gaminio modelis   <b>RO</b> : Model de produs   <b>SL</b> : Model izdelka		SBS-LAB-115
<b>DE</b> Hersteller   <b>EN</b> Manufacturer   <b>PL</b> Producent   <b>CZ</b> Výrobce   <b>FR</b> Fabricant   <b>IT</b> Produttore   <b>ES</b> Fabricante   <b>HU</b> Termelő   <b>DA</b> Producent   <b>FI</b> Valmistaja   <b>NL</b> Producent   <b>NO</b> Produsent   <b>SE</b> Tillverkare   <b>PT</b> Fabricante   <b>SK</b> Výrobca   <b>BG</b> Производител   <b>EL</b> Κατασκευαστής   <b>HR</b> Proizvođač   <b>LT</b> Gamintojas   <b>RO</b> Producător   <b>SL</b> Proizvajalec		expondo Polska sp. z o.o. sp. k.
<b>DE</b> Anschrift des Herstellers   <b>EN</b> Manufacturer Address   <b>PL</b> Adres producenta   <b>CZ</b> Adresa výrobce   <b>FR</b> Adresse du fabricant   <b>IT</b> Indirizzo del produttore   <b>ES</b> Dirección del fabricante   <b>HU</b> A gyártó címe   <b>DA</b> Producentens adresse   <b>FI</b> Valmistajan osoite   <b>NL</b> Adres producent   <b>NO</b> Produsentens adresse   <b>SE</b> Tillverkarens adress   <b>PT</b> Endereço do fabricante   <b>SK</b> Adresa výrobcu   <b>BG</b> Адрес на производителя   <b>EL</b> : Διεύθυνση κατασκευαστή   <b>HR</b> Adresa proizvođača   <b>LT</b> Gamintojo adresas   <b>RO</b> Adresa producătorului   <b>SL</b> Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra   Poland, EU

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This User Manual has been translated for your convenience using machine translation. Reasonable efforts have been made to provide an accurate translation; however, no automated translation is perfect nor is it intended to replace human translators. The official User Manual is the English version. Any discrepancies or differences created in the translation are not binding and have no legal effect for compliance or enforcement purposes. If any questions arise related to the accuracy of the information contained in the User Manual, please refer to the English version of those contents which is the official version.

## Technical data

Parameter description	Parameter value
Product name	Mini Magnetic Stirrer
Model	SBS-LAB-115
Voltage/frequency	230 V~/50 Hz
Power	5 W
Max rotation speed	350-1800 rpm
Material	ABS plastic, PET
Stir bar size	20-30 mm
Protection class	IP 42
Dimensions (LxWxH)	16.8 x 16 x 4.7 cm
Platform dimensions	Ø 135 mm
Volume of agitated substance	2000 cm <sup>3</sup>
Weight	0.55 kg

# 1. Introduction

The magnetic stirrer is designed with advanced international standards and manufacturing techniques. It features a smooth appearance, a compact and durable structure, and is highly versatile. The device is primarily intended for stirring low-viscosity liquids or solid-liquid mixtures and is straightforward to operate.

## 2. Operating the device

### Application

The device is suitable for use in educational institutions, environmental protection, sanitation, epidemic prevention, the chemical industry, water treatment, and medical applications.

### Step-by-Step Operation

1. Power On: Plug the device into a power source and turn the control knob to the right to switch on the power.
2. Adjust Speed: Rotate the control knob to set the desired speed.
3. Power Off: Turn the control knob to the left to switch off the power.

### Operating environment

**Indoor Use:** The instrument is designed for indoor use only.

**Altitude:** Operates at altitudes up to 2000 meters above sea level.

**Temperature Range:** The working temperature range is +5°C to +40°C.

**Humidity:** Suitable for environments with relative humidity up to 80%.

**Vibration and Airflow:** Ensure there is no vibration or airflow that could affect performance.


**Air Quality:** The surrounding air should be free from conductive dust, explosive gases, or corrosive gases.

**Safety Distance:** Maintain a safety distance of at least 10 cm in front of and behind the instrument during operation.

## 3. Safety precautions

Do not unplug or plug in the power connector or toggle the power button with wet hands. Never unplug the power cord while the device is still powered on. Do not perform maintenance or cleaning while the device is powered on. Ensure the device is installed on a stable, level surface free from swaying or vibration.

## 4. Product overview

1. Stirring platform
2. Stirring speed control knob
3. Support feet
4. DC socket

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