USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

ORBITAL SHAKER

DE	Produktname	Rollenmischer
EN	Product name	Orbital Shaker
PL	Nazwa produktu	Wyrzasarka orbitalna
CZ	Nazev vyrobku	Orbitaln třepaka
FR	Nom du produit	Agitateur orbital
IT	Nome del prodotto	Agitatore orbitale
ES	Nombre del producto	Agitador orbital
HU	Termek neve	Orbitalis razogep
DA	Produktnavn	Orbital ryster
FI	Tuotteen nimi	Orbiittisekoitin
NL	Productnaam	Orbitale schudder
NO	Produktnavn	Orbital rister
SE	Produktnamn	Orbitalskakare
PT	Nome do produto	Agitador orbital
SK	Nazov produktu	Orbitalna trepaka
BG	Име на продуктa	Орбитален шейкr
EL	Όνομα προϊόντος	Ταλαντωτης τροχας
HR	Naziv proizvoda	Orbitalni tresilnik
LT	Produktas pavadinimas	Orbitin maiykl
RO	Numele produsului	Agitator orbital
SL	Ime izdelka	Orbitalni stresalnik
DE Modell EN Product model PL Model produktu CZ Model vyrobku FR Modle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LAB-125
DE Hersteller EN Manufacturer PL Producent CZ Vyrobce FR Fabricant IT Produttore ES Fabricante HU Termel DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Vyrobca BG Производител EL Κατασκευαστης HR Proizvoa LT Gamintojas RO Productor SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa vyrobce FR Adresse du fabricant IT Indirizzo del produttore ES Direccin del fabricante HU A gyrt cme DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereo do fabricante SK Adresa vyrobci BG Адрес на производителя EL : Διευθυνση κατασκευαστη HR Adresa proizvoaa LT Gamintojo adresas RO Adresa productorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Gra Poland, EU

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This User Manual has been translated for your convenience using machine translation. Reasonable efforts have been made to provide an accurate translation; however, no automated translation is perfect nor is it intended to replace human translators. The official User Manual is the English version. Any discrepancies or differences created in the translation are not binding and have no legal effect for compliance or enforcement purposes. If any questions arise related to the accuracy of the information contained in the User Manual, please refer to the English version of those contents which is the official version.

FOR YOUR SAFETY AND CORRECT OPERATION, PLEASE CAREFULLY READ THIS INSTRUCTION BEFORE USING IT.


NOTE: All of the information reported herein is based on data available at the moment of printing. The factory reserves the right to modify its own products at any moment without notice and incur any sanction. So, it is suggested to always verify possible updates.

1. TECHNICAL DATA

Parameter description	Parameter value
Product name	Orbital shaker
Model	SBS-LAB-125
Voltage/frequency	230 V / 50 Hz
Rated power [W]	24 / 10
Rotation speed range	20 – 200 rpm
Max load	2 kg
Roll diameter	20 mm
Timer	Yes
Dimensions [Width x Depth x Height; mm]	360 x 250 x 140
Weight [kg]	5.15

2. PRODUCT OVERVIEW

EN

1. Speed control knob - rotate the knob to set the speed

2. Start/Stop button - press the knob to start or stop the program

General product information

This shaker is mainly used in research applications such as bacterial culture, fermentation, hybridization and biochemical reactions, enzymes, cell tissue research, etc., which place high demands on the oscillation frequency.

It is an indispensable laboratory equipment for the precision cultivation and preparation of plants, biology, microorganisms, genetics, viruses, environmental protection, medicine and other scientific research, teaching and production departments.

3. USAGE

Widely used in the experiments of various universities, environmental protection, sanitation, epidemic prevention, chemical industry, tap water, medical treatment and other units; it is the scientific research of plants, biology, microorganisms, genetics, viruses, environmental protection, medicine, education and production departments for precision training and preparation. Essential laboratory equipment.

4. INSTRUMENT USING ENVIRONMENT

- Indoor use
- Altitude: $\leq 2000\text{m}$
- The working temperature range of the instrument is $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$
- The applicable temperature range of the instrument is $\leq 80\%$
- There is no vibration and airflow that affect performance around
- There is no conductive dust, explosive gas or corrosive gas in the surrounding air
- Keep a safety distance of 10cm before and after the instrument is working

5. STEPS FOR USAGE

- Turn on
- Left knob sets the time, Right knob sets the speed
- Press the right knob to start/Stop.

6. Safety

It is strictly forbidden to unplug and plug the power connector and toggle the power button when the hands are wet with liquid!

It is strictly forbidden to unplug the power plug when the instrument is live!

It is strictly forbidden to maintain and clean the instrument when it is powered on!

It is strictly forbidden to install the instrument on an uneven, swaying and vibrating work surface!

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