USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso Manual
de instrucciones Használati
útmutató Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning Instruções
de utilização Používateľská
príručka Ръководство за
употреба Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

VACUUM FILTRATION DEVICE

DE	Produktname	Vakuum-filtrationsgerät
EN	Product name	Vacuum filtration device
PL	Nazwa produktu	Urządzenie do filtracji próżniowej
CZ	Název výrobku	Vakuové filtrační zařízení
FR	Nom du produit	Dispositif de filtration sous vide
IT	Nome del prodotto	Dispositivo di filtrazione sotto vuoto
ES	Nombre del producto	Dispositivo de filtración al vacío
HU	Termék neve	Vákuum szűrőberendezés
DA	Produktnavn	Vakuumfiltrationsanordning
FI	Tuotteen nimi	Tyhjiösuodatuslaite
NL	Productnaam	Vacuümfiltratie-apparaat
NO	Produktnavn	Vakuumfiltreringsanordning
SE	Produktnamn	Vakuumfiltreringsanordning
PT	Nome do produto	Dispositivo de filtração a vácuo
SK	Názov produktu	Vakuové filtračné zariadenie
BG	Име на продукта	Вакуумно филтрационно устройство
EL	Όνομα προϊόντος	Συσκευή φιλτραρίσματος κενού
HR	Naziv proizvoda	Uređaj za vakuum filtraciju
LT	Produkto pavadinimas	Vakuuminio filtravimo įrenginys
RO	Numele produsului	Dispozitiv de filtrare în vid
SL	Ime izdelka	Naprava za vakuumsko filtracijo
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-FD-130
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

DE

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info@expondo.com

www.expondo.de

This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Safety and operation instructions

Product name: Vacuum Filtration Device

Model name: SBS-FD-130

DO NOT USE THE DEVICE UNLESS YOU HAVE THOROUGHLY READ AND UNDERSTOOD THIS USER MANUAL.

Safety Instructions

Warning: Read all safety warnings and instructions. Failure to follow the warnings and instructions may result in serious injury.

1. General Safety

- This product is intended for intended use only.
- Do not use the product for purposes other than those described in this manual.
- Keep this manual for future reference.
- A Vacuum Filtration Device is designed for laboratory use to separate solids from liquids through the application of a vacuum. It is commonly used in chemical, biological, and environmental laboratories for tasks such as sample preparation, purification, and analysis.

2. Operational Safety

- Use the product only on a stable, level surface.
- Do not operate the product if it is damaged or malfunctioning.
- Keep away from children and pets.

3. Maintenance Safety

- Only qualified personnel should perform repairs.
- Use only original spare parts.

4. Warning Symbols

- **Caution:** Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

EN

- **Warning:** Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
 - **Potential risks and hazards:** Lists all risks and hazards directly or indirectly connected with the use or misuse of the product:
 - **Mechanical Hazards:** Improper handling or assembly of the device can lead to breakage of glass components, which may cause cuts or other injuries. Always handle glassware with care and follow the manufacturer's assembly instructions.
 - **Vacuum Hazards:** The use of vacuum pressure can pose risks such as implosion of glassware or sudden release of vacuum, which can cause injury. Ensure that all components are rated for vacuum use and inspect them regularly for damage.
 - **Chemical Hazards:** The device is often used with hazardous chemicals that can pose risks such as burns, poisoning, or respiratory issues. Use appropriate personal protective equipment (PPE) and work in a well-ventilated area or fume hood.
 - **Electrical Hazards:** If the vacuum pump is electrically powered, there is a risk of electrical shock or fire. Ensure that the pump is properly grounded and regularly maintained.
 - **Contamination Hazards:** Improper cleaning and maintenance of the device can lead to cross-contamination of samples, affecting the accuracy of results. Follow proper cleaning protocols and regularly maintain the device.
-

Maintenance and Cleaning

1. Regular Maintenance

- Inspect the product regularly for signs of wear or damage.
- Replace any damaged parts immediately.

2. Cleaning

- Use a soft, damp cloth to clean the exterior.
- Do not use abrasive cleaners or solvents.

3. Storage

- Store in a cool, dry place.
 - Protect from dust and moisture.
-

Disposal

Environmental Protection

- Do not dispose of this product as unsorted municipal waste.

EN

- By ensuring proper disposal, you help prevent potential negative consequences for the environment and human health.
-

Declaration of Compliance

We, expondo Polska sp. z o.o. sp. k., declare under our sole responsibility that the Vacuum Filtration Device, Model SBS-FD-130, to which this declaration relates, is in conformity with the relevant Union harmonization legislation:

- Regulation (EU) 2023/988 (General Product Safety Regulation)

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