USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

LAB WATER BATH

DE	Produktname	Laborwasserbad
EN	Product name	Laboratory water bath
PL	Nazwa produktu	łaźnia wodna laboratoryjna
CZ	Název výrobku	Laboratorní vodní lázeň
FR	Nom du produit	Bain-marie de laboratoire
IT	Nome del prodotto	Bagno d'acqua da laboratorio
ES	Nombre del producto	Baño de agua de laboratorio
HU	Termék neve	Laboratóriumi vízfürdő
DA	Produktnavn	Laboratorievandbad
FI	Tuotteen nimi	Laboratoriovesikylpy
NL	Productnaam	Laboratoriumwaterbad
NO	Produktnavn	Laboratorievannbad
SE	Produktnamn	Laboratorievattenbad
PT	Nome do produto	Banho-maria de laboratório
SK	Názov produktu	Laboratórny vodný kúpeľ
BG	Име на продукта	Лабораторна водна баня
EL	Όνομα προϊόντος	Εργαστηριακό λουτρό νερού
HR	Naziv proizvoda	Laboratorijska vodena kupelj
LT	Produkto pavadinimas	Laboratorinis vandens vonia
RO	Numele produsului	Baie de apă de laborator
SL	Ime izdelka	Laboratorijska vodna kopel
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-WB-01 SBS-WB-02
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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
This User Manual has been translated for your convenience using machine translation. Reasonable efforts have been made to provide an accurate translation; however, no automated translation is perfect nor is it intended to replace human translators. The official User Manual is the English version. Any discrepancies or differences created in the translation are not binding and have no legal effect for compliance or enforcement purposes. If any questions arise related to the accuracy of the information contained in the User Manual, please refer to the English version of those contents which is the official version.

Technical data

Parameter description	Parameter value	
Product name	Lab Water bath	
Model	SBS-WB-01	SBS-WB-02
Rated voltage [V~] / frequency [Hz]	230/50	230/50
Rated power [W]	800	1000
Temperature range [°C]	-5 ~ 100	-20 ~ 100
Tank volume [l]	6	6
Dimensions [width x depth x height; mm]	375 x 320 x 665	375 x 320 x 665
Weight [kg]	24.5	24.5


- 1 – ON/OFF button
- 2 – Setting button
- 3 – Control box
- 4 – The shift button
- 5 – Button for increasing the input values
- 6 – Button for decreasing the input values
- 7 – Button for switching on the circulation
- 8 – Button for switching on the refrigeration

- 1 – ON/OFF switch
- 2 – Fuse
- 3 – Power cord
- 4 – Water inlet
- 5 – Water outlet

Operating Steps:

1. Connect the circulating pump before startup (if there is no water inlet and outlet, this step can be omitted):
 - 1.1. For the connection of circulating pump during internal circulation, connect the water outlet and water inlet directly with both ends of the hose supplied with the goods.
 - 1.2. For the connection of circulating pump during external circulation, connect the water outlet to the inlet of container / equipment outside the tank with hose, and connect the water inlet to the outlet of container / equipment outside the tank.
2. Add the corresponding liquid medium into the tank. The liquid level of the liquid medium shall not be 20 mm lower than the worktable to avoid damage caused by exposed dry burning of the heater.
 - 2.1. When the operating temperature is below 8°C, industrial alcohol or antifreeze is generally used for liquid media.
 - 2.2. At an operating temperature of 8°C-75°C, pure water is generally selected for liquid media.
 - 2.3. At an operating temperature of 75°C-100°C, dimethyl silicone oil with a viscosity of 5cs is generally selected for liquid media.
 - 2.4. At an operating temperature of 101-200°C, dimethyl silicone oil with a viscosity of 10cs is generally selected for liquid media.
 - 2.5. When the operating temperature is above 200 °C, dimethyl silicone oil with a viscosity of 50cs is generally selected for liquid media.
3. Plug in the power, turn on the "main power" switch, and then turn on the "power" button on the operation panel. (Note: "main power" switch is on the back of the instrument).
4. Perform the instrument operation as follows:
 - 4.1. Display description:

- | | |
|---|---|
| 1. Main interface | 2. Settings |
| 3. Self-diagnose | 4. Time unit |
| 5. Time | 6. Display |
| 7. Settings | 8. Heating enabled |
| 9. Cooling enabled | 10. Working status |
| 11. Stop | 12. Circulation |
| 13. Heating | 14. Cooling |
| 15. Alarm | 16. Settings button |
| 17. Shift button | 18. Button for decreasing the input values |
| 19. Button for increasing the input values | 20. ON/OFF button |
| 21. Button for switching on the refrigeration | 22. Button for switching on the circulation |

4.2. Description of the buttons:

- ◀: shift key
- ▼: decrease key
- ▲: increase key

4.3. Setting of temperature parameters:

Press the setting function button "set" to enter the setting state of the temperature setting value. At this time, the last value of the setting temperature "SV" on the display screen flashes.

First press the shift button, then press the increase or decrease buttons to set the working temperature, and then press the setting function button "set" once. At this time, the last value of the timed shutdown time (min) "0000" on the upper right corner of the display screen flashes.

If it is not necessary to keep the default "0000" regularly, press the setting function key "set" again to exit the parameter settings and save the parameter settings values. At this time, the display "PV" shows the temperature of the liquid medium in the current tank.

5. After the temperature parameter setting is completed, first open the "circulation" button on the operation panel, and then open the "refrigeration" button (Note: when the working temperature exceeds 40°C, the compressor automatically enters the protection, and the refrigeration does not start), and then the microcomputer enters the automatic control state.
6. Description of correction between quantity temperature and actual temperature:

In the normal display state, long press the setting function key "set" for 3 seconds, the password prompt "LC" is displayed in the upper row of the display window, and the password value is displayed in the lower row. At this time, press the increase button "▲" to enter the password "00000 3", then click the setting function button "set" to enter the internal parameter setting state, and then click the setting function button "set" to select the parameters to be modified. Select the indicator in the upper window to display the word "PB1". At this time, it is the temperature correction interface. If the actual temperature is higher than the displayed temperature, press the increase button "▲" to increase, press the decrease button "▼" to decrease, and then long press the setting function button "set" for 3 seconds. The parameter values can be automatically saved and exit this state.

Precautions for use:

1. Liquid medium shall be added into the tank before use, and air opening is strictly prohibited. (Note: the liquid shall not contain acid and alkali).
2. The power supply power shall be greater than the total power of the instrument, and the power supply must have a good "grounding" device. The voltage used shall be subject to the technical parameter table.
3. The instrument shall be placed in a ventilated and dry place, and the back and both sides shall be 300mm away from the obstacles to maintain good heat dissipation.

After use, set all switches to off state, unplug the power plug, open the drain valve, discharge the medium in the tank, and keep the tank dry and tidy.

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