USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

PRECISION BALANCE

DE	Produktname	Präzisionswaage
EN	Product name	Precision balance
PL	Nazwa produktu	Precyzyjna waga
CZ	Název výrobku	Přesné vyvážení
FR	Nom du produit	Balance de précision
IT	Nome del prodotto	Bilancia di precisione
ES	Nombre del producto	Balanza de precisión
HU	Termék neve	Precíziós mérleg
DA	Produktnavn	Præcisionsvægt
FI	Tuotteen nimi	Tarkka tasapaino
NL	Productnaam	Precisiebalans
NO	Produktnavn	Presisjonsbalanse
SE	Produktnamn	Precisionsbalans
PT	Nome do produto	Balança de precisão
SK	Názov produktu	Presné vyváženie
BG	Име на продукта	Прецизен баланс
EL	Όνομα προϊόντος	Ισορροπία ακριβείας
HR	Naziv proizvoda	Precizna vaga
LT	Produkto pavadinimas	Tikslus balansas
RO	Numele produsului	Balanta de precizie
SL	Ime izdelka	Natančna tehtnica
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LW-1010 SBS-LW-1020
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value	
Product name	Precision balance	
Model	SBS-LW-1010	SBS-LW-1020
Max load [g]	60	120
Min load [g]	0.1	
Division [g]	0.005	
Battery	6 x AAA Batteries	
Power supply [VDC / mA]	9 / 300	
Dimensions (Width x Length x Height) (mm)	250x170x130	520x350x300
Weight [kg]	1.08	10
Protection class	II	

Purpose

The product is used to measure the weight of objects with high accuracy and precision, typically within a milligram or even finer increments. Its purpose is essential in environments where precise measurements are crucial, such as in laboratories, research facilities, pharmaceuticals, and industrial applications.

Usage

Before weighing

Please always turn on the scale and allow it to warm up for 4-5 minutes before weighing to ensure proper functionality and accurate measurements.

NOTE

- Calibration may be required before weighing. Read "Calibration" first and if necessary, calibrate your scale for accurate weighing.
- Error messages

Err-0: Overload

: Low battery

Weighing procedures

1. Press [ON/OFF] to turn on the scale.

All display segments will briefly appear on the LCD. After a few seconds, "0.000" will be displayed, indicating that the scale is ready.

2. Select the weighing unit with [UNIT].

Press **[UNIT]** to select a weighing unit, eg. "g", "ct", "oz", "dwt", "ozt", "gn" and "mg". Once the unit has been selected, the selected unit will be displayed to the weight value.

3. Start weighing

- **Without a container:**
 - 1) Ensure the reading is "0". If not, press **[ZERO]** to reset the display to "0".
 - 2) Place the object to be weighed on the platform.
 - 3) Once the reading becomes stable, the final weight will be displayed.
- **With a container:**
 - 1) Place an empty container on the platform.
 - 2) Wait for the stable reading to appear and then press **[ZERO]** to reset the display to "0".
 - 3) Once the reading stabilizes, the accurate weight will be displayed.

Calibration

NOTE Calibration may be required when the scale is initially installed or before weighing.

1. Turn on the scale by pressing **[ON/OFF]** and allow it to warm up for 2-3 minutes.
2. Enter calibration mode
 - 1) Turn off the scale by pressing **[ON/OFF]**.
 - 2) Press and hold **[UNIT]**, then press **[ON/OFF]** to turn the scale back on. The display will show "CAL" followed by the A/D value.
3. Calibrate using a weight
 - 1) Wait for the stable indicator to appear on the display.
 - 2) Press **[UNIT]**, and the display will show the flashing calibration weight.
 - 3) Place the specified calibration weight on the platform.
 - 4) Once the stable indicator is displayed again, press **[UNIT]**. The display will show "-----" followed by the A/D value, indicating that calibration is complete.
 - 5) Turn off the scale by pressing **[ON/OFF]**, then turn it back on and place the calibration weight on the platform to verify accuracy.
 - 6) If the weight is not accurate, repeat steps 1-2.
4. Return to weighing mode
 - 1) Turn off the scale by pressing **[ON/OFF]**.
 - 2) Turn the scale on again by pressing **[ON/OFF]** to return to normal weighing mode.

Function Settings

1. **Enter function setting mode.**
 - 1) Turn off the scale by pressing **[ON/OFF]**.

- 2) Press and hold the **[ZERO]** button, then press **[ON/OFF]** to turn the scale back on. The display will show either "A_ON" or "A_OFF."

2. Selection of auto shut off mode.

- When the display shows "A_ON" or "A_OFF", press **[ZERO]** to toggle between the two options:
 - **A_ON**: Auto shut-off function is activated.
 - **A_OFF**: Auto shut-off function is deactivated.

3. Selection of auto backlight mode.

- After setting the auto shut-off mode, press the **[UNIT]** button. The display will show either "L_ON" or "L_OFF".
- Press **[ZERO]** to toggle between:
 - **L_ON**: Auto backlight function is activated.
 - **L_OFF**: Auto backlight function is deactivated.

4. Return to weighing mode.

- 1) Press **[ON/OFF]** to turn the scale off.
- 2) Press **[ON/OFF]** again to turn it back on and return to normal weighing mode.

Disposing of Used Devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

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