## USER MANUAL

Bedienungsanleitung  
Instrukcja obsługi  
Návod k použití  
Manuel d'utilisation  
Istruzioni per l'uso  
Manual de instrucciones  
Használati útmutató  
Brugsanvisning  
Käyttöohje  
Gebruiksaanwijzing  
Bruksanvisning  
Instruções de utilização  
Používateľská príručka  
Ръководство за употреба  
Οδηγίες χρήσης  
Upute za uporabu  
Naudojimo instrukcija  
Manual de utilizare  
Navodila za uporabo

---

---

# MAGNETIC STIRRER

|  |                      |   |
|--|----------------------|---|
| <b>DE</b>  | Produktname          | Magnetrührer  |
| <b>EN</b>  | Product name         | Magnetic Stirrer  |
| <b>PL</b>  | Nazwa produktu       | Mieszadło magnetyczne   |
| <b>CZ</b>  | Název výrobku        | Magnetická míchačka   |
| <b>FR</b>  | Nom du produit       | Agitateur magnétique  |
| <b>IT</b>  | Nome del prodotto    | Agitatore magnetico   |
| <b>ES</b>  | Nombre del producto  | Agitador magnético  |
| <b>HU</b>  | Termék neve          | Mágneses keverő   |
| <b>DA</b>  | Produktnavn          | Magnetomrører   |
| <b>FI</b>  | Tuotteen nimi        | Magneettisekoitin   |
| <b>NL</b>  | Productnaam          | Magnetische roerder   |
| <b>NO</b>  | Produktnavn          | Magnetrører   |
| <b>SE</b>  | Produktnamn          | Magnetomrörare  |
| <b>PT</b>  | Nome do produto      | Agitador magnético  |
| <b>SK</b>  | Názov produktu       | Magnetické miešadlo   |
| <b>BG</b>  | Име на продукта      | Магнитна бъркалка   |
| <b>EL</b>  | Όνομα προϊόντος      | Μαγνητικός αναδευτήρας  |
| <b>HR</b>  | Naziv proizvoda      | Magnetska miješalica  |
| <b>LT</b>  | Produktu pavadinimas | Magnetinė maišyklė  |
| <b>RO</b>  | Numele produsului    | Agitator magnetic   |
| <b>SL</b>  | Ime izdelka          | Magnetno mešalo   |
| <b>DE</b> Modell   <b>EN</b> Product model  <br><b>PL</b> Model produktu  <br><b>CZ</b> Model výrobku   <b>FR</b> Modèle  <br><b>IT</b> Modello   <b>ES</b> Modelo   <b>HU</b> Modell  <br><b>DA</b> Model   <b>FI</b> Tuotteen malli  <br><b>NL</b> Productmodel   <b>NO</b> Produktmodell  <br><b>SE</b> Produktmodell   <b>PT</b> Modelo do produto<br>  <b>SK</b> Model   <b>BG</b> Модел на продукт  <br><b>EL</b> Μοντέλο προϊόντος   <b>HR</b> Model<br>proizvoda   <b>LT</b> : Gaminio modelis   <b>RO</b> :<br>Model de produs   <b>SL</b> : Model izdelka  |                      | SBS-MR-1500-3/1<br>SBS-MR-1500-6/1<br>SBS-MR-1500-1/10              |
| <b>DE</b> Hersteller   <b>EN</b> Manufacturer  <br><b>PL</b> Producent   <b>CZ</b> Výrobce  <br><b>FR</b> Fabricant   <b>IT</b> Produttore  <br><b>ES</b> Fabricante   <b>HU</b> Termelő  <br><b>DA</b> Producent   <b>FI</b> Valmistaja  <br><b>NL</b> Producent   <b>NO</b> Produsent  <br><b>SE</b> Tillverkare   <b>PT</b> Fabricante  <br><b>SK</b> Výrobca   <b>BG</b> Производител  <br><b>EL</b> Κατασκευαστής   <b>HR</b> Proizvođač  <br><b>LT</b> Gamintojas   <b>RO</b> Producător  <br><b>SL</b> Proizvajalec   |                      | expondo Polska sp. z o.o. sp. k.                                    |
| <b>DE</b> Anschrift des Herstellers  <br><b>EN</b> Manufacturer Address   <b>PL</b> Adres<br>producenta   <b>CZ</b> Adresa výrobce  <br><b>FR</b> Adresse du fabricant   <b>IT</b> Indirizzo del<br>produttore   <b>ES</b> Dirección del fabricante  <br><b>HU</b> A gyártó címe   <b>DA</b> Producentens<br>adresse   <b>FI</b> Valmistajan osoite  <br><b>NL</b> Adres producent   <b>NO</b> Produsentens<br>adresse   <b>SE</b> Tillverkarens adress  <br><b>PT</b> Endereço do fabricante   <b>SK</b> Adresa<br>výrobcu   <b>BG</b> Адрес на производителя  <br><b>EL</b> : Διεύθυνση κατασκευαστή   <b>HR</b> Adresa<br>proizvođača   <b>LT</b> Gamintojo adresas  <br><b>RO</b> Adresa producătorului   <b>SL</b> Naslov<br>proizvajalca |                      | ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra   Poland, EU |

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference.

## Technical data

| Parameter description                   | Parameter value  |                 |                  |
|---|------------------|-----------------|------------------|
| Product name                            | Magnetic Stirrer |                 |                  |
| Model                                   | SBS-MR-1500-3/1  | SBS-MR-1500-6/1 | SBS-MR-1500-1/10 |
| Rated voltage [V~] / frequency [Hz]     | 100-240 / 50-60  |                 |                  |
| Rated power [W]                         | 45               | 90              | 30               |
| Display                                 | LCD              |                 |                  |
| Max stirring volume [L]                 | 1*3              | 1*6             | 10               |
| Speed range [rpm]                       | 100-1500         | 100-1500        | 50-1500          |
| Working plate size [mm]                 | N/A              | N/A             | 205*205          |
| Protection class                        | II               |                 |                  |
| Dimensions [width * depth * height; mm] | 540*200*95       | 545*275*93      | 370*223*95       |
| Weight [kg]                             | 3.55             | 4.45            | 2.8              |

## Description

### SBS-MR-1500-3/1 & SBS-MR-1500-6/1

This product features a compact design with high versatility and simple operation. It is an instrument for scientific research and laboratory experiments. Primarily designed for liquid mixing, this product is suitable for stirring or simultaneously heating and stirring low-viscosity liquids or solid-liquid mixtures. It is commonly used as an experimental tool in universities, environmental protection, sanitation, epidemic prevention, chemical tap water, medical treatment and other units.

### SBS-MR-1500-1/10

This product works with test tubes, centrifuge tubes, and funnels, and is commonly used in biochemistry, genetic engineering, and medical labs for various experiments.

**The user is liable for any damage resulting from unintended use of the device.**

## Operation


### SBS-MR-1500-3/1 & SBS-MR-1500-6/1

#### 1. Operation steps

- a) Plug in the power cord and turn on the power switch.
- b) Rotate the left knob to set the time and use the right knob to adjust the speed.
- c) Press the knob to confirm the selected value.
- d) Long press the right knob to start or stop the operation.



SBS-MR-1500-1/10

**1. Operation interface**

- 1- Display of the set speed value
- 2- Display of the status icon
- 3- Display of the running direction icon
- 4- Display of the current speed
- 5- Display of set time/run time
- 6- Display of the program-controlled sequence number

**2. Key functions**

|   |   |
|---|---|
| <p>Prog.</p>       | <p>"Prog" key<br/>Press the key to enter the program control mode.</p>  |
| <p>Switch</p>      | <p>"Switch" key<br/>Press this key to set the magneton clockwise and counterclockwise operation mode.</p>                           |
| <p>Start/Stop</p>  | <p>"Start/Stop" button<br/>Press this button to Start/Stop the operation of the instrument.</p>                                     |
|                    | <p>Time adjustment knob<br/>Short press to set the running time, counterclockwise value decreases, clockwise value increases.</p>   |
|                    | <p>Speed adjustment knob<br/>Short press to set the running speed, counterclockwise value decreases, clockwise value increases.</p> |

**3. Operation steps**

- a) Place the magnetic stir bar in the solution container and position the container on the magnetic stirrer plate.

- b) Ensure the magnetic stir bar is centred at the bottom of the container.
- c) Short press the left knob briefly to make the time "SET" flash, then rotate the left knob to adjust the timing. If set to "00:00," the instrument will operate continuously in the forward timing mode. The "SET" flash will stop automatically after three seconds of inactivity.
- d) Short press the right knob briefly to make the speed "SET" flash, then rotate the right knob to adjust the speed. During operation, short press the right knob to adjust the speed in real-time.
- e) Press the "Switch" button to change the running direction of the stir bar. You can also short press the "Switch" button during operation to reverse the direction in real-time.
- f) Short press the "Start/Stop" button to begin the operation according to the set parameters. Press it again to stop the operation.

#### 4. Program control function

- a) Short press "Prog" to enter the program control mode. Short press "Prog" again to move to the next program setting. Up to 6 program sequences, labelled "P1-P6," can be set, with "P1" displayed by default.
- b) Short press the left knob to make the time "set" flash, then rotate the left knob to adjust the time. If the time is set to 00:00, the program will run continuously.
- c) Short press the right knob to make the speed "set" flash, then rotate the right knob to adjust the desired speed.
- d) Short press the "Switch" button to change the running direction of the stirrer in program-controlled mode.
- e) After completing the settings for all 6 program sections, press the "Start/Stop" button to begin the operation according to the set program. Press the button again to stop the operation.

## Operation environment

SBS-MR-1500-3/1 & SBS-MR-1500-6/1

1. For indoor use only.
2. Maximum operating altitude:  $\leq 2000\text{m}$ .
3. The instrument is designed to operate within an ambient temperature range of  $+5^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$ .
4. The allowable humidity for the instrument's environment is  $\leq 80\%$ .
5. Ensure the surroundings have no vibration and airflow that could affect performance.
6. The air around the instrument must be free from conductive dust, explosive gases, or corrosive gases.
7. Maintain a safety clearance of at least 10 cm in front and behind the instrument during operation.
8. Suitable stir bar sizes range from 20mm to 40mm.

SBS-MR-1500-1/10

1. The workspace should be kept clean and dry.
2. Place the product on a flat, stable, flame-retardant surface with adequate friction.
3. Position the workbench in an area with minimal vibration.

4. Keep the instrument away from objects or equipment that produce or are affected by magnetic fields.
5. Do not use the product in explosive or hazardous environments.
6. Avoid prolonged use in areas with high humidity or excessive dust.

**CAUTION**

1. NEVER plug in or unplug the power connector or operate the power switch while holding liquid in both hands.
2. DO NOT insert or remove the power plug while the instrument is powered on.
3. Maintenance and cleaning of the instrument must NOT be performed while it is connected to a power source.
4. The instrument must NOT be installed on an uneven, unstable, or vibrating work surface.

## Maintenance

1. Do not immerse electronic devices in cleaning solutions.
2. Avoid using hard objects when cleaning the machine to prevent surface damage.
3. If the machine will not be used for an extended period, unplug the power cord, store it in a cool, dry place, and cover it with a dust bag.
4. Before restarting the machine after a long period of inactivity, check that all functions are operating normally.


## Troubleshooting

SBS-MR-1500-1/10



| Problem   | Possible cause   | Solution  |
|---|--|---|
| Unable to start normally.                         | 1. The power cable is not tightly connected.   | 1. Check whether the power cable is tightly connected or connected to another power socket.   |
| The product does not stir/stir when jumpers.      | 1. Motor fault.<br>2. The bottom of the beaker is uneven and out of position.<br>3. Wrong speed adjustment step.<br>4. The actual voltage is within the rated voltage range. | 1. Check the motor wiring or replace the motor.<br>2. Lay the bottom of the container flat, reposition the beaker, and start again.<br>3. Gradually adjust from low speed to high speed.<br>4. Adjust voltage to the rated voltage range. |
| There is no vortex in the liquid after operation. | 1. The rotational speed is improperly set.<br>2. The viscosity of the stirred mixture was too high.  | 1. Resetting the speed.<br>2. Reduce the viscosity of the mixture.  |

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