USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

HEATING MAGNETIC STIRRER

DE	Produktname	Magnetrührer mit Heizplatte
EN	Product name	Heating Magnetic Stirrer
PL	Nazwa produktu	Mieszadło magnetyczne
CZ	Název výrobku	Vyhřívací magnetické míchadlo
FR	Nom du produit	Agitateur magnétique chauffant
IT	Nome del prodotto	Agitatore magnetico riscaldante
ES	Nombre del producto	Agitador magnético de calentamiento
HU	Termék neve	Fűtés mágneses keverő
DA	Produktnavn	Opvarmning af magnetisk omrører
FI	Tuotteen nimi	Lämmitys magneettisekoitin
NL	Productnaam	Verwarmende magnetische roerder
NO	Produktnavn	Oppvarming magnetrører
SE	Produktnamn	Värmemagnetomrörare
PT	Nome do produto	Agitador magnético de aquecimento
SK	Názov produktu	Vyhrievacie magnetické miešadlo
BG	Име на продукта	Загриваща магнитна бъркалка
EL	Όνομα προϊόντος	Μαγνητικός αναδευτήρας θέρμανσης
HR	Naziv proizvoda	Magnetska miješalica za grijanje
LT	Produkto pavadinimas	Šildymo magnetinė maišyklė
RO	Numele produsului	Agitator magnetic de încălzire
SL	Ime izdelka	Grelno magnetno mešalo
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-MR-170 SBS-MR-180
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value	
Product name	Heating Magnetic Stirrer	
Model	SBS-MR-170	SBS-MR-180
Rated voltage [V~] / frequency [Hz]	220 / 50	220 / 50/60
Motor power [W]	10	13 (x4)
Heating power [W]	400	420 (x4)
Max capacity [ml]	5000	2000 (x4)
Max speed [rpm]	1500	50-2000 (without load)
Max temperature [°C]	160	Range: RT+5-160
Dimensions [width * length * height; mm]	170*315*490	713*350*478
Weight [kg]	3.25	13.3

Description


This product is designed for stirring low-viscosity liquids or solid-liquid mixtures. It is commonly used in chemical, biotechnological, and pharmaceutical industries.

The user is liable for any damage resulting from unintended use of the device.

Operation

Operation panel

SBS-MR-180

- **Display 1:** Shows various indicators, including measured temperature and rotational speed.
- **Heating Symbol (4):** Illuminates when there is heating output after temperature control is activated.

- **Rotational Speed Symbol (10):** Lights up when the display shows rotational speed.
- **Heating Ratio (9):** Indicates the current percentage of heating output.
- **Setting Symbol (SV):** Appears during temperature, speed, time, or parameter setting (query).
- **Measurement Symbol (PV):** Displays when the measured temperature and speed are shown.
- **Operation Symbol (8):** Constantly displayed when temperature control is active.
- **Alarm Symbol (3):** Activates if the heating plate temperature exceeds the AP value, shutting off heating, disconnecting the circuit, and sounding the buzzer. When the temperature exceeds the AL value, the temperature control flashes, and the buzzer sounds.
- **Self-Tuning Symbol (Palm Symbol 6):** Flashes during self-tuning.
- **Temperature Key:** Rotate to adjust the temperature setting (setting and temperature lights on); rotate left to decrease, and right to increase. If no rotation occurs within 0.5 seconds, exit setting mode. Tap to activate temperature control, tap again to turn off, and hold to enter parameter settings (setting symbol on).
- **Speed Key:** Rotate to adjust the speed setting (setting and speed lights on); rotate left to decrease, and right to increase. Exits automatically if no rotation within 0.5 seconds. Tap to activate speed control, and tap again to deactivate.
- **Buzzer:** Sounds during abnormalities, and also sounds 10 times when the timing function completes.

Operation instruction

SBS-MR-170

1. Turn on the power switch.
2. Rotate the heating knob; the red indicator light will illuminate, signalling heating output. Adjust the knob to the desired temperature, and the microcomputer will automatically heat to that set temperature.
3. Rotate the stirring knob; the green indicator light will turn on, indicating output from the stirring motor. The microcomputer will automatically accelerate to the set speed.
4. After switching off, be cautious of the heating plate, as it may still be hot.

SBS-MR-180

Turn on the switch, and the display will show the temperature sensor's reading, with the speed indicator set to OFF. Rotate the temperature knob to set the desired temperature, then press it to begin heating. Pressing it again will stop heating. Rotate the speed knob to set the target speed, then press it to gradually increase speed until it reaches the set level. Pressing it again will stop mixing.

If the upper half of the display shows "□□□□", it indicates the sensor is disconnected. If the lower half shows "□□□□", it signals a short circuit in the sensor, triggering the buzzer. Press any button to silence it. In either case, inspect the sensor or replace it if needed. If the temperature sensor for the control beaker disconnects, the controller automatically switches to the heating plate sensor for temperature control.

When the timer function is active, timing starts once the temperature stabilizes and will automatically stop when complete, accompanied by 10 beeps from the buzzer.

In the standard mode, long press the temperature button. The upper display will show "LK" and the lower display will show a value. Change this value to 0 to enter the user parameter settings interface, where values can be adjusted as needed.

Sign	Meaning	Description
rT	Time running value	OFF for no time setting, always running
T2	Heating plate temperature	Heater temperature display
ST	Time setting	0-9999 running time setting, 0 is always running


Please note that each adjustment to a function parameter may impact control performance. If no button is pressed within one minute, the device will automatically revert to standard mode, and certain function parameters may remain unchanged.


ATTENTION! When installing the temperature probes, ensure that each probe is matched to its corresponding instrument station number—avoid random placement. Quality inspection personnel rigorously test and calibrate each station's probe before shipment to ensure accuracy. Following the installation order improves measurement precision. Temperature probes for all other stations should also be installed according to this guideline.

Wiring diagram

SBS-MR-170

Speed board terminals

Heating board terminals

A- Power supply

B- Motor main circuit

C- YLQD-101-100W low voltage, 24V DC brushless motor driver

Speed board terminals:


- 24V+ and 24V-: connect to a 24V DC power supply
- U, V, W: connect to the main circuit of the brushless motor
- HU, HV, HW: connect to the hall sensors of the brushless motor
- +5V, GND: connect to the power supply for the brushless motor hall sensors

Note: The number of pulses per revolution for this brushless driver is equal to the number of motor pole pairs * 2.

Heating board terminals:

- Terminals 1 and 2: PT100 sensor connection
- Terminals 3 and 4: output connection
- Terminals 5 and 6: 220VAC power supply connection

SBS-MR-180

A- Heating plate

B- Container

C- Motor


Disposing of used devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual, and packaging. The plastics used to construct the device can be recycled following their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.


Parts diagram

SBS-MR-170

Part number	Description
A	Heating plate
B	Temperature adjustment knob
C	Rotational speed adjustment knob

SBS-MR-180

Part number	Description
A	Temperature adjustment knob
B	Rotational speed adjustment knob
C	Heating Panel
D	Liquid Crystal Display (LCD)

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