USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

3D SHAKER

DE	Produktname	3D-Schüttler
EN	Product name	3D Shaker
PL	Nazwa produktu	Wstrząsarka 3D
CZ	Název výrobku	3D Shaker
FR	Nom du produit	Shaker 3D
IT	Nome del prodotto	Agitatore 3D
ES	Nombre del producto	Agitador 3D
HU	Termék neve	3D Shaker
DA	Produktnavn	3D Shaker
FI	Tuotteen nimi	3D Shaker
NL	Productnaam	3D-schudder
NO	Produktnavn	3D Shaker
SE	Produktnamn	3D Shaker
PT	Nome do produto	Agitador 3D
SK	Názov produktu	3D Shaker
BG	Име на продукта	3D шейкър
EL	Όνομα προϊόντος	3D Shaker
HR	Naziv proizvoda	3D Shaker
LT	Produkto pavadinimas	3D purtyklė
RO	Numele produsului	Agitator 3D
SL	Ime izdelka	3D stresalnik
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-LVS-200
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobci BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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
This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value
Product name	3D Shaker
Model	SBS-LVS-200
Rated voltage [V~] / frequency [Hz]	100-240 / 50-60
Rated power [W]	5
Tilt angle [°]	0-8
Max. load capacity [kg]	2
Protection class	II
IP	20
Dimensions [width * length * height; mm]	190*210*162
Weight [kg]	1.95

Description

The product is used to achieve the mixing results for different liquids.

The user is liable for any damage resulting from unintended use of the device.

Operation

Unpacking

Unpack the instrument carefully and check for any damages which may happen during transport. If any damage happens, please contact us for technical support.

NOTICE! If there is any apparent damage to the system, please do not connect to the power line.

Preparation

1. Unpack the device and check carefully for any damage during transportation. If yes, please contact the manufacturer/supplier for technical support.
2. Place the instrument on a stable, flat, and horizontal platform.
3. The ideal operation environment temperature is 0-40°C. Do not place the machine in a place where the temperature exceeds 40°C or below 0°C.
4. Make sure that the hands are dry when touching the power cord to avoid the risk of electric shock.
5. Ensure that the power supply of the machine is properly grounded.

Operation instructions

- Power on: Connect the power adapter to the DC power interface.
- Running: Adjust the speed control knob (B, Fig. 1) clockwise or counterclockwise to set the speed. Left to right: 0-100 rpm.
- Angle adjustment: Adjust the angle adjustment knob (C, Fig. 1) clockwise or counterclockwise in the stopped state to the required angle (0-8° adjustable), and adjust the speed control knob (B, Fig. 1) to run the machine.

Cleaning and maintenance

- Proper maintenance can keep the instrument working properly and lengthen its lifetime.
- Do not spray cleanser into the instrument when cleaning.
- Do not remove the power line when cleaning. Only use recommended cleansers:

Dyes	Isopropyl alcohol
Construction materials	Water containing tenside / Isopropyl alcohol
Cosmetics	Water containing tenside / Isopropyl alcohol
Foodstuffs	Water containing tenside
Fuels	Water containing tenside

- For cleansers not listed, please consult the manufacturer's after-sales service department first.
- Before using other methods for cleaning or decontamination, the user must ascertain with the manufacturer that this method will not damage the instrument.
- Wear the proper protective gloves when cleaning the instrument.

NOTICE!

- Electronic devices cannot be cleaned with cleanser.
- If maintenance service is required, the instrument must be cleaned first to avoid pollution of hazardous substances.
- If the instrument will not be used for a long time, please switch it off and place it in a dry, clean, room temperature, and stable location.

Disposing of used devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual, and packaging. The plastics used to construct the device can be recycled following their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

Troubleshooting

Problem	Possible cause	Solution
When switched ON, the instrument does not work.	Check whether the power cable is connected.	Connect the power cable.
	Check whether the fuse is broken.	Replace the broken fuse with a new one of the same type and rating.
The running speed can't reach the set value.	Check whether the load of the shaker exceeds the max. load capacity.	Reduce the load on the shaker to ensure it does not exceed the specified maximum load capacity.
Unbalanced forces on the working platform.	Check whether the load of the shaker exceeds the max. load capacity.	Reduce the load on the shaker to ensure it does not exceed the specified maximum load capacity.
	Check whether the loading is placed in the centre of the platform.	Placed in the centre of the platform.

Parts diagram


Fig. 1

Part number	Description
A	Platform
B	Speed control knob
C	Angle adjustment knob
D	DC power interface

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