USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Bruksanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

WATER BATH

DE	Produktname	Thermostatisches Wasserbad
EN	Product name	Water bath
PL	Nazwa produktu	Laboratoryjna łaźnia wodna
CZ	Název výrobku	Laboratoryjna łaźnia wodna
FR	Nom du produit	Bain thermostaté
IT	Nome del prodotto	Bagnomaria termostatico
ES	Nombre del producto	Baño termostático
HU	Termék neve	Termosztatikus vízfűdő
DA	Produktnavn	Termostatisk vandbad
FI	Tuotteen nimi	Termostaattinen vesihaude
NL	Productnaam	Thermostatisch waterbad
NO	Produktnavn	Termostatisk vannbad
SE	Produktnamn	Termostatiskt vattenbad
PT	Nome do produto	Banho-maria termostático
SK	Názov produktu	Termostatický vodný kúpeľ
BG	Име на продукта	Термостатна водна баня
EL	Όνομα προϊόντος	Θερμοστατικό μπεν μαρί
HR	Naziv proizvoda	Termostatska vodena kupelj
LT	Produkto pavadinimas	Termostatinis vandens vonelė
RO	Numele produsului	Baie de apă termostatică
SL	Ime izdelka	Termostatska vodna kopelc
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-TWB-300
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobcu BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value
Product name	Water bath
Model	SBS-TWB-300
Voltage [V] / Frequency [Hz]	200-240~ / 50-60
Power [W]	1000
Reservoir capacity [L]	13
Dimensions (Width x Length x Height) [mm]	350 x 465 x 425
Weight [kg]	8.45
Protection class	I

Product overview

The product is used to provide gentle, controlled heating for delicate processes. It ensures even heat distribution, preventing direct exposure to high temperatures that could cause overheating or damage. It is commonly used in cooking, laboratories, and industrial applications, a water bath helps maintain stable temperatures for tasks that require precision and care.

The user is liable for any damage resulting from unintended use of the product.

Installation

Assembly of the Equipment and Components

The water bath consists of three main components: the bath, the lid, and the base tray. The base tray is positioned inside the bath with its feet downward to create a gap between the tray and the bottom of the chamber.

**NOTE**

- For optimal temperature stability, make sure the base tray and sample container do not contact the sides of the chamber during operation.
- Always use the handle to lift the lid, as other parts may become hot.
- Keep the vent/thermometer hole unsealed to prevent pressure buildup inside the bath.

Instrument Placement

Position the water bath on a level, non-combustible surface, and make sure that both the mains plug and switch are easily accessible.

Electrical Supply

- Verify that the supply voltage and the type of mains plug match your mains outlet, which must be grounded.
- To disconnect the equipment, remove the mains plug from the outlet.

Proper Usage**Avoid Operating the Bath Without Water**

- Always ensure there is water in the chamber before operating the bath. Running the bath dry can shorten the lifespan of key components.
- The bath is equipped with a dry-start protection mechanism that typically detects this condition and prevents further heating.

In case of a dry start, the bath will display a warning and sound an alarm. If this protection program fails to detect the dry heating, a built-in safety temperature cut-out switch will activate and stop the heating when the bath reaches a certain temperature.

**NOTE**

- The dry-start protection program is enabled by default. To turn it off, refer to section *Setting the Dry-Start Protection Alarm*.
- If a dry start occurs, the chamber's internal surface may become very hot, even if the safety cut-out has worked.
- Avoid touching the chamber until it has cooled down.

Filling the Baths


For safe operation, do not fill the bath with water while it is heating. If additional water is needed, switch off the power before filling. If adding water triggers a safety warning during heating, restart the machine to clear the warning.

Emptying the Baths

Before emptying the bath, disconnect it from the power supply and remove the plug. Allow the water to cool to a safe temperature and take precautions to prevent accidental spillage.


The drainage operation is as follows:

1. Insert the drain pipe into the drain port, following the direction indicated in the figure below.

- 1- Drain pipe
- 2- Drain port

2. After draining, loosen the lock of the drain port by turning it in the direction of the arrow, as shown in the figure below. Then, quickly pull out the drain pipe.

Using Transparent Lid


- When the set temperature is above 60°C, use the lid to maintain accurate temperature control and help the water reach the set point.
- The lid also reduces excessive evaporation, minimizing the need to refill the bath and saving energy.

Using Base Tray

- The base tray improves temperature stability for containers placed inside the bath.
- Do not place containers or other items directly in the stainless-steel chamber and avoid direct contact between containers and the chamber's inner walls.

Operating instructions

Bath

- 1- Preset program
- 2- Set
- 3- °C indicator
- 4- Heating indicator
- 5- Alarm indicator
- 6- °F indicator
- 7- Run/Stop
- 8- Countdown
- 9- Over-temperature alarm key
- 10- °C/°F switch

Run/Stop

When the display shows " OFF," the water bath is in an unheated status. When the actual temperature is shown on the display, it indicates that the water bath is in heating mode.

Press  to switch between heating and stopping the heating process.

Example:

- Current temperature: **50.0°C**.
- To start heating, press the **Run/Stop** key. The display will show the - - - - alongside the current temperature.

Operation	Display	
	<i>OFF</i>	The bath is in unheated status
	→ - - - - &50.0	Press the run/stop key to start heating. The current temperature will be displayed, countdown is off.
	→ OFF	Press the run/stop key again to stop heating. The screen will display "OFF."

Heating temperature setting

Method 1

Directly press   to enter the temperature setting mode. Follow the steps below:

Example:

The current temperature is 37.0°C, and set the desired temperature to 50.0°C.

Operation	Display	
 	→ 50.0	Directly set the desired heating temperature.
	→ 37.0	Press the set key to save the value and exit the setting mode. The display will then show the current temperature or "OFF."

Method 2

Press **SET** to enter the temperature setting sub-menu and adjust the temperature by following these steps:

The current temperature is 37.0°C, and set the desired temperature to 50.0°C.

Operation	→	Display	
SET	→	0 C	Press the set key to enter the first-level menu.
▲ ▼	→	TEMP	Navigate to the temperature setting submenu.
SET	→	50.0	Press the set key again to enter the temperature setting option.
▲ ▼	→	55.0	Set the desired temperature.
SET	→	37.0	Press the set key to save the setting and exit. The display will then show the current temperature or "OFF."

Method 3:

Run the preset **P1** **P2** **P3**. Please refer to the detailed steps outlined in the section "Configuring and Running Temperature Presets" for instructions on how to configure and run the temperature presets.

Configuring and Running Temperature Presets

The bath includes three temperature presets **P1** **P2** **P3**, allowing you to conveniently store and access settings for routine use. Before using a preset, you need to save your desired temperature value in one of the presets.

To set a preset program, follow these steps. For example, to configure preset P1:

1. Ensure the current temperature is displayed (e.g., 37.0°C).
2. Set the desired preset temperature (e.g., 50.0°C).
3. Save the setting to preset P1 for future use.

These steps allow you to quickly recall preset temperatures for efficient operation.

Operation	→	Display	
P1	→	P _r 1	Select the preset you wish to set.
▲ ▼	→	50.0	Set the value you wish to use.
P1	→	37.0	Press the preset button again to save the value. The bath will then automatically return to displaying the current temperature or "OFF."

Only use distilled or deionized water to fill the bath. The use of tap water is not recommended.

To run a preset, follow these steps. For example, to run preset P1:

Operation	→	Display	
P1	→	P _r 1	Press the preset you want to use;
SET	→	50.0	Press "Set" to confirm the preset. The display will then show the current temperature or "OFF."

Switching Between Fahrenheit (°F) and Celsius (°C)

When the °C indicator is on, the temperature is displayed in Celsius (°C). When the °F indicator is on, the temperature is displayed in Fahrenheit (°F).

To switch between units, follow these steps:

For example, to change from °C to °F:


Setting the Over-Temperature Alarm

Protect your sample by setting a maximum temperature limit for the bath. This limit prevents the bath from heating beyond a specified temperature.

To set the over-temperature alarm, follow these steps:

For example, if the current temperature is 37.0°C, you might set the over-temperature alarm limit to 90.0°C.

If the liquid temperature exceeds the set over-temperature alarm limit, the water bath will stop heating, display a warning message, and sound an alarm. You can turn off the alarm by following one of these methods:

- **Power Off:** Turn off the power to stop the alarm. Restart the machine once the water temperature falls below the over-temperature limit.
- **Add Cold Water:** Add cold water to lower the temperature below the alarm limit, which will turn off the alarm.
- **Adjust Alarm Limit:** Increase the over-temperature alarm limit to above the current water temperature to stop the alarm.

Setting the Countdown Timer

The countdown timer can be set from 0 to 99 hours and 59 minutes. When the countdown ends, an alarm will sound three times, heating will stop, and "OFF" will be displayed.


NOTE

- The countdown can only start when heating is active. If heating stops, the countdown will pause.
- If heating is active when the timer value is saved, the countdown will start immediately.
- If heating is inactive when the timer value is saved, press "Run/Stop" to start the countdown.

- The countdown settings are independent of temperature settings. Consider the time needed to reach the set temperature when configuring the countdown.

To set the countdown timer, follow these steps:

For example: set the countdown time from 0 to 1 hour 10 minutes.

When the countdown reaches 0, the display shows "OFF" along with the current temperature, and the countdown is deactivated. If the countdown is active (not at 0), the display shows the remaining countdown time along with the current temperature.

The first two digits flash, indicating that you can set the countdown hours.

Set countdown hours.

The last two digits flash, indicating that you can set the countdown minutes.

Adjust the minutes as needed.

Set countdown minutes.

Press the designated key  to save the countdown setting. If you do not save, the setting will be invalid. The countdown timer will begin as soon as the setting is saved. If there is no operation for 10 seconds, the setting mode will exit automatically.

Setting the Dry-Start Protection Alarm

The dry-start protection alarm is automatically enabled when the device is powered on. If the system detects that the inner pot is dry, it will stop heating and sound an alarm.

Calibration

Calibration Conditions

To ensure accurate calibration results, the following conditions are recommended:

1. The thermometer used should be at least 10 times more accurate than the water bath display.
2. Perform calibration in a stable ambient temperature (+/-1°C) with minimal air movement.
3. Allow the liquid to reach the set calibration temperature and stabilize for 30 minutes before beginning calibration.
4. Fill the water bath with water to 50 mm above the tray and place the thermometer in the center of the bath, 40 mm above the tray.

Dual Point Calibration

Due to variations in ambient temperature, humidity, and bath design, there may be slight deviations in temperature measurement. To enhance accuracy, each unit is calibrated before shipment, using two calibration points—low and high temperature:


- **Low Temperature Point (LCP):** 5–50°C
- **High Temperature Point (HCP):** Above 55°C

It is recommended to use 40°C for the low temperature calibration point and 80°C for the high temperature calibration point.

Low Temperature Point Calibration Example:

Before calibration: The display shows 40.0°C, while the actual liquid temperature is 42.0°C.


After calibration: The display adjusts to 41.8°C, matching the actual liquid temperature of 42.0°C.

High Temperature Point Calibration Example:

Before calibration: The display shows 95.0°C, while the actual liquid temperature is 98.0°C.

After calibration: The display adjusts to 97.7°C, aligning closely with the actual liquid temperature of 98.0°C.

After completing both high and low temperature point calibrations, restart the device. Allow the bath temperature to stabilize, then verify the calibration accuracy.

Troubleshooting

Error code	Cause	Solution
<i>Err1Sh1</i>	Temperature sensor 1 short-circuit	Check the temperature sensor 1 or motherboard for a short circuit fault.
<i>Err2OPn</i>	Temperature sensor 1 open-circuit	Check if the sensor 1 is well connected.
<i>Err3Sh2</i>	Temperature sensor 2 short-circuit	Check the temperature sensor 2 or motherboard for a short circuit fault.
<i>Err4OPn</i>	Temperature sensor 2 open-circuit	Check if the sensor 2 is well connected.
<i>Err5dry</i>	Dry-start alarm, no water in the chamber or the water level is below the min.	Switch the bath off and refill the bath with water
<i>Err6OtP</i>	Over-temperature Alarm 1. The initial water temperature is above the <i>OtP</i> value 2. The <i>OtP</i> setting value is too low	1. Let water cool. 2. Reset the <i>OtP</i> value
<i>Err7rOn</i>	Temperature sensors malfunction	Contact supplier or manufacturer.
<i>Err8Out</i>	Calibration temperature not within the required range	1. Check if the low calibration point or high calibration point is within the required range. 2. Check if the temperature difference between the one measured by a thermometer and the one measured on the water bath is more than 10°C (If it is, there is issue on the temperature probe on the water bath)

Maintenance and Cleaning

Proper use of the instrument and regular cleaning will help maintain its condition and extend its service life.

1. **Power Disconnection:** Disconnect the power supply before performing any maintenance or cleaning.
2. **Cleaning:** Routine maintenance is not required, aside from regular cleaning. Wipe the exterior of the equipment with a damp cloth moistened with soap, dishwashing liquid, or alcohol. Allow it to dry completely before resuming use.
3. **Internal Components:** There are no user-serviceable parts inside the unit.


Disposing of Used Devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to

construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

Parts

- ① Transparent lid
- ② Stainless-steel chamber
- ③ LED display
- ④ Water bath
- ⑤ Base tray
- ⑥ Drain port
- ⑦ Handle
- ⑧ USB port
- ⑨ Power switch
- ⑩ Power socket

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