## USER MANUAL

Bedienungsanleitung  
Instrukcja obsługi  
Návod k použití  
Manuel d'utilisation  
Istruzioni per l'uso Manual  
de instrucciones Használati  
útmutató Brugsanvisning  
Käyttöohje  
Gebruiksaanwijzing  
Bruksanvisning Instruções  
de utilização Používateľská  
príručka Ръководство за  
употреба Οδηγίες χρήσης  
Upute za uporabu  
Naudojimo instrukcija  
Manual de utilizare  
Navodila za uporabo

---

---

# AUTOCLAVE

<b>DE</b>	Produktname	Autoklav
<b>EN</b>	Product name	Autoclave
<b>PL</b>	Nazwa produktu	Autoklaw
<b>CZ</b>	Název výrobku	Autokláv
<b>FR</b>	Nom du produit	Autoclave
<b>IT</b>	Nome del prodotto	Autoclave
<b>ES</b>	Nombre del producto	Autoclave
<b>HU</b>	Termék neve	Autokláv
<b>DA</b>	Produktnavn	Autoklav
<b>FI</b>	Tuotteen nimi	Autoklaavi
<b>NL</b>	Productnaam	Autoclaaf
<b>NO</b>	Produktnavn	Autoklav
<b>SE</b>	Produktnamn	Autoklav
<b>PT</b>	Nome do produto	Autoclave
<b>SK</b>	Názov produktu	Autokláv
<b>BG</b>	Име на продукта	Автоклав
<b>EL</b>	Όνομα προϊόντος	Αυτόκλαβο
<b>HR</b>	Naziv proizvoda	Autoklav
<b>LT</b>	Produkto pavadinimas	Autoklavas
<b>RO</b>	Numele produsului	Autoclavă
<b>SL</b>	Ime izdelka	Avtoklav
<b>DE</b> Modell   <b>EN</b> Product model   <b>PL</b> Model produktu   <b>CZ</b> Model výrobku   <b>FR</b> Modèle   <b>IT</b> Modello   <b>ES</b> Modelo   <b>HU</b> Modell   <b>DA</b> Model   <b>FI</b> Tuotteen malli   <b>NL</b> Productmodel   <b>NO</b> Produktmodell   <b>SE</b> Produktmodell   <b>PT</b> Modelo do produto   <b>SK</b> Model   <b>BG</b> Модел на продукт   <b>EL</b> Μοντέλο προϊόντος   <b>HR</b> Model proizvoda   <b>LT</b> : Gaminio modelis   <b>RO</b> : Model de produs   <b>SL</b> : Model izdelka		SBS-AC-2300 SBS-AC-1800 SBS-AC-1200
<b>DE</b> Hersteller   <b>EN</b> Manufacturer   <b>PL</b> Producent   <b>CZ</b> Výrobce   <b>FR</b> Fabricant   <b>IT</b> Produttore   <b>ES</b> Fabricante   <b>HU</b> Termelő   <b>DA</b> Producent   <b>FI</b> Valmistaja   <b>NL</b> Producent   <b>NO</b> Produsent   <b>SE</b> Tillverkare   <b>PT</b> Fabricante   <b>SK</b> Výrobca   <b>BG</b> Производител   <b>EL</b> Κατασκευαστής   <b>HR</b> Proizvođač   <b>LT</b> Gamintojas   <b>RO</b> Producător   <b>SL</b> Proizvajalec		expondo Polska sp. z o.o. sp. k.
<b>DE</b> Anschrift des Herstellers   <b>EN</b> Manufacturer Address   <b>PL</b> Adres producenta   <b>CZ</b> Adresa výrobce   <b>FR</b> Adresse du fabricant   <b>IT</b> Indirizzo del produttore   <b>ES</b> Dirección del fabricante   <b>HU</b> A gyártó címe   <b>DA</b> Producentens adresse   <b>FI</b> Valmistajan osoite   <b>NL</b> Adres producent   <b>NO</b> Produsentens adresse   <b>SE</b> Tillverkarens adress   <b>PT</b> Endereço do fabricante   <b>SK</b> Adresa výrobcu   <b>BG</b> Адрес на производителя   <b>EL</b> : Διεύθυνση κατασκευαστή   <b>HR</b> Adresa proizvođača   <b>LT</b> Gamintojo adresas   <b>RO</b> Adresa producătorului   <b>SL</b> Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra   Poland, EU

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This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via [info@expondo.com](mailto:info@expondo.com).

## Technical Data

Description of the parameter	Value of the parameter		
Product name	Autoclave		
Model	<b>SBS-AC-2300</b>	<b>SBS-AC-1800</b>	<b>SBS-AC-1200</b>
Supply voltage [V~] / Frequency [Hz]	230/50		
Rated power [W].	2300	1800	1800
Current intensity [A]	10	8	8
Safety class	I		
IP code	IPX0		
Dimensions [Width x Depth x Height; mm]	700 x 440 x 410	600 x 440 x 410	600 x 440 x 410
Tray dimensions [width x length; mm]	225 x 395 + 210 x 395 + 200 x 395	225 x 295 + 210 x 295 + 200 x 295	175 x 305 + 150 x 305 (x2)
Chamber dimensions [diameter x depth; mm]	250 x 450	250 x 340	200 x 360
Maximum pressure on trays [kg/cm <sup>2</sup> ]	6	5	4
Weight [kg]	46	42	40
Capacity [l]	22.1	16.7	11.3
Minimum / maximum working pressure [Bar]	-0.9 – 2.3		

Maximum temperature generated [°C]	145	
Sterilization temperature [°C]	121/134	
Residual humidity [%]	0.26	0.36
Water tank capacity [L]	3 (distilled water) + 1.5 (waste water)	
Application temperature [°C]	0-55	
Load capacity of a single tray [g]	1000	
Type of fuse	T15AL 250V	

## 1. General Description

The instruction manual is intended to assist in safe and reliable use. The product is designed and manufactured strictly according to technical specifications using the latest technology and components and maintaining the highest quality standards.

**PLEASE CAREFULLY READ AND UNDERSTAND THIS  
INSTRUCTION MANUAL BEFORE OPERATION**

To ensure long and reliable operation of the unit, make sure to operate and maintain it properly in accordance with the guidelines in this instruction manual. The technical data and specifications contained in this instruction manual are up to date. The manufacturer reserves the right to make changes in order to improve the quality. Taking the technical progress and the possibility of reducing noise into account, the unit is designed and built in such a way so that risks resulting from noise emissions are reduced to the lowest possible level.

### Explanation of symbols

	The product complies with applicable safety standards.
---	--

	Please read the instructions before use.
	Recyclable product.
	CAUTION! or WARNING! or REMINDER! describing a situation. (general warning sign).
	Wear protective gloves.
	CAUTION! Warning of electric shock!
	Caution! Hot surface can cause burns!
	For indoor use only.

**CAUTION! The illustrations in this instruction manual are for reference only and may differ from the actual product in some details.**

## 2. Safety of use

 **CAUTION!** Read all safety warnings and all instructions. Failure to follow the warnings and instructions may result in electric shock, fire and/or severe personal injury or death.

The term "unit" or "product" in the warnings and in the description of the instructions refers to the Autoclave.

---

The product is out of the scope of Medical Device Regulation (EU) 2017/745. It is not a medical device and should not be used for examinations and tests in hospitals and medical centers.

## 2.1. Electrical safety

- a) The plug of this unit must fit into the outlet. Do not modify the plug in any way. Original plugs and matching outlets reduce the risk of electric shock.
- b) Do not touch the unit with wet or damp hands.
- c) Do not use the cord in an unintended manner. Never use it to carry the unit or to pull the plug out of the socket. Keep the cord away from heat sources, oil, sharp edges or moving parts. Damaged or tangled cords increase the risk of electric shock.
- d) If you cannot avoid using the unit in a wet environment, use a residual current unit (RCD). Using an RCD reduces the risk of electric shock.
- e) Do not use the unit if the power cord is damaged or shows signs of wear. A damaged power cord should be replaced by a qualified electrician or the manufacturer's service department.
- f) To avoid electric shock, do not immerse the cable, plug, or unit itself in water or other liquid. Do not use the unit on wet surfaces.
- g) Do not use the unit in rooms with very high humidity / in the immediate vicinity of water tanks!

## 2.2. Safety in the workplace

- a) Keep the work area tidy and well lit. Disorder or poor lighting can lead to accidents. Be foresighted, watch what you are doing and use common sense when using the unit.
- b) Do not use the unit in an explosive area, for example in the presence of flammable liquids, gases or dust. The unit produces sparks that can ignite dust or fumes.
- c) If you find any damage or irregularities in the operation of the unit, immediately turn it off and report it to an authorized person.

- 
- d) If you have any doubts as to whether the unit is working properly or if it is damaged, contact the manufacturer's service department.
  - e) Only the manufacturer's service department can repair the unit. Do not carry out repairs yourself!
  - f) In case of open flames or fire, use only dry powder or snow (CO<sub>2</sub>) fire extinguishers to extinguish the live equipment.
  - g) No children or unauthorized persons are allowed in the work area. (Inattention may result in loss of control of the unit.)
  - h) Use the unit in a well-ventilated area.
  - i) Check the condition of the safety stickers regularly. Replace them if they are illegible.
  - j) Keep these instructions for use for future reference. If the unit is to be passed on to a third party, the operating instructions must also be handed over together with the unit.
  - k) Keep the packaging and small assembly parts out of the reach of children.
  - l) Keep the unit away from children and animals.
  - m) When using this unit together with other units, also follow the other instructions for use.

Please note! Keep children and other bystanders safe while operating the unit.

### 2.3. Personal safety

- a) Do not operate this unit if you are tired, ill or under the influence of alcohol, drugs or medication that could impair your ability to operate the unit.
- b) The unit may be operated by persons who are physically fit, capable of operating it and appropriately trained, and who have read this instruction manual and have been trained in occupational safety and health.
- c) Use caution and common sense when operating this unit. A moment's inattention during operation may result in serious personal injury.
- d) Use personal protective equipment as required when operating the unit as specified in Section 1 of the explanation of symbols.

---

The use of appropriate, approved personal protective equipment reduces the risk of injury.

- e) To prevent accidental start-up, make sure the switch is in the off position before connecting to a power source.
- f) Do not overestimate your capabilities. Maintain body balance and equilibrium at all times during operation. This allows for better control of the machine in unexpected situations.
- g) Pressurized air can cause serious injuries.
- h) The unit is not a toy. Children should be watched to ensure that they do not play with the unit.
- i) Do not place your hands or any objects inside the running unit!
- j) Do not open the dispenser door while the dispenser is in operation or the display shows "6Kpa".
- k) Do not lean over or put your hands on the water tank while the unit is in operation.
- l) The trays must only be removed after sterilization using the special handle.

## 2.4. Safe use of the unit

- a) Do not use the unit if the ON/OFF switch does not function properly (does not turn on and off). Units that cannot be controlled by the switch are unsafe, cannot operate, and must be repaired.
- b) Unplug the unit before making adjustments, changing accessories, or putting it away. This precaution reduces the risk of accidental start-up.
- c) Disconnect the unit from the power supply before adjusting, cleaning, or servicing. This precaution reduces the risk of accidental start-up.
- d) Keep unused equipment out of the reach of children and out of the reach of anyone unfamiliar with the unit or this instruction manual. These units are dangerous in the hands of inexperienced users.
- e) Keep the unit in good working condition. Inspect before each use for general damage or damage to moving parts. If damaged, have the unit repaired before use.
- f) Keep the unit out of the reach of children.

- 
- g) Repairs and maintenance should be carried out by qualified personnel using only original spare parts. This will ensure the safety of use.
  - h) To ensure the designed operational integrity of the unit, do not remove factory-installed covers or loosen screws.
  - i) When transporting or moving the unit from storage to the place of use, observe the health and safety rules for manual handling applicable in the country where the unit is used.
  - j) Do not move, shift, or rotate the machine while in operation.
  - k) Clean the unit regularly to prevent permanent dirt build-up.
  - l) The unit is not a toy. Cleaning and maintenance must not be performed by children without adult supervision.
  - m) Do not tamper with the unit to alter its performance or design.
  - n) Keep the unit away from sources of fire and heat.
  - o) Do not exceed the maximum permissible operating pressure of the unit!
  - p) Do not block the ventilation openings of the unit!
  - q) Use only ultra-pure distilled / demineralized water for operation!  
Use of other water may affect the efficiency of the unit or even lead to its failure.

**CAUTION!** Although the product has been designed to be safe, with adequate safeguards, and despite the additional safety features provided to the user, there is still a slight risk of accident or injury when handling the unit. You are advised to use caution and common sense when using this product.

### 3. Rules of use

The unit is designed to sterilize packaged or unpackaged, solid and empty objects. Toxic or volatile substances must not be sterilized!

**IMPORTANT!** When sterilizing unpackaged items, sterilization only makes sense if they are used immediately after sterilization or to prevent cross-infection.


**NOTICE!** The device is:

- 
- NOT intended for disease diagnosis, prevention, monitoring, treatment, or alleviation.
  - NOT for diagnosis or compensation of injury or disability.
  - NOT intended for anatomical investigations.
  - NOT for in vitro examination of human specimens.
  - NOT for conception control or support.
  - NOT for sterilizing reusable items like surgical and dental instruments, dressings, and glassware.

**The user is responsible for any damage resulting from misuse.**

### 3.1. Description

Explanatory drawings (models SBS-AC-2300, SBS-AC-1800, SBS-AC-1200 similar to each other)

1. Display
2. Distilled water tank
3. Door with lock
4. Printer
5. Power switch
6. Safety pneumatic valve (240 KPa)
7. Air filter
8. Rating plate
9. Distilled water drainage coupling
10. Fuses
11. Socket for power cord

---

12. Compartment with tray holders

13. Trays

14. Door gasket

15. Waste water drainage connector

16. Graduated cup

17. Door gasket

18. Drain hose

19. Tray holder



## 3.2. Preparation for operation

### **POSITIONING OF THE UNIT**

The ambient temperature must not exceed the range of 0-55°C and ambient humidity should not exceed 85%. Place the unit in a way that ensures good air circulation. Maintain a minimum clearance of 10 cm from any wall of the unit, and 20 cm from the top. Keep the unit away from any hot surfaces. Always operate the unit on a level, stable, clean, fireproof and dry surface and out of the reach of children and persons of impaired mental, sensory and intellectual functions. Place the unit in such a way that the mains plug can be reached at any time. Ensure that the power supply to the unit corresponds to that specified on the identification plate!

### **ASSEMBLY OF THE UNIT**

After setting the unit in the designated location, connect the supplied power plug to the unit. Plug the other end of the plug into a grounded power outlet. Place the paper roll in the printer (use a flat, blunt tool, e.g. a screwdriver, to open it) - see the pictures below:

---


Place the paper roll so that it unrolls from the top of the roll laid horizontally. The paper should pass over the black guide in the printer dock cover.

If you need to use the machine immediately after installation, pour high quality distilled/demineralized water into the tank (see further in the manual for water quality requirements).

**IMPORTANT!** The temperature of the poured distilled / demineralized water must not exceed 40°C!

### 3.3. Working with the unit


#### 3.3.1 Control panel

**Display:**

- Clock
- Date
- Pressure in the chamber (PRESSURE)
- Temperature in the chamber (TEMPERATURE)
  - Current operating mode (MODE)

- 
- Current status (STATE)
    - Total number of running cycles (CYCLE)
    - Heating ring / evaporator temperature (OUT-T / STE-T)
  - Distilled water tank and main tank level (WATER1 / WATER2)
    - Door lock status (DOOR)
    - Error code (ERROR)

**Buttons:**

- a) Menu / parameter setting button
- b) Button to scroll up an item/parameter
- c) Button to scroll down the item/parameter
- d) Button for accepting the parameter and starting or stopping the unit

**3.3.2 Switching on the unit**

After connecting the unit to the power supply turn on the power by setting the power switch (5) to the "I" position (the switch will light up).

When the power is turned on, the unit will indicate the last sterilization mode performed and the alarm "Please open the door before working" will sound if the door is closed before the power is turned on.

Open the door and you can operate the unit further, i.e. go to the desired mode.

---

### 3.3.3 Changing the settings

- A) Turning the unit on Then press the parameter setting menu button (a) on the control panel to access the following mode setting submenu:

S E T T I N G	
MODE:	
UNWRAPPED,	121°C
UNWRAPPED,	134°C
WRAPPED,	121°C
WRAPPED (B&D)	134°C
COTTON,	121°C
PLASTIC,	121°C

- "UNWRAPPED..." - disposable vacuum sterilization of unwrapped items; selectable at 120°C or 134°C; operating pressure 110 KPa or 210 KPa; duration of sterilization cycle at 121°C - 20 minutes, while at 134°C - 4 minutes; drying time at 121°C - 5 minutes, while at 134°C - 12 minutes.
- "WRAPPED..." - triple vacuum sterilization of packaged items; choice of temperature 121°C or 134°C; working pressure 110 KPa or 210 KPa; duration of sterilization cycle at 121°C - 20 minutes and at 134°C - 4 minutes; drying time at 121°C - 5 minutes, at 134°C - 12 minutes.
- "WRAPPED (B&D) - control sterilization according to the Bowie and Dick method. This is an automatic program testing the efficiency of the unit, conducted at 134°C.

NOTE: for this program a special test package according to Bowie and Dick method is needed. Such a package is not included in the unit!

- 
- "COTTON" - triple sterilization of cotton materials; working temperature 121°C; working pressure 110 KPa; duration of the cycles 15 minutes; drying time 15 minutes.
  - "PLASTIC" - single vacuum sterilization of plastic items; working temperature 121°C ; working pressure 110 Kpa, duration of the sterilization cycle 20 minutes; drying time 15 minutes.

Select a program by pressing the scroll buttons up / ( b or c) down, then press the accept button (d) to confirm and start the program.

The unit has a built-in mini printer which prints out the sterilization result.

- B) Turn on the unit by setting the power switch (5) to position "I". Then, press the parameter settings menu button (a) on the control panel to select and enter the following unit settings submenu:

Ster. T:	04
Dry. T:	10
V-COUNT:	03
BACKLIGHT:	M
LANGUAGE:	ENG
Clock:	2019-01-20 13:11:06
SETPRINT:	printer

Use the up/down scroll buttons (b / c) to select the parameter in question, then confirm it by pressing the parameter acceptance button (d) to change its value using the up/down scroll buttons (b / c) and after selecting the desired value, confirm it again by pressing the acceptance button (d).

- "Ster.T" - time of the sterilization process

- "Dry.T" - time of the drying process
- "V-COUNT" - number of vacuum processes
- "BACKLIGHT" - intensity of the backlight
- "LANGUAGE" - menu language selection
- "Clock: ..." - date and clock setting
- "SETPRINT" - printer on/off

### 3.3.4 Adding distilled/demineralized water

If after switching the unit on and selecting the sterilization mode and starting it up the message "EMPTY" appears on the display next to the WATER2 position it means that the distilled/demineralized water in the tank is low i.e. there is not enough water to carry out sterilization and therefore it should be added.

If the water has not been topped up the sterilization process will not be executed - for safety reasons the program will not start.

Pour distilled/ demineralized water into the water tank located on top of the unit (2).

To do so, remove the tank lid (2) and pour in water until the message "FULL" appears on the display next to the WATER2 position.

**ATTENTION!** If during the ongoing sterilization process the display shows "EMPTY", don't worry because there is enough water for the ongoing process.

It is only before the next sterilization process that it is necessary to top it up if there is not enough water for the next sterilization process. Requirements for distilled/demineralized water:

	Permissible level in water supplied to the generator/chamber	Allowable level in the condensed steam at sterilizer inlet
Evaporated residues	≤10 mg/l	-
Silicon dioxide (SiO <sub>2</sub> )	≤1 mg/l	≤0.1 mg/l
Iron (Fe)	≤0.2 mg/l	0.1 mg/l
Cadmium (Cd)	≤0.005 mg/l	≤0.005 mg/l

<b>Lead (Pb)</b>	≤0.05 mg/l	≤0.05 mg/l
<b>Heavy metal residues except iron, cadmium and lead</b>	≤0.1 mg/l	≤0.1 mg/l
<b>Chlorine (Cl)</b>	≤2 mg/l	≤0.1 mg/l
<b>Phosphorus pentoxide (P<sub>2</sub>O<sub>5</sub>)</b>	≤0.5 mg/l	≤0.1 mg/l
<b>Conductivity (at 25°C)</b>	3 to 15 µs/cm	≤3 µs/cm
<b>pH-value (degree of acidification)</b>	5-7.5	5-7
<b>Hardness (Σ alkaline earth ions)</b>	≤0.02 mmol/l	≤0.02 mmol/l
<b>Appearance</b>	Colourless, clear, without impurities	

**CAUTION!** Using distilled water which does not meet the above parameters will cause malfunction of the unit and may even lead to its damage!

Use only deionized water with a maximum conductivity of 15 µs/cm - higher values may lead to unit failure.

There is special "ultra pure" water available on the market dedicated to sterilization applications - the use of such water is recommended.

### 3.3.5 Emptying the waste water tank

If during the operation of the unit the message "FULL" appears on the display next to the WATER1 position, this means that the waste water tank is full and should be emptied. To do so, connect a hose to the waste water drainage connector (15). The waste water will be drained when it reaches a maximum temperature of 70°C - if it is higher than this check that the fan is working correctly.

**CAUTION!** Hot water - risk of scalding.

### 3.3.6 Placing items on the trays


Open the door (3) and using the tray holder (19) remove the respective tray to place the items to be sterilized on it.

**CAUTION!** Only pull out/place trays using the supplied tray handle as they may be hot - risk of burns!

---

**IMPORTANT!**

- When sterilizing items of different materials at the same time, they should be grouped on a given tray according to the same material i.e. do not place different materials on the same tray.
- There should be free space between items in the same tray - they must not touch each other.
- Any folding instrument, etc. should be in the open position (e.g. scissors)
- Do not overload the trays.
- Items to be sterilized should be clean and dry.
- When the ambient temperature is below 10°C the unit should be switched on 5-10 minutes before sterilization to warm up.
- Method of placing items on trays:

✓ - correctly

X - incorrect

After placing the trays with the items in the unit, close the door tightly (with the lock knob in the "OPEN" position) and lock it by turning the lock knob all the way in the direction of the "CLOSE" arrow.

**NOTE!** When the chamber inside is still warm, the door may be harder to close - this is normal and is caused by the accumulation of residual steam inside.

Alternatively, you can wait a little longer for the steam to escape from the chamber.

**IMPORTANT!** For safety reasons the unit will not start the sterilization process if the door is not closed properly - the message "OPEN" will appear and an alarm will sound.

---

To deactivate the alarm turn off the power (switch 5), close the door and restart the unit.

### 3.3.7 Start up and sterilization

After making sure that there is enough water in the unit and that the door is tightly closed, press the menu button (a), select the desired program and then press the accept/activate button (d) - the unit will start and the current status of the process will be shown on the display.

Depending on the selected program, objects and temperature, the whole process takes about 40-60 minutes.

The process is as follows:

- Heating (display shows "HEAT") - after starting the program, the unit heats the chamber to the appropriate temperature and maintains it automatically.
- Vacuum generation ("VACUUM" displayed) - air is pumped out of the chamber and steam is injected.

This stage lasts 3 cycles.

- Depressurizing (on the display: "DOWN") - under high pressure, air is evacuated from the chamber.
- Heating (on display: "UP") - the unit generates the correct temperature and pressure.
- Sterilization (on display: "STER") - sterilization takes place and the display indicates the current sterilization temperature and counts down its time.
- Pressure release (on the display: "DRY") - The pressure is released from the chamber. The display indicates the temperature and the time to complete this step. Once the pressure and temperature have been released, the unit enters the vacuum drying stage.
- Completion (on the display: "END") - an alarm will sound to indicate that the sterilization process is complete.

**CAUTION!** Once sterilization is complete do not open the door until the display shows a pressure below 6 Kpa! Remove trays only using the handle provided. It is best to let sterilized items cool down completely before storing or reusing them.

---

After sterilization there will always be some hot steam left in the chamber, so take care when opening the door - do not lean directly over it or hold your hand over the ajar door.

When sterilization is complete the unit will print out a receipt showing a summary of the sterilization process (pressure, temperature, duration, etc. for a given sterilization step, etc.), then it will switch to standby mode, i.e. holding the temperature until the next sterilization process begins.

### **3.3.8 Switching off**

When sterilization is finished and the unit is no longer used, switch it off by setting the power switch (5) to the off position - "O" - its illumination will switch off.

Leave the door open. In case of not using the unit for a longer period of time, additionally disconnect the unit plug from the power outlet.

### **3.3.9 Interruption of unit operation and emergency shutdown**

When the unit is in operation, pressing and holding the accept/start/stop button (d) for 5 seconds will cause the error message "E99" to appear in the display and the unit will switch to a non-standard operating mode. Press the accept/activate/stop button (d) and the vacuum pump will then start a 3 minute drying mode to remove the water from the chamber.

The unit will remind you to open the door. After that the door can be opened, but remember not to open it until the pressure on the display has fallen below 6Kpa!

To start a new sterilization program open and close the door again and select the desired program.

NOTE: if the unit has to be switched off suddenly, turn the power switch (5) to the off "O" position.

### **3.3.10 Interruption of power supply**

In the course of sterilization if there is a break in power supply, i.e. a temporary power failure, then after the power supply is restored the display will show the last performed sterilization stage and the error code "E98".

---

The acceptance/activation/stop button (d) must then be pressed on the control panel for the vacuum pump to start a 3-minute drying mode to remove water from the chamber.

The unit will remind you to open the door. After that the door can be opened, but remember not to open it until the pressure on the display has fallen below 6Kpa!

To start a new sterilization program open and close the door again and select the desired program.

### 3.3.11 Transport/Storage

For transport the unit must be completely cooled down and drained of liquids. Close the door to prevent it from opening unintentionally. The door should be ajar when storing the unit when it is not in use.

#### Error codes

The following error codes can appear in the display of the unit under ERROR in the event of malfunctions:

Error code	Meaning
E0	Problem with the temperature sensor at the chamber surface
E1	Temperature inside the chamber is exceeded
E2	Evaporator temperature is too high
E3	Pressure and temperature cannot be matched
E4	Problem maintaining the sterilization temperature and pressure.
E5	Door open during sterilization
E6	Working an overly long time
E7	Pressure cannot be released


## 3.4. Cleaning and maintenance

- a) Pull the mains plug and let the unit cool down completely before cleaning, adjusting or replacing accessories and when the unit is not in use.
- b) Use only non-corrosive cleaning agents for cleaning the surfaces.
- c) After each cleaning, all the parts should be dried well before the unit is used again.
- d) Store the unit in a dry and cool place protected from moisture and direct sunlight.


- e) Do not spray the unit with a stream of water or immerse it in water.
- f) Make sure that no water enters through the ventilation openings in the casing.
- g) Clean the ventilation openings with a brush and compressed air.
- h) Perform regular inspections of the unit checking technical fitness and any damages.
- i) Use a soft cloth for cleaning.
- j) Do not use sharp and/or metal objects (e.g. a wire brush or metal spatula) for cleaning, as these may damage the surface of the material from which the unit is made.
- k) Do not clean the unit with acidic substances, medical products, thinners, fuel, oil or other chemicals as this may damage the unit.
- l) To transport the unit, cool the unit completely and empty both water tanks.

#### 3.4.1 Maintenance interval

Type of service	Tips:
<b>Daily</b>	
Cleaning the door seal	Water with a little soap. After washing, wipe thoroughly with a damp, soft lint-free cloth.
Cleaning the chamber	When the chamber has cooled down completely, wipe it and the door from the inside with a damp, soft lint-free cloth.
<b>Weekly</b>	
Cleaning the chamber, trays and rack	When removed from the machine, they should be cleaned with a special medical disinfectant and then wiped thoroughly with a damp, lint-free cloth.
Cleaning the clean water drain filter	The drain filter is located at the bottom of the chamber (see photos below). Unscrew it and take apart the strainer and adapter. Cleaning it in an ultrasonic cleaner is recommended.

<b>Every month</b>	
Tank cleaning	Empty the distilled/ demineralized water tank and wipe the inside thoroughly dry and sanitize.
<b>Every year</b>	
Perform a bump test	Bowie and Dick (B&D) method - special test kit required.
<b>If needed</b>	
Cleaning/replacing the door seal	When the unit is unplugged and cold: insert a pointed but not sharp object (e.g. clean screwdriver or dull knife) between the gasket and the door panel and gently lever and with the other hand grasp and carefully remove the gasket. If it shows signs of wear (e.g. cracks/tears), replace with a new one; if not, wash it thoroughly, dry it and reattach it by sliding it evenly over the ring. The same tool can be used to slide it on as when removing it.

	 <p>The images show a person's hands cleaning a blue circular component from a washing machine drum. Step 1: The component is held in the hand. Step 2: The component is being scrubbed with a brush. Step 3: The component is being rinsed with water. Step 4: The component is being dried with a cloth.</p>
Additional cleaning	

### DISPOSAL OF USED UNITS

At the end of its useful life, this product should not be disposed of with normal household waste but should be taken to a collection point for the recycling of electrical and electronic equipment. This is indicated by the symbol on the product, operating instructions or packaging. The materials used in this unit are recyclable according to their marking. You will be making an important contribution to protecting our environment by reusing, recycling or otherwise disposing of used units.

Your local administration will provide you with information about the appropriate disposal point for used units.

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