## USER MANUAL

Bedienungsanleitung  
Instrukcja obsługi  
Návod k použití  
Manuel d'utilisation  
Istruzioni per l'uso  
Manual de instrucciones  
Használati útmutató  
Brugsanvisning  
Käyttöohje  
Gebruiksaanwijzing  
Bruksanvisning  
Instruções de utilização  
Používateľská príručka  
Ръководство за употреба  
Οδηγίες χρήσης  
Upute za uporabu  
Naudojimo instrukcija  
Manual de utilizare  
Navodila za uporabo

---

---

# HOMOGENIZER

<b>DE</b>	Produktname	Homogenisierer
<b>EN</b>	Product name	Homogenizer
<b>PL</b>	Nazwa produktu	Homogenizator
<b>CZ</b>	Název výrobku	Homogenizér
<b>FR</b>	Nom du produit	Homogénéisateur
<b>IT</b>	Nome del prodotto	Omogeneizzatore
<b>ES</b>	Nombre del producto	Homogeneizador
<b>HU</b>	Termék neve	Homogenizátor
<b>DA</b>	Produktnavn	Homogenisator
<b>FI</b>	Tuotteen nimi	Homogenisaattori
<b>NL</b>	Productnaam	Homogenisator
<b>NO</b>	Produktnavn	Homogenisator
<b>SE</b>	Produktnamn	Homogenisator
<b>PT</b>	Nome do produto	Homogeneizador
<b>SK</b>	Názov produktu	Homogenizátor
<b>BG</b>	Име на продукта	Хомогенизатор
<b>EL</b>	Όνομα προϊόντος	Ομογενοποιητής
<b>HR</b>	Naziv proizvoda	Homogenizator
<b>LT</b>	Produktu pavadinimas	Homogenizatorius
<b>RO</b>	Numele produsului	Homogenizator
<b>SL</b>	Ime izdelka	Homogenizator
<b>DE</b> Modell   <b>EN</b> Product model   <b>PL</b> Model produktu   <b>CZ</b> Model výrobku   <b>FR</b> Modèle   <b>IT</b> Modello   <b>ES</b> Modelo   <b>HU</b> Modell   <b>DA</b> Model   <b>FI</b> Tuotteen malli   <b>NL</b> Productmodel   <b>NO</b> Produktmodell   <b>SE</b> Produktmodell   <b>PT</b> Modelo do produto   <b>SK</b> Model   <b>BG</b> Модел на продукт   <b>EL</b> Μοντέλο προϊόντος   <b>HR</b> Model proizvoda   <b>LT</b> : Gaminio modelis   <b>RO</b> : Model de produs   <b>SL</b> : Model izdelka		SBS-UH-160 SBS-UH-170
<b>DE</b> Hersteller   <b>EN</b> Manufacturer   <b>PL</b> Producent   <b>CZ</b> Výrobce   <b>FR</b> Fabricant   <b>IT</b> Produttore   <b>ES</b> Fabricante   <b>HU</b> Termelő   <b>DA</b> Producent   <b>FI</b> Valmistaja   <b>NL</b> Producent   <b>NO</b> Produsent   <b>SE</b> Tillverkare   <b>PT</b> Fabricante   <b>SK</b> Výrobca   <b>BG</b> Производител   <b>EL</b> Κατασκευαστής   <b>HR</b> Proizvođač   <b>LT</b> Gamintojas   <b>RO</b> Producător   <b>SL</b> Proizvajalec		expondo Polska sp. z o.o. sp. k.
<b>DE</b> Anschrift des Herstellers   <b>EN</b> Manufacturer Address   <b>PL</b> Adres producenta   <b>CZ</b> Adresa výrobce   <b>FR</b> Adresse du fabricant   <b>IT</b> Indirizzo del produttore   <b>ES</b> Dirección del fabricante   <b>HU</b> A gyártó címe   <b>DA</b> Producentens adresse   <b>FI</b> Valmistajan osoite   <b>NL</b> Adres producent   <b>NO</b> Produsentens adresse   <b>SE</b> Tillverkarens adress   <b>PT</b> Endereço do fabricante   <b>SK</b> Adresa výrobci   <b>BG</b> Адрес на производителя   <b>EL</b> : Διεύθυνση κατασκευαστή   <b>HR</b> Adresa proizvođača   <b>LT</b> Gamintojo adresas   <b>RO</b> Adresa producătorului   <b>SL</b> Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra   Poland, EU

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*This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via [info@expondo.com](mailto:info@expondo.com).*

## Technical data

Parameter description	Parameter value	
Product name	Homogenizer	
Model	SBS-UH-160	SBS-UH-170
Rated voltage [V~]   Frequency [Hz]	220   50/60	
Rated power [W]	510	510
Speed range [rpm]	2000-11000	2000-11000
Working head diameter [mm]	70	700
Mixing Capacity [ml]	500-40000	500-40000
Max. processing viscosity [mPa.s]	8000	8000
Dimensions [Width x Depth x Height; mm]	400 x 355 x 710	350 x 215 x 735
Weight [kg]	13.6	10.05

**Product overview**

The device is driven by a compact high-speed micro-electric motor and is designed for laboratory use. It is ideal for mixing low to medium viscosity liquids, high-viscosity liquids, and for dispersing particles in shear emulsifying liquid media.

The machine consists of a high-speed motor, shear emulsification working head, speed controller, and other key components. The drive unit is constructed from high-quality, die-cast aluminum to ensure stable operation. The working head, made from refined stainless steel, is designed for a wide range of experimental applications.

The operating speed is regulated by a microprocessor, which compensates for sudden speed changes due to overload or voltage fluctuations. This ensures consistent performance. A built-in soft-start function enhances safety, and the stepless speed regulation allows you to adjust the speed as needed.

The digital speed display on the SBS-UH-170 model helps you collect accurate experimental data, while the integrated timer function enables precise control of operating time during experiments.

**Installation**

1. Place the base on the platform. Connect the two supporting stands by screwing them together, then fix them into the mounting hole on the base. Insert the supporting and fixing frame at the upper end of the supporting stand and turn the handle to secure it.
2. Attach the short beam to the emulsifying electric motor through the mounting bolt hole at the back of the motor. Insert the beam into the supporting and fixing frame, then turn the handle to fasten the motor to the supporting stand. The host motor can be freely adjusted using the supporting and fixing frame.

**Note**

- Before testing or operating the product, carefully read the "Use and Maintenance" manual.
- The machine must be installed on a bearing during operation. Handheld use is prohibited.

**Caution**

- Do not use the product in flammable or explosive environments.
- Ensure proper grounding of the input power supply socket. Do not connect the grounding wire to telephone lines, metal pipes, or lightning arresters.
- Never start the machine before placing the shear emulsifying head in the liquid medium. Ensure that the working head (feed port) is kept at a safe distance from the bottom of the container.

---

## Use and maintenance

1. The machine uses a single-phase, three-wire system, so the power socket must be a triple outlet (220V, 10A).
2. Never run the shear emulsification working head without liquid medium to prevent mechanical damage to the components.
3. Input power to the unit:
  - Turn on the power switch.
  - The built-in switch indicator will light up, and the Speed Digital display will show [0000].
  - The "SBS-UH-160" dynamic indicators confirm the machine is ready for operation.
4. **Speed setting:**
  - Rotate the speed control knob to enter the preset speed mode.
  - You can pre-set the desired stirring speed or dynamically adjust the speed during operation.
5. **Start control:**
  - Press the speed control knob to start the machine.
  - For the "SBS-UH-170" model, if no pre-set speed is configured, the machine will operate at the default speed.
6. **Pause/Restart:**
  - Press the speed control knob to pause the machine during operation.
  - Press it again to resume operation at the previously set speed (speed memory function).
7. **Time control:**
  - For the "SBS-UH-170" model, press and hold the speed control knob for 3 seconds until the display shows "SET."
  - Rotate the knob to set the experimental time, then press it again to exit timer mode and return to speed setting mode.
  - Rotate the knob to select the desired speed, then press it to start the machine.
  - The machine will operate at the selected speed, and the display will alternate between the running speed and the countdown timer.
  - If the machine is suspended, the preset operation will need to be reset.
8. **Operating state control:**
  - If power is interrupted or the machine is accidentally switched off, it will return to the preset initial state for safety when restarted.
9. **Operational failure:**
  - If a command cannot be executed, the display will show "Err1" and the "CAUTION" indicator will light up, signaling an alarm.
10. **After operation:**

- 
- Turn off the power, unplug the machine, and clean the shear emulsification working head before the next use.
- 

### **11. Stator selection for experiments:**

**11-1.** For mixing and emulsifying low-viscosity liquids, or for smashing and dissolving small particles, use the small-hole stator.

**11-2.** For medium-viscosity fluids and secondary particle materials, use the large-hole stator.

**11-3.** For high-viscosity fluids and large-particle emulsifying experiments, use the straight-hole stator.

---

### **12. How to change the stators:**

1. Unscrew the three nuts using the enclosed wrench.
  2. Remove the screws, then detach the stator coil from the stator base on the head.
  3. Replace the stator with the correct type for your experiment.
  4. Reassemble the stator coil and screws securely, then tighten the three nuts.
  5. Rotate the drive shaft by hand to ensure smooth operation.
- 

### **13. Maintenance:**

- After prolonged use or inactivity, a maintenance inspection is required before resuming operation.
- Use a megger to measure insulation resistance to ensure electrical safety.

### **14. Brush replacement:**

- The mixer brush is a consumable part and requires regular inspection.
  - Turn off the power, unplug the machine, unscrew the brush cap, and remove the brush.
  - If the brush is less than 6 mm, replace it with a new one of the same specification.
  - Ensure the new brush moves freely in the tube to avoid electrical sparks or motor failure.
- 

### **15. Motor failure troubleshooting:**

- If the motor fails, check the power supply, ensure the plug is secure, and confirm the brush is in good condition.
- If these issues are resolved but the motor still fails, contact after-sales support for assistance.
- Do not disassemble the machine without proper authorization to avoid accidents.

## **Precautions**

---

**⚠ Caution**

- Store the product in a dry and clean environment. The working temperature should not exceed 40°C, and prevent foreign objects from entering the device.
- If operating the product in wet conditions, ensure it is equipped with electric leakage protection.
- Voltage fluctuations may cause unstable speeds. Use a regulated power supply device to ensure consistent performance.
- Avoid using the machine in highly corrosive environments. If use in such conditions is unavoidable, ensure proper protective measures are taken.
- If the agitating medium is volatile or flammable, be cautious of potential electric sparks, which could lead to explosions or fires.

**Shear emulsifying principle**

The experimental device is driven by a high-speed motor with a shear emulsifying head. The rotor, working in alignment with the first stator, spins at high speed to create a circular tangent of high-frequency, strong velocity, and angular motion. This generates integrated kinetic energy, promoting the movement of experimental material from the bottom of the container into the rotor area.

With the resulting centrifugal force, the material is expelled through the stator. In the narrow gap between the stator and rotor, strong, reciprocating hydraulic shear and friction are produced. The experimental material undergoes iterative work, with up to tens of thousands of cuts, tears, and impacts, achieving the effect of shear emulsification. This process is about 1,000 times more effective than standard mixing methods.

The laboratory shear emulsifying machine is suitable for biotechnology organizations, sample preparation in pharmaceutical research, enzymatic treatments in the food industry, cosmetics manufacturing, and various applications in the paint and petrochemical industries. It can handle liquid-liquid and liquid-solid phases for different experimental needs.

Different head and stator configurations can be used to meet a wide variety of experimental requirements, making the machine versatile for emulsification, homogenization, suspension of solids, aggregation, and dissolution.

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