USER MANUAL

expondo.de

BEDIENUNGSANLEITUNG
INSTRUKCJA OBSŁUGI
NÁVOD K POUŽITÍ
MANUEL D'UTILISATION
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES

HFI-10

HIGH FREQUENCY INSTRUMENT

DE	3
EN	5
PL	7
CZ	9
FR	11
IT	13
ES	15

PRODUKTNAME	HOCHFREQUENZSTAB
PRODUCT NAME	HIGH FREQUENCY INSTRUMENT
NAZWA PRODUKTU	PRZYRZĄD WYSOKIEJ CZĘSTOTLIWOŚCI
NÁZEV VÝROBKU	VYSOKOFREKVENČNÍ PŘÍSTROJ
NOM DU PRODUIT	ÉLECTROSTIMULATEUR À HAUTES FRÉQUENCES
NOME DEL PRODOTTO	BACCHETTA AD ALTA FREQUENZA PER ELETTROTHERAPIA
NOMBRE DEL PRODUCTO	CABEZAL DE ALTA FRECUENCIA PARA ELECTROTHERAPIA
MODELL	
PRODUCT MODEL	
MODEL PRODUKTU	
MODEL VÝROBKU	HFI-10
MODÈLE	
MODELLO	
MODELO	
IMPORTEUR	
IMPORTER	
IMPORTER	
DOVOZCE	EXPONDO POLSKA SP. Z O.O. SP. K.
IMPORTATEUR	
IMPORTATORE	
IMPORTADOR	
ADRESSE VON IMPORTEUR	
IMPORTER ADDRESS	
ADRES IMPORTERA	
ADRESA DOVOZCE	UL. NOWY KISIELIN-INNOWACYJNA 7, 66-002 ZIELONA GÓRA POLAND, EU
ADRESSE DE L'IMPORTATEUR	
INDIRIZZO DELL'IMPORTATORE	
DIRECCIÓN DEL IMPORTADOR	

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LEGEND

Read instructions before use.

GENERAL SAFETY INFORMATION FOR THE USE OF ELECTRICAL DEVICES

To avoid injury from fire or electric shock, please make sure to follow these safety instructions when using the device. Please read the instructions carefully and make sure that you have understood them well. Keep the manual at hand to be able to read it again at any time. Always use grounded power outlets providing an adequate voltage (as indicated on the device). If you have any doubt, let an electrician check that your outlet is properly grounded. Never use a damaged power cable. Do not use the device in damp or wet environment, or if your hands or body are damp or wet. Protect the device from sunlight. Use the device in a protected environment to avoid damaging the equipment or endangering others. Make sure the device can cool down properly and avoid placing it too close to other heat-producing devices. Before cleaning, always unplug the device. Use a soft damp cloth for cleaning. Avoid using detergents and make sure that no liquid enters the device. The inner elements of this device don't require any maintenance from the user. An unauthorized opening of the device may void the guarantee.

SAFETY GUIDELINES

1. Please read and understand these instructions carefully before operating the device! Please observe the safety guidelines carefully in order to prevent any damage due to improper use!
2. Please keep this manual at hand for future reference. If the device is handed over to third parties, this manual must be handed over too.
3. Only use this device indoors and for its intended use.
4. The guarantee does not cover damages resulting from improper use or operation.
5. Before first use, ensure that the voltage and current requirements of the device match the available power supply.
6. This device should not be used by persons (including children) with limited physical, sensory or mental aptitude, or lacking the necessary experience and/or knowledge – unless they are supervised by a qualified person or have received safety instructions from this person.
7. **ELECTRICAL SHOCK HAZARD!** Do not attempt to repair the device yourself. In case of device failure, repair must be done by a qualified professional.
8. Please check the main plug and power cable regularly. If the power cable is damaged, it has to be

replaced by the manufacturer or its customer service or another qualified person, in order to prevent hazards.

9. Prevent damages to the power cable by avoiding squeezing, bending or rubbing on sharp edges and keep it away from hot surfaces and open flames.
10. **CAUTION DANGER TO LIFE!** During cleaning, never immerse the device in water or other liquids. Never put the device under running water or pour other liquids on it. Never pour / dispense liquids over the lamp. If any liquid should get inside the device, immediately unplug it and contact the service team.
11. Don't open the device under any circumstances.
12. **WARNING! USING THIS ELECTROSTIMULATION DEVICE HFI-10 IS NOT RECOMMENDED WITH THE FOLLOWING CONDITIONS:** rosacea, telangiectasia, allergic conditions of the skin, pregnancy (any surgery and treatment), hypertension, swellings, anxiety, heart problems, skin wounds, fever, asthma, nervous disorders, epilepsy, metal implants.
13. **WARNING!** Always use the device on dry skin and hair! The glass electrodes on the device must also be kept dry during use!
14. Never use the device with a damaged glass electrode.
15. Never leave the device and its packaging unattended.

PURPOSE

The HIGH FREQUENCY ELECTROSTIMULATION DEVICE HFI-10 is a cosmetic device used in SPA salons, beauty salons or fitness clubs. It is designed to generate high-frequency therapeutic micro-currents, in order to improve cell metabolism, clean and disinfect the skin, enhance microcirculation, improve the absorption of active principles from cosmetics, as well as to tone up the skin and reduce wrinkles.

In case of damage resulting from an improper use of the device, the user shall assume all liability.

BEFORE FIRST USE

On receipt of the parcel, make sure the packaging is undamaged before opening it. If the packaging is damaged, please contact your transport company and distributor within 3 days. Be as specific as possible when describing the damages. Do not turn the package upside down! When transporting the package, please ensure that it is kept horizontal and stable. The packaging materials (cardboard, plastic straps, polyurethane foam) should be kept at hand, so that the device can be sent back to the service centre, if needed, in the best conditions!

ENVIRONMENT OF USE

The room temperature must not be higher than 45°C and the relative humidity should be under 85%. Keep the device away from hot surfaces. Operate the device on flat, stable, clean, refractory and dry surfaces, out of the reach of children and people with mental disabilities. Make sure that the power outlet is always within reach. The power outlet must be properly grounded and match the technical requirements of the device.

OPERATING PRINCIPLE

1. main switch
2. „spoon” electrode – to treat hard to reach areas
3. „mushroom” electrode – to treat cheeks, forehead, neck, chest, back.
4. “spot” electrode – to treat hard to reach areas
5. “comb” electrode – to treat the scalp

APPLICATION

Before using the device, prepare the skin properly – it should be dry (wet skin may damage the glass electrode). Remove the device from the box WITH DRY HANDS. Put the switch in the off position.

Select the appropriate electrode and insert its metal tip in the device output. Please note that both the electrode and your hands must be completely dry.

Turn on the device by rotating the main knob clockwise, until you reach the desired intensity. Start with the lowest intensity - it can be increased gradually by turning the knob clockwise. Never touch metal objects with the device. Before use, the person on whom the treatment is to be performed should remove all metal items (such as jewellery, etc.).

Place the glass electrode to the skin and adjust electrical intensity. Please note! It is forbidden to use the electrodes on eyelids! Any treatment around the eyes should be done with extreme caution, only at the lowest intensity. You may want to combine the use of this device with skin cosmetics (such as masks): if so, make sure they are oil-based (not water-based) and alcohol-free. The device can produce sparks. To avoid them, always use the lowest electrical intensity before putting the electrode on the skin. Then increase gradually the intensity. Before removing the electrode from the skin, set the intensity back to its minimum and turn off the device.

WARNING! Never insert or remove glass electrodes when the device is on!

TRANSPORTATION AND STORAGE

When transporting the device, avoid shaking it, dropping it and turning it upside down. Store it in a cool, dry place, free of any corrosive gas.

CLEANING AND MAINTENANCE

- Unplug the device before cleaning.
- Use a mild, non-corrosive detergent.
- Let the device dry before using it again.
- Store the device in a dry, cool location, free from moisture and direct sunlight.
- The quartz balls can be washed in a sieve under running water with a mild detergent.

REGULAR CONTROL OF THE DEVICE

Regularly ensure that the device is in good condition. Should it be damaged, please stop using it immediately and contact your customer service to solve the problem.

What to do in case of problem?

Please contact your customer service and prepare the following information:

- Invoice number and serial number (the latter is to be found on the technical plate on the device).
- If relevant, a picture of the damaged, broken or defective part.
- It will be easier for your customer service advisor to identify the source of the problem if you give a detailed and precise description of the matter. The more details you give, the quicker the customer service will be able to solve your problem!

CAUTION: Never open the device without the authorization of your customer service. This may void the guarantee!

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